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EN | Instruction for use
Citralysat K2
Medical device IIB

1. Treatment indication
Citralysat K2 is a dialysis solution designed for continuous methods of blood purification (CRRT – Continuous Renal Replacement Therapy) with regional citrate anticoagulation provided by 4% trisodium citrate and concurrent calcium administration. Citralysat K2 can be used only on devices designated for continuous hemodialysis equipped with a module for citrate anticoagulation and with pumps for citrate and calcium infusion next to standard pumps for blood flow operation, dialysis solution and filtrate.

2. Composition
Citralysat K2 is delivered as a double-chamber bag made of SA (Sealed Air) with a total volume of 5000 ml. Dialysis solution for an immediate use is made by mixing contents of both bag compartments just before use. The smaller chamber contains 250 ml of acid solution with 15 mmol/l Mg²⁺, 40 mmol/l K⁺, 72 mmol/l Cl⁻, with excess of 2 mmol/l (corresponding to pH around 2.7) and 111 mmol/l of glucose. The larger chamber contains 4750 ml of solution with 140 mmol/l Na⁺, 118.84 mmol/l Cl⁻ and 21.16 mmol/l HCO₃⁻.

Dialysis solution for an immediate use contains:	Ions and glucose concentration in the dialysis fluid for immediate use are:
Sodium chloride 6,604 g/l	Na ⁺ 133 mmol/l
Sodium hydrogencarbonate 1,680 g/l	K ⁺ 2,00 mmol/l
Potassium chloride 0,1491 g/l	Mg ²⁺ 0,75 mmol/l
Magnesium chloride hexahydrate 0,1525 g/l	Cl ⁻ 116,50 mmol/l
Glucose monohydrate 1,100 g/l	HCO ₃ ⁻ 20 mmol/l
	Anhydrous glucose 1,0 g/l
	Theoretical osmolality ≥ 278 mosmol/l
	Hydrochloric acid 25%

Other components:
Water for injection
Hydrochloric acid 25%

3. Contraindications
Citralysat K2 must not be given intravenously or infused into extracorporeal circuit. Citrate impaired metabolism like in patients with impaired liver functions may lead to acidosis, hypocalcemia or an in-potentialy need for calcium substitution. In such a case citrate CRRT must be discontinued and another therapy of renal replacement must be applied.

4. Adverse reactions
Efficacious CRRT therapies may bring about hypophosphatemia. Phosphate levels in serum should be therefore monitored before and during the treatment and phosphate depletion should be compensated.

5. Directions for use
User shall be duly instructed on the proper application of Citralysat K2 before its use. Citralysat K2 may only be used by a doctor competent to provide continual replacement of the renal function or under supervision of such doctor. During its use other additional devices and solutions are needed and the solution of trisodium citrate 4% shall be administered via infusion in pre-dilution. The ratio of this 4% solution to that of blood flow under the standard conditions should be 1:34, which corresponds to an infusion of 4 mmol of citrate per 1 liter of blood treated.


Calcium with its concentration within 50 a 500 mmol/l shall be administered via system infusion or via venous CRRT blood set just before connecting to venous catheter's lumens. Calcium amount administered via infusion shall be set appropriately for the concentration of system ionized calcium to be controlled. Appropriate initial rate is routinely 1.7 mmol Ca per 1 liter of applied dialysis solution. If Citralysat K2 is injected at a rate of 2 l/h, it corresponds to 3.4 mmol Ca²⁺ infused in an hour. Metabolic acidobasic state of the patient can be controlled by the ratio between buffer infusion depended on the blood and citrate flows and buffer removals via a flow of dialysis solution. Increased flow of dialysis solution which may lead to acidosis shall be taken into account. This is what distinguishes it from other dialysis solutions used in CRRT combined with system anticoagulation (e.g. 35 mmol/l bicarbonate). If Citralysat K2 is injected at a rate of 2 l/h together with a blood flow rate of 100 ml/min, a routine balanced metabolic state can be achieved. This holds if an adequately big dialyzer not restricting transport of buffers has been used.

6. Continuous hemodialysis dosage
If it is not clinically contraindicated, the required efficacy of hemodialysis treatment is achieved by administering 1.5 - 2.5 l/h of dialysis solution in adults, depending on the body weight (20-25 ml/kg). No experience with treatment in children has been obtained to date.

7. Bag connection to extracorporeal circuit
After mixing up contents of both chambers the bag with the dialysis solution is connected to the dialysate circuit according to instructions for devices designated for continuous renal replacement therapy. It is necessary to avoid contaminations of the dialysis solution and any parts being in contact with the dialysis solution. The bag with Citralysat K2 disconnected from the dialysate circuit have to be disposed of. Before its use Citralysat K2 shall be warmed up closely to the patient's body temperature for it not to be substantially decreased.

- 8. Notice**
Dialysis solution for citrate CRRT in a double-chamber bag.
- Do not use before mixing up both the chambers.
 - Solution is not for intravenous infusion
 - Use only when fluid is clear and package undamaged
 - Dialysis solution must be used only combined with citrate infusion (citrate anticoagulation).
 - Dialysis calcium-free solution: a separate calcium infusion is necessary.
 - Connection of system citrate CRRT must be carefully checked up before the treatment initiation.
 - Incorrect connection of citrate and calcium infusion must be avoided above all. Correct connection should be approved by measuring a decrease of ionized calcium concentration induced with citrate in extracorporeal circuit within 20-30 min after the therapy initiation, at the latest. If such a decrease does not take place, it is necessary to check up the connection, because a changeover of citrate and calcium infusions may lead to an electrolytes serious disbalance.
 - Make use of high-flux dialyzer with minimum area of 1.4 m² a exchange it each 72 hours, at least.
 - High-Flux Dialyzer with minimum 1.4 m² active surface and change it after 72 hours.
 - Free of bacterial endotoxins.
 - Only for single use.


9. Manipulation
9.1 Pulling out of the cardboard box
The double-chamber bag containing the dialyzing solution must be pulled out of the cardboard box immediately before use.

9.2 Opening the second overpouch
Open the second overpouch at the tear notches located on the edges, see the picture. Do not use sharp objects and be careful not to damage the primary packaging.

WHEN TAKING THE BAG OUT OF THE SECONDARY PACKAGING, DO NOT GRIP THE TUBES WITH CONNECTORS.

9.3 Instructions for mixing both chambers
The contents of both chambers must be mixed immediately before use of the dialyzing solution. Before mixing both chambers, check the composition, production batch number and the use by date. Check that the solution is clear. Check the bag and safety outlet stopper are not damaged. The solution must not be used if the packaging or the safety outlet stopper is damaged.

9.4 Tilt the small chamber back. Grasp the small chamber with both hands (preferably with thumbs positioned on the top) and squeeze until the chamber divider (peel seal) rips apart.

9.5 Grasp the corner of the big chamber and bend it upwards. Apply sufficient pressure to rip the divider completely apart.

9.6 Mix the contents of the bag thoroughly.

When the contents of both chambers have been mixed, check that the seam is completely open, the solution is clear and colorless and no fluid is leaking from the bag.

9.4 Tilt the small chamber back. Grasp the small chamber with both hands (preferably with thumbs positioned on the top) and squeeze until the chamber divider (peel seal) rips apart.

9.5 Grasp the corner of the big chamber and bend it upwards. Apply sufficient pressure to rip the divider completely apart.

9.6 Mix the contents of the bag thoroughly.

When the contents of both chambers have been mixed, check that the seam is completely open, the solution is clear and colorless and no fluid is leaking from the bag.


Sample handling with two-chamber bags can be found at www.bio-medica.eu/roztoky_navody

Citralysat K2 must be used within 48 hours after mixing. Any unused product residue must be disposed of.

10. Expiry date and way of storage:
EXP 24 months. Expiry date: see information on the label. Store at the temperature within +4 °C and +30 °C.

11. Packaging disposal information: We recommend to dispose of the packaging as potentially hazardous waste.

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