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### Natrium Citricum 0.5%

(Trisodium Citrate 18 mM)

An anticoagulant solution for continual replacement of the renal function (CRRT - Continuous Renal Replacement Therapy). Medical device IIb

#### 1. Intended use

Natrium Citricum 0.5% is a solution for regional anticoagulation with citrate in the extracorporeal circulation during continual replacement of the renal function.

#### 2. Caution

Strictly follow the instructions for use according to the Package Leaflet of the product.

This solution must be used in predilution on a dialyser designated for continual replacement of the renal function, while the given dialyser must be suitable for anticoagulation with citrate. During the entire treatment the hemodynamic condition of the patient, the liquid balance, as well as the electrolyte and acid-base balance must be duly monitored.

Administration of any liquids, both related to CRRT and out of the CRRT prescription, requires careful assessment. In relation to the composition of the other solutions used during the treatment, a separate infusion containing a sufficient amount of calcium, magnesium or sodium may also be needed.

#### 3. Warning

Do not use for direct intravenous infusion.

Natrium Citricum 0.5% is designated for a single use only. It contains a sterile solution. Repeated use of the product may lead to an infection. To prevent microbiological contamination, it is necessary to use the solution immediately after opening. Do not use Natrium Citricum 0.5% after expiry date printed on the box label and on the back side of the bag.

Use only when the solution is clear and free of any particles. Use only when the overwrap and the bag containing the solution are undamaged.

Any unused solution must be disposed of.

Store at the temperature from +4° C to +30° C.

Provided that the solution needs to be warmed up to the body temperature (+37° C), this procedure must be carefully monitored and it is necessary to check that the solution is clear and free of any particles.

#### 4. Description

Natrium Citricum 0.5% is a single-chamber PP bag containing a sterile anticoagulation ready-for-use solution that is free of any bacterial endotoxins.

#### 5. Composition

One litre of the Natrium Citricum 0,5% solution contains sodium citrate dihydrate Ph. Eur. 5.29 g (0.5% weight by volume) and sodium chloride Ph. Eur. 5.03 g/l

(pH has been set with use of acidum citricum monohydrate Ph. Eur.). Water for injection Ph. Eur. ad 1,000 ml. pH 6.4-7.5

Composition in mmol/l:

|                                                                       |     |        |
|-----------------------------------------------------------------------|-----|--------|
| Citrate (C <sub>6</sub> H <sub>5</sub> O <sub>7</sub> <sup>3-</sup> ) | 18  | mmol/l |
| Natrium (Na <sup>+</sup> )                                            | 140 | mmol/l |
| Chloride (Cl <sup>-</sup> )                                           | 86  | mmol/l |
| Theoretical osmolarity                                                | 244 | mOsm/l |

#### 6. Indications

The anticoagulation solution of Natrium Citricum 0.5% is indicated for patients treated with CRRT in a citrate anticoagulation regime and mainly for patients with a high risk of bleeding or those who suffer from heparin-induced thrombocytopenia.

#### 7. Contraindications

No absolute contraindications exist for Natrium Citricum 0.5%. In clinical situations when the citrate cycle may change, for example upon severe liver failure or in the conditions related with severely reduced perfusion in a muscle, e.g. upon a shock, the benefits and risks of using citrate should be duly assessed by the doctor.

#### 8. Adverse reaction

The solution contains sodium citrate. It does not contain any calcium or magnesium. Using of the solution may cause severe disturbance of the acid-base balance (alkalosis) or severe decrease of the calcium levels (hypocalcaemia) or the magnesium levels (hypomagnesaemia) in blood.

#### 9. Instructions for use (manipulation)

Natrium Citricum 0.5% may only be used by a doctor competent to provide continual replacement of the renal function or under supervision of such doctor.

Citrate anticoagulation is considered as a suitable method of anticoagulation, particularly in patients with the risk of bleeding, and is recommended for all critically ill patients hospitalised at intensive care units.

- When removing the product from the secondary package do not grab it by the hoses with connectors.
- Remove the bag from the package only shortly before use.
- Check the composition, the batch number and the expiry date.
- Check whether the solution is clear and without any visible particles.
- Check the bag and the outlet safety lock whether they are not damaged. The solution must not be used if the package or the outlet safety lock is damaged.
- According to the instructions of the attending doctor add all necessary additives in the bag with use of an injection port. Assessment of compatibility of any additives with the Natrium Citricum 0.5% solution is in charge of the doctor.
- Proceed in compliance with the detailed instructions provided by the dialyser manufacturer.

#### 10. Expiry date and way of storage:

The expiry date of the Natrium Citricum 0.5% (Trisodium Citrate 18 mM) anticoagulation solution is indicated on the bag. The product must not be used after this date.

EXP 24 months

The anticoagulation solution of 0.5% Natrium Citrate (Trisodium Citrate 18 mM) must be protected from light and stored under the temperatures from + 4° C to + 30° C.


#### 11. Packaging:

PP bag of 5,000 ml.

#### 12. Information for package disposal:

We recommend to dispose the package as potentially hazardous waste.

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