

# StrokeSENS

---

v.2.1.0

(300)

Release Notes

## What's New in Version 2.1.0

### New StrokeSENS ASPECTS clinical module

- Alberta Stroke Program Early CT Scoring (ASPECTS) for automated detection and characterization of early ischemic changes on head CT

### New StrokeSENS mCTA clinical module

- Multi-phase CTA (mCTA) tools to produce visualizations that indicate dynamic changes in contrast-enhanced CTA related to parameters of perfusion and cerebral blood flow.

### Platform updates

- Support for Linux Docker container deployment
- Additional controls for licensing and distribution
- REST API service for triggering study processing and alerting host system of processing status
- StrokeSENS is now supported in the following languages:
  - o English (en\_US)
  - o French (fr\_FR)
  - o German (de\_DE)
  - o Italian (it\_IT)
  - o Portuguese (pt\_PT)
  - o Spanish (es\_ES)

## Known Issues

Number	Summary
<b>NV-3144</b>	JSON for ASPECTS has an extra blank entry <i>Workaround: User can still use all other keys/values in the JSON</i>

### In the Build Package

- o StrokeSEN\_2.1.0\_(300).tar.gz
- o 2.1.0
  - Software Usage Terms and Conditions.pdf
  - StrokeSENS v2.1.0 Software Integration Guide.pdf
  - StrokeSENS v2.1.0 DICOM Conformance Statement.pdf
  - UserManuals
    - StrokeSENS\_User\_Manual\_v2.1\_de\_DE.pdf
    - StrokeSENS\_User\_Manual\_v2.1\_en\_US.pdf
    - StrokeSENS\_User\_Manual\_v2.1\_es\_ES.pdf
    - StrokeSENS\_User\_Manual\_v2.1\_fr\_FR.pdf
    - StrokeSENS\_User\_Manual\_v2.1\_it\_IT.pdf
    - StrokeSENS\_User\_Manual\_v2.1\_pt\_PT.pdf

## External Libraries – Copyright and Licensing Requirements

Library name	version	Licensing agreement	Download link
Tensorflow	2.7.0	Apache License 2.0	<a href="https://www.tensorflow.org/">https://www.tensorflow.org/</a>
Qt	5.15.14	<a href="https://www.qt.io/terms-conditions/">https://www.qt.io/terms-conditions/</a>	<a href="https://login.qt.io/login">https://login.qt.io/login</a>
Eigen	3.3.7	Mozilla Public License, version 2.0	<a href="https://bitbucket.org/eigen/eigen/get/3.3.7.zip">https://bitbucket.org/eigen/eigen/get/3.3.7.zip</a>
DCMTK	3.6.7.1	<a href="http://dicom.offis.de/dcmj2k.php">http://dicom.offis.de/dcmj2k.php</a>	<a href="http://dicom.offis.de/dcmTk.php.en">http://dicom.offis.de/dcmTk.php.en</a>
SimpleITK	2.0.2	<a href="https://github.com/SimpleITK/SimpleITK/blob/master/LICENSE">https://github.com/SimpleITK/SimpleITK/blob/master/LICENSE</a>	<a href="https://github.com/SimpleITK/SimpleITK/releases/tag/v2.0.2">https://github.com/SimpleITK/SimpleITK/releases/tag/v2.0.2</a>
Dcm2niix	1.0.20171215	<a href="https://github.com/rordenlab/dcm2niix/blob/master/license.txt">https://github.com/rordenlab/dcm2niix/blob/master/license.txt</a>	<a href="https://github.com/rordenlab/dcm2niix">https://github.com/rordenlab/dcm2niix</a>
GDCM	3.0.21	<a href="https://gdcm.sourceforge.net/Copyright.html">https://gdcm.sourceforge.net/Copyright.html</a>	<a href="https://gdcm.sourceforge.net/wiki/index.php/Downloads">https://gdcm.sourceforge.net/wiki/index.php/Downloads</a>
Boost	1.66.0	<a href="https://www.boost.org/users/license.html">https://www.boost.org/users/license.html</a>	<a href="https://www.boost.org/users/history/version_1_66_0.htm">https://www.boost.org/users/history/version_1_66_0.htm</a>
PyTorch	1.10	<a href="https://pytorch.org/">https://pytorch.org/</a>	<a href="https://github.com/pytorch/pytorch/blob/master/LICENSE">https://github.com/pytorch/pytorch/blob/master/LICENSE</a>
RLM	15.1BL1	<a href="http://www.reprisesoftware.com/admin/software-licensing-downloads.php">http://www.reprisesoftware.com/admin/software-licensing-downloads.php</a>	<a href="http://www.reprisesoftware.com/admin/software-licensing-download.php?&amp;euagree=yes">http://www.reprisesoftware.com/admin/software-licensing-download.php?&amp;euagree=yes</a>
ICU	60.2	<a href="https://github.com/unicode-org/icu/blob/main/icu4c/LICENSE">https://github.com/unicode-org/icu/blob/main/icu4c/LICENSE</a>	<a href="https://icu.unicode.org/download/56">https://icu.unicode.org/download/56</a>
OpenCV	4.5.5	<a href="https://github.com/opencv/opencv/blob/master/LICENSE">https://github.com/opencv/opencv/blob/master/LICENSE</a>	<a href="https://github.com/opencv/opencv/tree/2.4">https://github.com/opencv/opencv/tree/2.4</a>
Go	1.20	<a href="https://golang.org/LICENSE">https://golang.org/LICENSE</a>	<a href="https://golang.org/dl">https://golang.org/dl</a>

## Support

For further information on installing/exchanging the dynamic libraries please contact Circle Cardiovascular Imaging support.

**Phone:** +1 (403) 338-1870

**E-mail:** [support@circlecvi.com](mailto:support@circlecvi.com)

**Website:** [circlecvi.com](http://circlecvi.com)

**Circle Cardiovascular Imaging Inc.**

1800, 707 - 8th Avenue SW

Calgary, Alberta

Canada T2P 1H5