

# StrokeSENS

---

v.1.4.2 (325)

Release Notes

## What's New in Version 1.4.2

### Intended Use

StrokeSENS is a decision-aid software package to be used by clinicians to perform image processing, analysis, viewing and communication of computed tomography (CT) scans of the brain in patients with suspected acute stroke. Data and images are acquired through DICOM-compliant imaging devices prior to processing and analysis in StrokeSENS.

The StrokeSENS software provides analysis capabilities for imaging datasets acquired with standard CT imaging and contrast-enhanced CT Angiography (CTA) modalities. Analysis of non-contrast CT images is provided by the StrokeSENS ASPECTS module, which includes assessment of regions with suspected acute ischemic tissue (ie. ASPECTS scoring). Analysis of contrast-enhanced CT images is provided by the StrokeSENS LVO module, which includes automated detection of suspected large vessel occlusion (LVO).

StrokeSENS ASPECTS is a computer-aided diagnosis software device used to assist the clinician in the assessment and characterization of brain tissue abnormalities using CT image data. The Software automatically registers images and segments and analyzes ASPECTS Regions of Interest (ROIs). StrokeSENS ASPECTS extracts image data for the ROI(s) to provide analysis and computer analytics based on morphological characteristics. The imaging features are then synthesized by an artificial intelligence algorithm into a single ASPECT (Alberta Stroke Program Early CT) Score. StrokeSENS ASPECTS is indicated for evaluation of patients presenting for diagnostic imaging workup with anterior circulation Large Vessel Occlusion (LVO), for evaluation of extent of disease. Extent of disease refers to the number of ASPECTS regions affected which is reflected in the total score. This device provides information that may be useful in the characterization of early ischemic brain tissue injury during image interpretation (within 6 hours). StrokeSENS ASPECTS provides a comparative analysis to the ASPECTS standard of care radiologist assessment using the ASPECTS atlas definitions and atlas display including highlighted ROIs and numerical scoring.

StrokeSENS LVO is a radiological computer-aided triage and notification software indicated for use in the analysis of CTA head images. The device is intended to assist hospital networks and trained clinicians in workflow triage by flagging and communication of suspected positive findings of Large Vessel Occlusion (LVO) in head CTA images.

StrokeSENS LVO uses a software algorithm to identify suspected LVO findings. In the case of a suspected LVO, the system will send a notification to a pre-configured destination(s), notifying the clinicians of the existence of a suspected LVO that requires review. The notification system is intended to be used in parallel to the standard of care workflow to notify clinicians of the existence of the case earlier that they may have been notified as part of the standard of care workflow.

Notifications may include a compressed preview of images. Notifications are meant for informational purposes only and are not intended for diagnostic use beyond notification. The StrokeSENS LVO device does not alter the original medical image and is not intended to be used as a diagnostic device.

The results of StrokeSENS LVO are intended to be used in conjunction with other patient information and based on professional judgement, to assist with triage / prioritization of medical images. Notified clinicians are responsible for viewing full images per standard of care.

### Device Description

StrokeSENS is a software package that provides for the visualization and study of changes in brain tissue using digital images captured by Computed Tomography (CT) diagnostic imaging systems. StrokeSENS is compliant with the DICOM standard and can be installed to run on an off-the-shelf physical or virtual server.

It provides viewing, quantification, analysis, and reporting capabilities which, when interpreted by a trained physician, may yield information useful in clinical decision making.

StrokeSENS works with the following types of (DICOM-compliant) medical image data:

- CT (Computed Tomography)

StrokeSENS acquires (DICOM-compliant) medical image data from the following sources:

- DICOM file
- Network using DICOM protocol

StrokeSENS provides tools for performing the following types of analysis:

- ASPECTS - Ischemic Tissue Changes (ASPECTS)
- LVO - Large Vessel Occlusion detection

StrokeSENS provides tools for performing the following types of analysis under evaluation (Not for clinical use):

- mCTAp - Maps related to cerebral blood flow

## What's New

### 1. Export to PACS

- This functionality enables users to configure StrokeSENS to automatically export processed outputs/results to PACS by configuring export DICOM node(s).




Figure 1.1: DICOM Secondary Captures prepared for export following LVO Detection. In the case that a suspected LVO was detected, it is clearly indicated by the addition of "Suspected LVO Detected!" to the output.






Figure 1.2: DICOM Secondary Captures prepared for export following ASPECTS results. DICOM export includes the original series CT with predicted ASPECTS masks written onto the CT DICOM image, as well as a report image as visualized above.

### 2. Licensing

- Licensing was updated to allow for licensing control of StrokeSENS installation at a institution/installation level as well as for each unique version

- Licenses will expire after a prescribed time
  - Licenses can be renewed/replaced by Circle support following update or upgrade
3. LDAP Configuration for Active Directory group integration
- Admin users can configure a connection to the hospital or health system's Lightweight Directory Access Protocol (LDAP) for Active Directory group integration
  - This enables configuration to allow for end-users to log in to StrokeSENS using the same credentials they would use throughout the hospital or health system
  - LDAP configuration is the responsibility of an ADMIN IT user upon installation
4. ASPECTS Module Updates
- ASPECTS processing module was updated to enable processing of thin slice non-contrast CT's (0.6-5mm)
  - ASPECTS results no longer require confirmation of LVO
5. LVO Module Updates
- The LVO detection algorithm was updated for increased performance
    - Sensitivity of 90.3% (+0.9%) & Specificity of 95.1% (+7.7%) for detecting ICA and M1 occlusions
    - Mean processing time ~2 minutes
6. Multi-phase CT Angiography Perfusion (mCTAp) Updates
- This release includes updates to the mCTA Perfusion module for non-clinical use
  - The release includes updates to the algorithm to improve the performance and resolution of the perfusion maps.
  - The update includes calculation of relative cerebral blood flow (rCBF) and relative Delay (rDelay) maps
  - The update includes calculation and visualization of discrete thresholded maps with pre-defined thresholds and corresponding volumetric analysis
  - mCTAp maps are labelled not for clinical use






Figure 6.1: Multi-phase CTA perfusion outputs. Please user manual for detailed descriptions.

## Resolved Issues

The following list contains resolved usability related issues, in addition to all other internal issues resolved during the product development process.

Issue key	Summary
NV-1755	Heat map gradient showing up in non-brain region within pane
NV-1760	Some series found to separate Internal NCCT
NV-1773	Scaling on MPR incorrect with MIPS and switch to workstation
NV-1774	Switching between heatmaps and normal view 15 times blanks NCCT
NV-1775	mCTAp studies classify 2 times
NV-1779	Patient list doesn't update with new studies for users logged in
NV-1780	Processing Time Criteria for mCTAp not met
NV-1781	Important buttons and info cut off at min width
NV-1788	ML node occasionally crash when processing mCTAp algorithm on release and debug mode
NV-1791	Second row of viewers doesn't fit on screen in Safari
NV-1794	Institution Name tag not visible on study UI by default
NV-1795	Spoke user admin can't modify users to or from spoke admin
NV-1796	No hospital edit/delete confirmation dialog
NV-1797	No 'mCTA failed' messaging on front-end

## Known Issues

Issue key	Summary
NV-999	Swipe does not respond as expected on mobile <i>Workaround: No workaround available.</i>
NV-1267	Manufacturer is not listed on the worklist. <i>Workaround: No workaround available.</i>
NV-1299	Patient Protocol images not carried over using combine studies feature <i>Workaround: No workaround available.</i>
NV-1297	Inconsistent series slice thickness info after combining studies <i>Workaround: No workaround available.</i>
NV-1289	Password field from the LDAP Config view need to be encrypted. <i>Workaround: No workaround available.</i>
NV-1257	Patient List Name Sorting isn't Truly Alphabetical <i>Workaround: No workaround available.</i>
NV-1256	Patient Names Overlap Date in Worklist on Mobile devices. <i>Workaround: No workaround available.</i>
NV-1161	Presence of CTP ML loading icon in newly uploaded studies <i>Workaround: No workaround available.</i>
NV-1619	Older versions of Safari browser might show case high usage of system memory. <i>Workaround: Use updated browser version or other supported browsers.</i>
NV-1620	Switching to heatmaps on workflow page loads extra tiny image first. <i>Workaround: No workaround available.</i>
NV-1623	Loading MPR on 2 viewers sometimes causes the images to be stuck and sometime blank in the viewers. <i>Workaround: No workaround available.</i>
NV-1631	Layout frames column and rows have conflicting dimensions <i>Workaround: No workaround available.</i>
NV-1657	Header disclaimer "ASPECTS results auto-generated..." gets cut off when the browser is shrunk. <i>Workaround: No workaround available.</i>
NV-1665	Inconsistent Modal pop-ups between Hospital management and User Management/Network Management <i>Workaround: No workaround available.</i>
NV-1666	Inconsistent resets when editing hospitals with invalid info <i>Workaround: No workaround available.</i>

	Improvements to password error messaging
NV-1691	<i>Workaround: No workaround available.</i>
	Quick findings icon always highlighted
NV-1723	<i>Workaround: No workaround available.</i>
	In mobile width, using the windowing tool horizontally can swipe series
NV-1782	<i>Workaround: No workaround available.</i>
	Multiview doesn't resize correctly
NV-1785	<i>Workaround: No workaround available.</i>
	Workstation series tab previews don't resize correctly
NV-1784	<i>Workaround: No workaround available.</i>
	Default Zoom Opened Patient Study too far Zoomed In
NV-1871	<i>Workaround: No workaround available.</i>
	No way to access mCTAp when navbar buttons hide due to resize
NV-2016	<i>Workaround: No workaround available.</i>
	Not all series show up in viewer at all times
NV-2038	<i>Workaround: No workaround available.</i>
	Windows preset is not preserved between mCTA phases
NV-2069	<i>Workaround: No workaround available.</i>
	The 'OK' button to continue session does not always work on first click
NV-2123	<i>Workaround: No workaround available.</i>
	Sync Button not present on workstation viewer
NV-2125	<i>Workaround: No workaround available.</i>
	stroke.log_svc performance on Windows 10
NV-2144	<i>Workaround: No workaround available.</i>
	Slice Scrolling Blocked by Series Tab in Workstation
NV-2185	<i>Workaround: No workaround available.</i>
	mCTAp Brain Windowing Option Blocked by the Gradient Scroller for Heatmaps
NV-2211	<i>Workaround: No workaround available.</i>
	Difference in Size between NCCT and Reconstructed NCCT in mCTAp Page
NV-2243	<i>Workaround: No workaround available.</i>
	Zoom level is not retained on switching between normal view and heatmaps
NV-2251	<i>Workaround: No workaround available.</i>
	Dicom Info Button Inaccessible on NCCTs in 2x2 and 3x2 Workstation
NV-2252	<i>Workaround: No workaround available.</i>
NV-2253	MPR buttons style off while Mips is open on Edge browser

	<i>Workaround: User can use a different browser.</i>
NV-2256	Toggle Contours Button Still Accessible when using Mips/Mean <i>Workaround: No workaround available.</i>
NV-2257	Slicesync on mCTAp page doesn't disable for series when enabled on workflow <i>Workaround: No workaround available.</i>
NV-2279	Product name error in "Activation license not available" error message <i>Workaround: No workaround available.</i>
NV-2278	mCTAp heatmap gradient moves with arrow keys <i>Workaround: No workaround available.</i>
NV-2280	User manual switches pages too quickly while zoomed in <i>Workaround: No workaround available.</i>
NV-2287	Windowing Retained in Coronal and Sagittal Views when Switching from Heatmaps to Normal <i>Workaround: No workaround available.</i>
NV-2288	mctap link is disabled on reloading the browser while in patientlist page <i>Workaround: No workaround available.</i>
NV-2318	StrokeSENS performing ML on coronal image and not axial image for certain datasets. <i>Workaround: No workaround available.</i>
NV-2415	StrokeSENS Login / LDAP Solution on Mobile & Tablet is creating multiple (duplicate) users <i>Workaround: No workaround available.</i>
NV-2391	Series do not show in the Patient List according to their SeriesNumber <i>Workaround: No workaround available.</i>
NV-2488	ASPECTS has unhandled crash on 3mm NCCT study <i>Workaround: No workaround available.</i>
NV-2432	License server did not install after upgrade to version 1.4 <i>Workaround: No workaround available.</i>

## External Libraries – Copyright and Licensing Requirements

This is a list of free and open source software used in our products, including any modifications made by Circle Neurovascular Imaging Inc., throughout the course of our development.

**Qt v5.15.3**  
Project Website  
License  
Download

Qt Project  
Qt Commercial Enterprise License  
Source Code and Binaries

### **DCMTK v3.6.7.1**

Project Website  
License(s)

Download

DICOM Toolkit (DCMTK) Project  
DCMTK License  
DCMJP2K Commercial License  
DCMTK v3.6.7 Source Code (.zip)

### **Eigen v3.3.7**

Project Website  
License  
Download

Eigen Project  
Mozilla Public License v2.0  
Eigen v3.3.7 Source Code (.zip)

### **GLEW v1.9.0**

Project Website  
License(s)

Download

The OpenGL Extension Wrangler (GLEW) Project  
The OpenGL Extension Wrangler Library License  
Mesa 3-D Graphics Library License  
The Khronos Group Inc. License  
GLEW v1.9.0 Source Code

### **GLM v0.8.4.4**

Project Website  
Licenses  
Download

OpenGL Mathematics (GLM) Project  
The Happy Bunny and MIT Licenses  
GLM v0.9.3.4 Source Code

### **OpenCV v2.4.3**

Project Website  
License  
Download

Open Source Computer Vision (OpenCV) Project  
OpenCV (3-Clause BSD) License  
OpenCV v2.4.3 Source Code

### **OpenSSL v1.0.2u**

Project Website  
License  
Download

OpenSSL Project  
OpenSSL License  
OpenSSL Source Code

### **Boost v1.66.0**

Project Website  
License  
Download  
Modifications

Boost Org  
Boost Software License  
Boost v1.66.0 CPP Source Code (.zip)

### **QtService v2.8**

Project Website  
License  
Download

Qt Project  
BSD License  
QtService v2.8 CPP Source Code (.zip)

### **QtSingleapplication v2.8**

Project Website  
License  
Download

Qt Project  
BSD License  
QtSingleapplication v2.8 CPP Source Code (.zip)

### **OpenLDAP v2.4.35**

Project Website  
License  
Download

OpenLDAP Project  
The OpenLDAP Public License  
OpenLDAP 2.4.35 Source Code (.tgz)

### **Qt Cryptographic Architecture (QCA) v2.1.0**

Project Website  
License  
Modifications

QCA Project  
GNU Lesser General Public License v2.1  
Disable MD2 if not available in OpenSSL (patch)  
Fix Lookup in Templates (patch)  
Fix Outdated Link to Qt Documentation (patch)

### **Go v1.13.12**

Project Website

[Go Project](#)

License  
Download  
BSD License  
Source Code and Binaries [\[SB\]](#)

**Consul v1.7.3**  
Project Website  
License  
Download  
Consul Project  
Mozilla Public License  
Source Code and Binaries

**React JS v16.14.0**  
Project Website  
License  
[React Project](#)  
MIT License

**Redux JS v4.0.0**  
Project Website  
License  
[Redux JS Project](#)  
MIT License

**HammerJS v2.0.8**  
Project Website  
License  
Download  
HammerJS Project  
MIT License  
HammerJS Source Code

**Pako v1.0.5**  
Project Website  
License  
Download  
Pako Project  
MIT License  
Pako Source Code

**LRU-Cache v3.1.0**  
Project Website  
License  
LRU-Cache Project  
ISC License

**Idle-Session-Timeout v1.0.1**  
Project Website  
License  
Idle-Session-Timeout Project  
MIT License

**Network-Speed v1.2.1**  
Project Website  
License  
Network-speed Project  
MIT License

**service v1.0.0**  
Project Website  
License  
service Project  
Zlib License

**Scheduler v0.0.0**  
Project Website  
License  
Scheduler Project  
MIT License

**Gorilla WebSocket v1.4.1**  
Project Website  
License  
Gorilla WebSocket Project  
BSD 2-Clause License

**mgo v0.0.0**  
Project Website  
License  
mgo Project  
BSD 3-Clause License

**Amqp v0.0.0**  
Project Website  
License  
amqp Project  
BSD 2-Clause License

**Msgpack v2.9.1**  
Project Website  
License  
msgpack Project  
BSD 2-Clause License

**Nginx v1.18.0**  
Project Website  
Nginx

License

BSD 2-Clause License

**RabbitMQ v3.9.10**

Project Website  
License

RabbitMQ  
Mozilla Public License

**MongoDB v3.4**

Project Website  
License

MongoDB Inc.  
Licensed under the AGPL/SSPL

**Metabase v0.36.0**

Project Website  
License

Metabase Project  
GNU Affero General Public License

**Erlang v24.1.17**

Project Website  
License

Erlang  
Apache License 2.0

**Tensorflow v1.13.2**

Project Website  
License

Tensorflow Project  
Apache License 2.0

**Simple ITK v2.0.2**

Project Website  
License

Simple ITK  
Apache 2.0 License

**ANTs v2.3.4**

Project Website  
License

ANT  
ANT License

**dcm2niix v1.0.20171215**

Project Website  
License

dcm2niix  
dcm2niix License

**PyTorch 1.10.0**

Project Website  
License

pytorch  
pytorch License

## Support

For further information on installing/exchanging the dynamic libraries please contact Circle Neurovascular Imaging support.

**Phone:** +1 (403) 338-1870

**E-mail:** [support@circlenvi.com](mailto:support@circlenvi.com)

**Website:** [www.circlenvi.com](http://www.circlenvi.com)

### **Circle Cardiovascular Imaging Inc**

Doing business as,

### **Circle Neurovascular Imaging**

1100, 800 5<sup>th</sup> Avenue SW

Calgary, AB T2P 3T6

Canada