Ultherapy®

P R I M E

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REFERENZENZEN

- Araco A. Prospective Study on Clinical Efficacy and Safety of a Single Session of Microfocused Ultrasound With Visualization for Collagen Regeneration. *Aesthet Surg J.* 2020;40(10):1124-1132. doi:10.1093/asj/sjz363
- Suh DH, Choi JH, Lee SJ, Jeong KH, Song KY, Shin MK. Comparative histometric analysis of the effects of high-intensity focused ultrasound and radiofrequency on skin. *J Cosmet Laser Ther.* 2015;17(5):230-236. doi:10.3109/14764172.2015.1022189
- Fabi SG, Joseph J, Sevi J, Green JB, Peterson JD. Optimizing Patient Outcomes by Customizing Treatment With Microfocused Ultrasound With Visualization: Gold Standard Consensus Guidelines from an Expert Panel. *J Drugs Dermatol.* 2019;18(5):426-432.
- Park JY, Lin F, Suwanchinda A, et al. Customized Treatment Using Microfocused Ultrasound with Visualization for Optimized Patient Outcomes: A Review of Skin-tightening Energy Technologies and a Pan-Asian Adaptation of the Expert Panel's Gold Standard Consensus. *J Clin Aesthet Dermatol.* 2021;14(5):E70-E79.
- Pavicic T, Ballard JR, Bykovskaya T, et al. Microfocused ultrasound with visualization: Consensus on safety and review of energy-based devices. *J Cosmet Dermatol.* 2022;21(2):636-647. doi:10.1111/jocd.14666
- Werschler WP, Werschler PS. Long-term Efficacy of Micro-focused Ultrasound with Visualization for Lifting and Tightening Lax Facial and Neck Skin Using a Customized Vectoring Treatment Method. *J Clin Aesthet Dermatol.* 2016;9(2):27-33.
- White WM, Makin IR, Slayton MH, Barthe PG, Gliklich R. Selective transcutaneous delivery of energy to porcine soft tissues using Intense Ultrasound (IUS). *Lasers Surg Med.* 2008;40(2):67-75. doi:10.1002/lsm.20613
- White WM, Makin IR, Barthe PG, Slayton MH, Gliklich RE. Selective creation of thermal injury zones in the superficial musculoaponeurotic system using intense ultrasound therapy: a new target for noninvasive facial rejuvenation. *Arch Facial Plast Surg.* 2007;9(1):22-29. doi:10.1001/archfaci.9.1.22
- Fabi SG, Massaki A, Eimpunth S, Pogoda J, Goldman MP. Evaluation of micro-focused ultrasound with visualization for lifting, tightening, and wrinkle reduction of the décolletage. *J Am Acad Dermatol.* 2013;69(6):965-971. doi:10.1016/j.jaad.2013.06.045
- Data were collected in a survey of 9 clinics with 68 patients in 2024. Data on file.
- Chang YC, Croix J, Jawaji S, Hernandez S, Chapas AM, MacGregor J. Patient satisfaction and our clinical experience with 459 microfocused ultrasound treatments. *Lasers Surg Med.* 2019;51(6):495-499. doi:10.1002/lsm.23058
- Merz Aesthetics®. Data on File. 2023.
- Dobke MK, Hitchcock T, Misell L, Sasaki GH. Tissue restructuring by energy-based surgical tools. *Clin Plast Surg.* 2012;39(4):399-408. doi:10.1016/j.cps.2012.07.008
- Fabi SG, Goldman MP, Dayan SH, Gold MH, Kilmer SL, Hornfeldt CS. A prospective multicenter pilot study of the safety and efficacy of microfocused ultrasound with visualization for improving lines and wrinkles of the décolleté. *Dermatol Surg.* 2015;41(3):327-335. doi:10.1097/DSS.0000000000000322
- Oni G, Hoxworth R, Teotia S, Brown S, Kenkel JM. Evaluation of a microfocused ultrasound system for improving skin laxity and tightening in the lower face. *Aesthet Surg J.* 2014;34(7):1099-1110. doi:10.1177/1090820X14541956
- Sasaki GH, Abelev N, Papadopoulos L. A Split Face Study to Determine the Significance of Adding Increased Energy and Treatment Levels at the Marionette Folds. *Aesthet Surg J.* 2017;37(8):947-960. doi:10.1093/asj/sjx055
- Guo S, Dipietro LA. Factors affecting wound healing. *J Dent Res.* 2010 Mar;89(3):219-29. doi: 10.1177/0022034509359125. Epub 2010 Feb 5. PMID: 20139336; PMCID: PMC2903966
- Hirayama D, Iida T, Nakase H. The Phagocytic Function of Macrophage-Enforcing Innate Immunity and Tissue Homeostasis. *Int J Mol Sci.* 2017 Dec 29;19(1):92. doi: 10.3390/ijms19010092. PMID: 29286292; PMCID: PMC5796042.
- Kerscher M, Nurrisyanti AT, Eiben-Nielson C, Hartmann S, Lambert-Baumann J. Skin physiology and safety of microfocused ultrasound with visualization for improving skin laxity. *Clin Cosmet Investig Dermatol.* 2019;12:71-79. Published 2019 Jan 14. doi:10.2147/CCID.S188586
- Moretti M. The Global Aesthetic Market Study: XXI. Medical Insight, Inc; 2024.
- Sasaki G, Grossman J, Fabi SG, et al. Stimulation of Collagen Synthesis in Human Skin Following Microfocused Ultrasound Therapy. Presented at the ASDS 2018 Annual Meeting.

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© 2024, Merz Aesthetics (Schweiz) AG. Alle Rechte vorbehalten.
Merz Aesthetics (Schweiz) AG