


# Luer Lock Syringes: A Simple Guide in Easy Steps

**STEP 1:** Syringes are delivered as such; air bubbles are typical and shouldn't be expelled to avoid spillage. Needles are provided separately within this package.



**STEP 2:** Immediately before administration, remove the round cap from the syringe and unwrap the needle. Hold the needle by its plastic cover; avoid touching the open end.



**STEP 3:** Twist the needle securely onto the syringe body. Ensure a tight seal. The pharmacy is not liable for leaks between the syringe and needle interface.



**STEP 4:** This is how your assembled needle and syringe combo should look.



**STEP 5:** Prepare the injection site with the provided alcohol swab. Remove the clear plastic needle cover by pulling directly away and administer the injection. Do not recap used needles. Dispose of the used syringe immediately in a "sharps" or similar container.

