

## WRIST & THUMB BRACE



### BREG (USA) CLASSIC WRIST BRACE

The Classic Wrist Brace is constructed of lightweight, durable material, with a thin web space area for improved patient comfort.

#### Common Examples of Use

- Wrist sprains
- Tendonitis
- Soft tissue injuries
- Post cast support

LEFT#	RIGHT#	SIZE	WRIST CIRC.
100637-110	100637-210	XSmall	5.25" – 6.25" (13-15.8 cm)
100637-120	100637-220	Small	6.25" – 7.25" (15.8-18 cm)
100637-130	100637-230	Medium	7.25" – 8" (18-20 cm)
100637-140	100637-240	Large	8" to 9" (20-22.8 cm)
100637-150	100637-250	X-Large	9" to 10" (22.8-25.4 cm)



### MEDIROYAL (SWEDEN) VENTUS, SHORT WRIST SUPPORT

This model is made from a three-piece design to provide an optimal fit around the wrist. The model has straighter cut pattern to compensate for swollen joints. There is a dynamic seam over the wrist to provide extra stability when the splint is loaded. The splint pocket is made from double woven nylon for maximum durability. The thumb grip is adjustable. This model is supplied with the standard sky-blue short wrist splint. Right or left model.

COLOR	PART#	SIZE	WRIST CIRC.
Black	MR2220	XSmall	13-14 cm
		Small	14-16 cm
Beige	MR2223	Medium	16-18 cm
		Large	18-20 cm
		X-Large	20-22 cm

\* Circumference measurement, wrist  
Indication:  
Wrist instability - Post-fracture treatment - Tendonitis - Carpal tunnel syndrome or rheumatic problems like arthritis



### MEDIROYAL (SWEDEN) NRX, MANEX RADIAL 2.0

This orthosis has been specially designed for radial nerve palsy or drop-hand syndrome but can also be used for other diagnostic groups. It can replace and support weak muscles and support a functional grip. The kit consists of a dorsal wrist part with a thermoplastic insert, finger platform, a thumb piece and a foam strap. The foam strap can replace the finger platform when used as a night splint. The thermoplastic piece can be molded individually over the dorsal part of the hand by using hot water (75-80 degrees C). That can be done either directly in the pocket of the brace or externally and then inserted into the pocket once it's stable. The dorsal splint creates a platform where the finger and thumb parts can be attached for individual extension assist.

#### The new NRX® Manex Radial 2.0 has been improved on several points:

- A thinner and more flexible NRX® material provides a lower profile
- Injection molded hook with excellent performance
- Improved thumb design with a slightly less stretch for more extension support and a full hook receivable surface where the new micro hook can be attached anywhere.
- An improved finger platform design, also with a hook receivable surface allowing attachment of the hook anywhere.

PART	SIZE	WRIST CIRC.
MR2285	X-Small	13-14 cm
	Small	14-16 cm
	Medium	16-18 cm
	Large	18.20 cm

\* Circumference measurement, wrist

#### Indication:

Radial nerve palsy - Drop hand syndrome after CVA -  
Carpal tunnel syndrome - Post traumatic after distal radius fracture or ulnar deviation of the MCP-joints.