

WRIST & THUMB BRACE



MEDIROYAL (SWEDEN) VENTUS POLLUS, SHORT WRIST

This model is made from a three-piece design to provide an optimal fit around the wrist. The short thumb piece provides a good support for the CMC-joint and the thenar eminence. There is a dynamic seam over the wrist to provide extra stability when the splint is loaded. The splint pocket is made from double woven nylon for maximum durability. The thumb grip is adjustable. This model is supplied with the standard wrist splint.

Specification:

COLOR	PART#	SIZE	WRIST CIRC.
Black	MR2222	XSmall	13-14 cm
		Small	14-16 cm
Beige	MR2225	Medium	16-18 cm
		Large	18-20 cm
		X-Large	20-22 cm

* Circumference measurement, wrist

Indication:

Wrist instability and problems from the CMC-joint arthritis or tendonitis.



MEDIROYAL (SWEDEN) RESTO® BASIC RHP

This model has the thumb in a resting position instead of a functional position as for the E25/E26. The inner lining is made from a fixed Neptune™ material added comfort. It's soft and absorbs surface moisture and quickly dries to provide an optimal skin comfort. The complete splint can be machine washed in 40 degrees when using a laundry bag. E27 is fitted with two straps around the wrist, one separate thumb strap and a strap for the finger platform.

Left or Right model.

Specifications

PART#	SIZE	Measurement MCP width	Length of splint
E 27	Small	60-75 mm	30 cm
	Medium	70-85 mm	34 cm
	Large	85-95 mm	37 cm

* Circumference in cm over wrist

Indications

Stroke Or Neurological Diseases Where The Hand Needs Support And Positioning. Can Be Used For Slack Paresis Or Mild To Moderate Spasticity