

FINGER SPLINT



MEDIROYAL (SWEDEN) 502 COIL-EXTENSION SPRING-PIP

This spring provides an extra high degree of extension to the PIP joint. The coil springs are placed on either side of the joint. Pads are placed distal and proximal to the joint to distribute the pressure evenly. The front strap is adjustable with hook and loop. The force over the joint is approx. 12,4 N in full extension.

SIZE	LENGTH
A	64 mm
B	70 mm
C	76 mm

Indication:
Flexionkontraktur i PIP-led

* Length in mm. The measurement is taken from the distal crease in the palm to the crease of the DIP joint. The finger should be extended during the measurement.

MEDIROYAL (SWEDEN) 504 COIL-EXTENSION SPRING-PIP SHORT


- Extends the PIP-Joint but does not affect the MCP-joint.
- Smaller sizes can be used for the DIP-joint as well.
- The spring tension can be adjusted individually.
- The standard force is approx. 2.3 N at 15° of flexion.

SIZE	LENGTH	WIDTH
A	48 mm	22 mm
B	51 mm	25 mm
C	57 mm	29 mm

Spring tension: Medium

* Length in mm. Take the measurement on the volar side of the finger, from the base of the joint.

Indication:
Extension defect of 45 degrees or less in the PIP joint, correction of Boutonniere defect.



PRIMO (M'SIA) FINGER SPLINT WITH STRAPS, Padded-12 per pack

Finger splint is used to immobilize a finger when it is injured. Its main job is to prevent finger from bending and protect it from second injury. It is constructed of malleable aluminum with foam padding. And there can be several kinds of splint to meet different requirements of customers.

- Material: Aluminum and Foam
- Color: Blue
- Application: First-aid Device for Immobilizing Fingers

Features

- Radiolucent
- Good immobilization
- Keep interphalangeal joints in correct position
- Malleable and pliable Good air permeability
- Easy to put on or put off with hook and loop strap

PART#	SIZE	LENGTH IN INCH / CM
AB002-1	Size 1	2.5 inch (5.5cm)
AB002-2	Size 2	3.5 inch (8.2cm)
AB002-3	Size 3	4.25 inch (9.2cm)