

## SHOULDER BRACING



### BREG (USA) SHOULDER ABDUCTION PILLOW

Shoulder Abduction Pillow is a shoulder immobilizer designed for varying degrees of abduction (10°-75°). When pillow is inflated, arm straps can be applied to limit posterior shift of the shoulder following rotator cuff repairs.

#### Features

- inflatable bladder for abduction of 10°-75°
- Limits posterior shift of the shoulder
- Universally sized abduction pillow
- Universal left or right

PART#	DESCRIPTION
01851	Shoulder Abduction Pillow (Universal)

#### Indication:

Large rotator cuff repairs



### BREG (USA) UNIVERSAL ESSENTIAL SHOULDER SLING

The Universal Essential Shoulder Sling is a one-size-fits-most envelope sling.

The material can be folded to fit most patients while maintaining the use of the thumb loop.

#### Features

- Roll-back universal design
- Ultra-soft straps
- Closed pocket design
- Front strap adjustments
- Unique design to prevent D-ring rotation
- Ultra-breathable mesh material

PART#	SIZE
VP20103-000	One Size

#### Indication:

For non-surgical or post-op support

- Mild strains and sprains of the shoulder, elbow and forearm
- Arthroscopic repair
- Shoulder injuries or instabilities



### MEDIROYAL (SWEDEN) ERIXTHREE SHOULDER SUPPORT

The Erix Three® Neuro shoulder is a shoulder brace that has been designed to provide positioning, support and proprioceptive feedback to patients with neurological disorders. The shoulder joint is the largest non weight bearing joint in the body that completely rely on proprioceptive feedback and balance between the muscles and ligaments. Any imbalance can affect the functionality of the joint negatively.

ErixThree® Neuro has been designed to provide the shoulder joint with a posterior safe position. It will retract the shoulder posteriorly in order to increase functionality and to create a better position for the joint. In many neurological disorders the anterior antagonists can become short or dominating when the posterior antagonists like Supraspinatus, Infraspinatus and Teres Minor are not functioning properly. This will cause the shoulder joint to slide anteriorly.

The ErixThree® Neuro shoulder has been specially designed for neurological patients with a smaller volume over the biceps and shoulder joint area due to lack of muscle function. The inside has a special friction pattern over selected muscles to provide better proprioception and biofeedback. These friction surfaces should not be worn on bare skin. A t-shirt under the brace is necessary to not cause any skin friction, especially on neurological patients with more sensitive skin. The posterior part incorporates an area that covers the Scapula along with a chest strap. These areas are flexible in order to retract the shoulder joint posteriorly. The inside of the chest strap has a friction surface with micro perforation that will run under the Axilla and attach on the front side of the brace.

#### Indication:

The product has been designed to provide support to the shoulder joint after neurological disorders, stroke or mild to moderate anterior shoulder subluxations.

PART	SIZE	
NRX940 Left & Right	XSmall	23 - 27 cm
	Small	27 - 32 cm
	Medium	32 - 35 cm
	Large	35 - 38 cm
	X-Large	38 - 42 cm

\* Measurement chart:

\* The size measurement is taken around the biceps.