

COLD THERAPY



BREG (USA) POLAR CARE KODIAK (INTERNATIONAL VERSION)

The Polar Care Kodiak Cold Therapy System is the most convenient and versatile offering in Breg's Polar Care line. Its easy to use, compact design makes it great for clinic, hospital, and home-use. With the addition of a little ice and water, your patients will enjoy 6-8 hours of effortless cold therapy. Offering the only battery powered option on the market, Breg ensures your patients can enjoy the benefits of the cold therapy from anywhere: on the sidelines or in the back yard; the Kodiak will be a reliable companion. Each battery pack comes with 4 replaceable AA batteries ready to power 10-14 hours of motorized cold therapy.

Product Features

- Physicians and patients benefit from reduced pain and swelling-follow-up visit focus on recovery, rather than pain management
- Whisper quiet pump technology allows patients to rest and recover in comfort
- Durable, powerful construction provides peace of mind
- Dual filtration pump system prevents clogs
- Battery Pack option increases compliance for patients on the go
- Intelli-Flo Pad technology meets the unique demands of the treatment area with fine tuned temperature control

PART#	DESCRIPTION
100291-000	Polar Care Kodiak International

IntelliFlo Pads (Intl)	
100293-000	Shoulder, X-Large
100295-000	Knee, Compression
100297-000	Multi-Use, Long Stem
100392-000	Ankle, Compression
100300-000	Back
100302-000	Hip
100303-000	Hand/Wrist

Accessories	
100316-000	INTL Power Supply
97050	Battery Pack
10699	Intelli Flo Pad Y Adaptor

Indication:

- Post Operative
- Arthroscopic Procedures
- Plastic Surgery
- General Surgery
- Post-trauma
- Chronic pain
- Physical therapy

Intelli-Flo Y Adaptor

WARNING

The Polar Care Kodiak can be cold enough to seriously injure the skin. Follow the Operating Instructions (on the back and side of the unit), Product Insert (in the pouch on side of unit), and the Intelli-Flo Pad Fitting Instructions (provided with each Intelli-Flo Pad).

This Y-adaptor allows two BREG Intelli-Flo pads to be applied simultaneously, in combination with a Polar Care Kodiak, for cold therapy applications in multiple treatment areas.

NOTE: The Kodiak system requires enough water to feed both pads. The addition of water may be necessary after initial start-up, due to water displaced into the pads. Once both pads are filled, fill unit with additional ice and water as necessary.

Application Instructions:

- 1 Connect the center coupling of the Intelli-Flo Y-adaptor to the Kodiak system.
- 2 Connect the remaining two couplings to each of the Intelli-Flo pads.
- 3 Ensure pads are cooling properly upon start-up. If necessary add additional ice and water to the Kodiak system.

BREG
an orthoCare Company

BREG Inc.
VIA, CA 92081 USA

Telephone (800) 521-0807
(760) 588-5000
Telex (800) 521-0734

EXP 11/04/05 REV 3