

Makes every bed a safe bed

Delivery, Installation & Training


CloudCuddles that are purchased through NHS, Council, Respite or Charity funding as a long term sleeping solution are delivered and installed with full training by Spacesaverbeds Ltd trained personnel. In some areas of the UK where installation training has been given, those trained personnel will install the CloudCuddle. In some cases where a private purchase is made, delivery is by Royal Mail next day service to all parts of the UK including Islands. An installation guide is downloadable on our website.

www.spacesaverbeds.co.uk/Resources



CE marked Class 1 Medical Device

	Junior	Maxx
In the carry bag dimensions	50cm x 40cm x 30cm	75cm x 53cm x 32cm
Outside Size (L) (W) (H)	210cm x 110cm x 105cm	200cm x 110cm x 126cm
Inside Size (L) (W) (H)	200cm x 90cm x 70cm	200cm x 90cm x 90cm
Maximum Sitting Height	70cm	90cm
Maximum User Weight	By Assessment	By Assessment
Optimum Mattress Size	200cm x 90cm x 16cm-26cm	200cm x 90cm x 16cm-20cm
Weight - Bed Tent	3.7Kg	7.4 Kg
Inflatable Pressure	0.45bar (6.5psi)	0.69bar (6.5psi)
Frame Material	Polyester & Polyurethane	Polyester & Polyurethane
Cover Material	Polyester	Polyester
Sheet Material	100% Cotton	100% Cotton
Pump	Handle - Body - Hose	Handle - Body - Hose
Machine Washable Tent	60 degrees	60 degrees
Machine Washable Sheets	60 degrees	60 degrees
Conformity	CE Marked	CE Marked
Classification	Class 1 Medical Device	Class 1 Medical Device
Overall Carrying Weight	6.8Kg	11Kg
Delivered Box Size	60cm x 43cm x 33cm	80cm x 53cm x 33cm



*The
CloudCuddle
is available
throughout
the UK*

2025/2026 EDITION

Cloud Cuddle®

Junior



Maxx

All New
Maxx



CloudCuddle UK Call us on 01905 347 538
info@spacesaverbeds.co.uk
www.spacesaverbeds.co.uk



Junior

500 in use
across the UK

A multi award winning
bed tent that turns every
bed into a safe bed

CloudCuddle – Safe & calming long term sleeping solution

The CloudCuddle Bed Tent has found a place in the UK as a long term sleeping solution for young children with Special Needs through Paediatric Occupational Therapists and Nurse recommendation. The CloudCuddle has been fitted to all types of domestic single beds as well as electric Hi/Lo beds or used just on a mattress on the floor. By doing this it has created a safe sleeping environment that has allowed a child to relax and unwind in a truly calming space as well as being safe and secure, for peace of mind of parents and carers.

Tough, Robust and Durable – CE Class 1 Medical Devices

Although lightweight to transport, the CloudCuddle is tough, robust and durable once fitted to a bed – this combined with the beds calming effect has created a successful sleeping place for some pretty robust characters. The CloudCuddle meets the very highest specification required and is a CE marked Class 1 Medical Device (documentation available on request).

Machine Washable/Disinfected

If it's in a private or care environment, washing the CloudCuddle couldn't be easier. Firstly the zipped in base sheet can be removed and washed/disinfected separately without taking down the CloudCuddle at 60°. Then if required all of the covers and sheets can be unzipped and washed/disinfected at 60°.



CloudCuddle Junior - Where it's used

For a child with special needs the transition from a traditional cot to their next bed can be complicated and expensive for those providing the bed. The need is to find a safe sleeping environment for an active child who isn't safe in a domestic single bed and often expensive high soft sided beds are prescribed. This is where the Cloud Cuddle had really found a place for children as a great cost saving option for the provider.

CloudCuddle Junior - The facts so far

80% of CloudCuddle Juniors are used by children at home as a long term sleeping solution, following recommendation and funding provisions by Paediatric Occupational Therapists and Nurses through NHS or Local Council funding.

15% are used as an option when visiting relatives, respite, hospital or holiday (often charity funded).

5% are used in Children's Respite Facilities or Children's Hospitals.

From April 2021 to June 2025, 56 NHS Trusts or Local Councils across the UK have provided CloudCuddle Juniors for childrens use at home.

Of the CloudCuddle Juniors prescribed for home use:

60% are for children in the age bracket 3-5 Years

30% are for children in the age bracket 6-10 Years

10% are for adolescents



Maxx

Introducing the CloudCuddle
Maxx - A larger and even more
sturdy version



Maxx - A great option to a high sided cot bed

The great success of the CloudCuddle Junior for children from 3 years to 10 years has identified the need for a larger version of the CloudCuddle for children and adolescents.

The higher internal height of 90cm is enough for taller children and adolescents to comfortably sit up if they want.

Maxx - Fits to virtually all community care beds

In response to this incredible demand, we are proud to introduce the CloudCuddle Maxx – a larger, even more robust version. Designed for use with a variable height bed as an alternative to a cot bed, the Maxx offers greater

flexibility. It can also be fitted to compatible domestic beds, providing a more spacious option than the Junior.

Maxx - Fits to virtually all domestic beds

Like the Junior, the CloudCuddle Maxx has been engineered without any entrapment points, ensuring maximum safety. The Maxx features an increased internal height of 90cm (compared to the Juniors 70cm), allowing an adolescent to sit comfortably when needed.

Maxx - Opening from the inside can be enabled or disabled (DOLS)

HOME ASSESSMENTS are available across all of the UK

