

Northside Hospital Pharmacy: Surface Sampling Log: Biological Safety Cabinet (BSC)

Sample Area: NSF Main Pharmacy

Method: ☐ Paddle ☐ Plate ☐ Swab Media: Manufacturer / Type : _____

Lot#: _____ Expiration Date: _____

Certificate of Analysis (COA) for Media, printed and attached: ☐

Employee Collecting Samples: _____

Date / Time Samples Collected: _____

Employee Reading Results: _____

Employee Documenting Results: _____

Action	Date/Time
Interval #1 30 - 35 deg C - Incubate 48 hours	Start:
	End:
Interval #2 20 - 25 deg C - Incubate additional 5 days	Start:
	End:

ISO 5 Area	Specific Location 1 of each (if applicable)	Interval #1 Growth	CFU Count	Interval #2 Growth	CFU Count	Total CFU Count	Test Results
N/A	Control	Y / N		Y / N			N/A
BSC Hood #1	Top	Y / N		Y / N			Pass / Fail
	Bottom	Y / N		Y / N			Pass / Fail
	Side	Y / N		Y / N			Pass / Fail
	Back	Y / N		Y / N			Pass / Fail
	Front	Y / N		Y / N			Pass / Fail
	Equipment	Y / N		Y / N			Pass / Fail

Corrective Action Summary / Comments:

Surface Sampling Action Levels, per 2023 USP <797>	
ISO Class	Surface Sampling Action Levels (CFU/Device or Swab)
5	>3
7	>5
8	>50

**note all sample retests will be completed using new sampling log*

Revised 1-2025