



DECONTAMINATION
CONSUMABLES
CATALOGUE



THE DECONTAMINATION SPECIALISTS

01223 440 475

iM-Med.com



Equipment

Choice of Steelco CSSD and endoscopy cutting-edge decontamination solutions.



Academy

Delivering decontamination training and compliance for continuous patient safety.



Technical

Specialist CSSD and endoscopy services: guaranteeing compliance, optimising performance.




Consumables

Comprehensive range of high-quality endoscopy and CSSD products and accessories.

CONTENTS

The Complete Decontamination Strategy	>
iM Clean	>
iM Consumables	>
iM Clean Stage 1: Procedure Room	>
iM Consumables Team	>
Single-Use Disposable Valves	>
Air/Water Suction Buttons	>
Process Cleaning Button	>
Bite Blocks	>
Endoscopy Biopsy Valves	>
Air/Water Tubing	>
Irrigation Tubing	>
Auxiliary Water Connectors	>
Consumables	>
iM Clean Stage 2: Cleaning	>
Consumables	>
EPW 100S Manual Cleaning	>
Endoscope Cleaning and Flushing	>
iM Clean Stage 3: Reprocessing	>
Endoscopy Chemistries	>
CSSD Chemistries	>
PAA and H2O2 Monitoring	>
Consumables	>
Connectors	>
Leak Test Connectors	>
Valve House Connectors and Channel Separators	>
Air, Water, Suction, Biopsy, Auxiliary and Elevator	>
Channel Connectors	>
Scope Channel Connectors	>
Steelco One-Time Connection System	>
iM Clean Stage 4: Drying, Storage, Transportation	>
Endoscope Trolleys	>
iM Endopax™ System	>
Consumables	>
iM Equipment	>
Steelco Endoscopy and CSSD Equipment	>
iM Academy	>
iM Technical	>
Orders and Enquiries	>





The Complete Decontamination Strategy

Our focus is on enabling patient safety for healthcare providers, whilst providing real choice throughout the endoscope and surgical instrument decontamination journey.

This is why we provide a complete decontamination strategy.

Time to talk to iM Med.





iM Clean

Comprehensive range of high-quality products and accessories for the decontamination journey of medical devices.

The iM Clean Range has been developed to support all the stages of the safe and compliant procedures required for the decontamination journey of the endoscope.

From use, through the decontamination stages, right back to the procedure room, ensuring the endoscope has been reprocessed to its optimum state for the safe use on the next patient.

iM-Med.com/iM-clean





Comprehensive range of
high-quality endoscopy
and CSSD products and
accessories.



CONSUMABLES

Prioritising infection prevention in endoscopy with a portfolio of solutions for procedural patient care.



CLEAN STAGE 1

iM Med is passionate about patient safety.

With a portfolio of solutions along the patient pathway, the iM Clean Range has been developed to enable healthcare providers the choice to procure innovative products of the highest quality, to support safe patient care.

Arrange a visit by a member of the iM Med decontamination specialist team for a discussion on how we can support your individual requirements.

We can also provide departmental audits, product trials and departmental training.

Call **01223 440 475**
Email **enquiries@iM-Med.com**

1



Eliminate the patient risk from contaminated valves

The intricate design and shape of valves means that ensuring safe reprocessing is extremely difficult and problematic. That, coupled with growing risks such as multiple-resistant germs, increases the demand of changing from re-usable components to safe disposable products.

Single packed disposable valves:

- Reliability
- Traceability
- Cost-effective alternative
- Guidance supported

The iM Med sterile single-use valves are manufactured to the highest quality and at a cost-effective price, enabling departments to transition to this patient safety essential at a price point to support local cost improvement initiatives.

Minimise cross-contamination risks

- Sterile single-use pack
- Elimination of risk during reprocessing
- Reduce human error

Overcome limitations of tracking and traceability

- Single-use valves come with Lot No. for documentation

Offer consistent performance and improve unit efficiency

- Not subject to wear and tear
- Less expensive
- Reduce reprocessing time

Time to talk to iM Med.



AIR/WATER SUCTION BUTTONS



FOUR BUTTON VALVE SET

Contains suction, air/water, biopsy valves and auxiliary water connector fitting for Olympus endoscopes series 140, 160, 165, 180, 190, 240, 260, 290

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	200	830017	FAL19566



THREE BUTTON VALVE SET

Contains suction, air/water, biopsy valves fitting for Olympus endoscopes series 140, 160, 165, 180, 190, 240, 260, 290

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	200	830018	FAL19569



TWO BUTTON VALVES

Suction, air/water fitting for Olympus endoscopes series 140, 160, 165, 180, 190, 240, 260, 290

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	200	830059	FAL19569



AUXILIARY WATER CONNECTOR

For Olympus endoscopes

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	500	830013	FAL19562



PROCESS CLEANING BUTTON

CLEANING ADAPTOR

Fitting for Olympus endoscopes series 140, 160, 165, 180, 190, 240, 260, 290

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	200	830019	FAL19575



ENDOSCOPE BUTTON VALVE SET

4-piece set
Air/water, suction, biopsy, auxiliary water connector

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fujinon	200	830020	FAL19553



ENDOSCOPE BUTTON VALVE SET

3-piece set
Air/water, suction, biopsy valve

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fujinon	200	830062	Enquire with iM Med



STERILE BIOPSY VALVE

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fujinon / Olympus	400	830022	FAL19561



SINGLE USE POLYP TRAP - TWO CHAMBERS

Reducing the need for specimen transfer during a procedure, maximising efficiency.
Designed to handle specimen safely and prevent loss of specimen.

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	80	830036	FAL19544



BITE BLOCKS



BITE BLOCK - SENIOR
Open Size: 18x19mm

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	500	830033	FAL19579



BITE BLOCK - ADULT
Open Size: 25x20mm

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	500	830032	FAL19545



BITE BLOCK - PEDIATRIC
Open Size: 16x19mm

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	500	830034	FAL19554



STERILE BIOPSY VALVE

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fujinon / Olympus	400	830022	FAL19561



STERILE BIOPSY VALVE

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Pentax	400	830024	FAL19577



STERILE BIOPSY VALVE

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Bronchoscope	400	830029	FAL19550



ENDOSCOPY BIOPSY VALVES



BIOPSY VALVE WITH IRRIGATION PORT

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	400	830027	FAL19537



BIOPSY VALVE WITH IRRIGATION PORT

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fujinon / Olympus	400	830041	Enquire with iM Med



AIR/WATER TUBING WITH CO2 EXTENSION

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	40	830002	FAL19563



AIR/WATER TUBING WITH CO2 EXTENSION

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Pentax	40	830004	FAL19568



AIR/WATER TUBING WITH LIGHT SOURCE

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Pentax	40	830005	FAL19549



AIR/WATER TUBING WITH CO2 EXTENSION

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fujinon (500/600)	40	830007	FAL19574



AIR/WATER TUBING WITH CO2 EXTENSION

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fujinon (700)	40	830009	FAL19559



IRRIGATION TUBING



IRRIGATION TUBING

Olympus OFP2

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	40	830012	Enquire with iM Med



IRRIGATION TUBING

Olympus OFP, Endo-Stratus EGA,
Endo-Start Unit II

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus and others	40	830050	Enquire with iM Med



IRRIGATION TUBING

Olympus AFU-100 and Endogator EGP100
irrigation pumps

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus and others	40	830010	Enquire with iM Med



EXTENSION TUBING FOR IRRIGATION TUBING SET

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Most luer lock tubing sets	100	830051	Enquire with iM Med

1



AUXILIARY WATER CONNECTOR

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Olympus	500	830013	FAL19562



AUXILIARY WATER CONNECTOR

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Pentax	500	830014	FAL19567



AUXILIARY WATER CONNECTOR

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fujinon	500	830015	FAL19543



CONSUMABLES



RED TRANSPORT BAGS

Moisture retaining used endoscope containment bags

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fibro-cleaner	250	121002	FAL18543
ED flow	250	125006	FAL18544



ENDOSCOPE BEDSIDE CLEANING KIT

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
500ml	24	110105	FAL18395
300ml	24	110103	FAL18394

1



DIGNITY SHORTS

Made from recycled materials

SIZE	QUANTITY	iM PRODUCT CODE	NHS CODE
Extra large	100	810001	Enquire with iM Med
Large	100	810002	Enquire with iM Med

To place orders, please call or email:

T. +44 (0)1223 440 475

E. enquiries@im-med.com

If you require smaller order quantities of any of our consumable products, please call or email us.



CLEANING



Optimising the critical
cleaning stage of device
decontamination with an
effective portfolio of
cleaning solutions.



CONSUMABLES



SINGLE USE ENDOSCOPE CLEANING BRUSHES - DOUBLE ENDED WITH HANDLES

Head 1 diameter 10mm; length 30mm
Head 2 diameter 5mm; length 20mm
Total length 160mm; Handle length 70mm

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	1000	830039	FAL19552



SINGLE USE ENDOSCOPE CLEANING BRUSHES - DOUBLE ENDED WITH HANDLES

Head 1 diameter 5mm; length 25mm
Head 2 diameter 10mm; length 45mm
Total length 2300mm; Handle length 45mm

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	400	830044	Enquire with iM Med



SINGLE USE ENDOSCOPE CLEANING BRUSHES - DOUBLE ENDED

Head 1 diameter 5mm (or 4mm)
length 20mm
Total length 2300mm

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	400	830045	Enquire with iM Med



CLEANING/PROTECTION SPONGES

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	100	830037	FAL19542



NEUTRAL NON-ENZYMATIC DETERGENT FOR
MANUALLY CLEANING INSTRUMENTS

2x 5 Litre bottles

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	2	110205	FAL18396



UNIVERSAL FACE SHIELD SETS

Three visors per set

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	400	810011	Enquire with iM Med



FACE SHIELDS/WATER REPELLENT PROTECTIVE
VISORS

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	200	810005	Enquire with iM Med



GAUNTLET SLEEVE

Water repellent sleeve cover in packs of 100

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	2000	810004	Enquire with iM Med



WATER REPELLENT GOWNS

COLOUR	QUANTITY	iM PRODUCT CODE	NHS CODE
Blue	200	810008	Enquire with iM Med



ENDOSCOPE MANUAL CLEANING



Steelco

Miele Group
Member

STEELCO EPW 100S

The EPW 100S allows the complete assistance of manual cleaning, taking care of the most stressing repeated actions (syringe flushing and rinsing) and drastically reducing the risk of errors thanks also to the OCS system, the One-time Connection System shared between pre-washing, ashing/disinfection and drying/storage devices.

Manual Cleaning automation

Leak test, powerful and effective channel flushing & rinsing

Fast cycle start up

Endoscope and operator are selected by the "built in" barcode reader or RFID System

Process traceability

Detailed print report of all the cycle phases transferred to a USB stick or integrated into SteelcoData


Suitable for organic residual sampling

Immediately after the channels rinsing phase

2



: one-time connection system



SCOPETECH DUAL FLUSH

4-in-1 Flexible Endoscope Cleaning System

The ScopeTech Dual Flush dramatically improves the speed of endoscope reprocessing, while providing step-by-step cleaning process management, cycle validation and reports.

- Built-in automated scope inflation: for water immersion leak test
- Integrated detergent pump: doses the right amount and tracks detergent conditions
- Two flushing / rinsing ports: with single or dual mode operation
- Cleaning sequence control: promotes consistent processing while recording
- Scope inflation and cleaning process reports: Windows compatible software and Excel reports
- Sequence and express processing modes: simple and easy to use
- Compact, watertight cabinet mounts: allows for mounting on a wall or countertop
- Eliminates technician muscle fatigue: quiet and efficient pumping through endoscope channels



FLUSHTECH GI

Automated Endoscope Flushing - No Manual Syringing

The FlushTech GI is a powerful, easy to operate endoscope flushing aid that uses single and dual flushing modes to move chemistry and water through all the channels of an endoscope.

- Simultaneous port flushing
- Pulsating flush technology
- Flushing sensor with alarm
- Manual or auto flush settings
- High-capacity flushing
- Reporting capabilities
- No syringing required
- System decontamination cycle



Disinfectant for self disinfection of flushing systems available, enquire for further information



REPROCESSING



Guaranteeing effective
decontamination in
automated reprocessing
with a portfolio of UK
manufactured chemistries
and accessories.



CLEAN STAGE 3

BRITISH MANUFACTURED CHEMISTRY

Peracetic acid endoscopy disinfectants designed and manufactured within an accredited ISO 13485:2016 Quality Management System and CE marked.

With BSI as our notified body, we manufacture type tested disinfectants in our manufacturing facility in Yorkshire, England.

We supply decontamination facilities across the country with their AED approved chemistry. By choosing a chemistry manufactured in the UK, Trusts can add to their ongoing environmental targets and carbon footprint reduction by selecting UK manufactured rather than imported chemistry.



QUALITY ASSURED

A devotion to quality and patient safety forms the core of all we do and iM Med is proud to hold the ISO 13485: 2016 accreditation.

The ISO 13485 standard is an effective solution to meet the comprehensive requirements for a Quality Management System. Adopting ISO 13485 provides a practical foundation for manufacturers to address the Medical Device Directives, regulations and responsibilities, as well as demonstrating a commitment to the safety and quality of medical devices.



STEELCO APPROVED PERACETIC ACID

For customers selecting the Steelco endoscopy washer disinfector range, iM Med provides the Steelco manufacturer approved iM Peracept 15™ peracetic acid and iM EZcept™ enzymatic detergent. These UK manufactured chemistries are type tested and approved on the NHS Supply Chain Decontamination Accessories Framework and they are available for immediate delivery, removing the time risks that imported products can encounter.



iM Peracept15™

PERACETIC ACID (PAA) DISINFECTANT

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	110515	FYY011

iM EZCcept™

ENZYMATIC DETERGENT

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	110050	FYY066

ENDOSCOPY CHEMISTRIES

PERACETIC ACID (PAA) DISINFECTANT

For the disinfection of endoscopes in an endoscopy washer disinfectant compatible with a single bottle 2% PAA

 PeraceptFC™

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	110502	FAL18390

PERACETIC ACID (PAA) DISINFECTANT

For the disinfection of endoscopes in an endoscopy washer disinfectant (EWD) compatible with a single bottle 15% PAA. Steelco manufacturer approved for use in the Steelco EWD range

 Peracept15™

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	110515	FYY011

PERACETIC ACID (PAA) DISINFECTANT

For the disinfection of endoscopes in an endoscopy washer disinfectant compatible with a two part 5% PAA. Must be used in combination with iM Peracept DUO B™

 PeraceptDuoA™

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	110505	FYY014

ALKALINE CHEMISTRY

For the disinfection of endoscopes in an endoscopy washer disinfectant compatible with a two part 5% PAA. Must be used in combination with iM Peracept DUO A™

 PeraceptDuoB™

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	110500	FYY007

ENZYMATIC DETERGENT

For use in endoscopy washer disinfectors utilising iM Peracept 15™. Steelco manufacturer approved for use in the Steelco EWD range

 EZCept™

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	110050	FYY066



iM KartWash™

ACIDIC AUTOMATED TROLLEY WASHING DETERGENT

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
10 Litres	1	113210	Enquire with iM Med
20 Litres	1	113220	Enquire with iM Med



iM EnzWash™

pH BALANCED ENZYMATIC DETERGENT

For use in washer disinfectors and ultrasonic machines

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	113005	Enquire with iM Med
10 Litres	1	113010	Enquire with iM Med

iM AlkWash™

ALKALINE HOSPITAL INSTRUMENT CLEANER AND PRE-SOAK CONCENTRATE

For use in washer disinfectors, ultrasonic baths and flexible endoscope processing

VOLUME	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Litres	2	113105	Enquire with iM Med
10 Litres	1	113110	Enquire with iM Med

PAA AND H2O2 MONITORING

ARE YOU MONITORING EXPOSURE TO PAA IN YOUR ENDOSCOPY DEPARTMENT?

This is a JAG/IHEEM annual audit requirement.

DO YOU MONITOR H2O2, SUCH AS FROM LOW TEMPERATURE STERILISERS?

This is a part of COSHH risk management.



ChemDAQ STERI-TRAC 2+ AREA GAS MONITOR

Continuously monitor and alert you to hazardous vapor levels in the air: customisable alerts, plug and play set up, full data collection

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	Single or multiple unit options	ST2P-MON-W	Enquire with iM Med



ChemDAQ SAFECIDE 2.0 PORTABLE GAS MONITOR

Rugged handheld device - fully portable gas monitoring system: fast, accurate, real-time measurements of chemical vapor

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	Single or multiple unit options	Enquire with iM Med	Enquire with iM Med



ChemDAQ DATA ACQUISITION SYSTEM

PC-based computer software system that continuously gathers, analyses, reports, displays and stores

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	Covers multiple units	Enquire with iM Med	Enquire with iM Med



PRINTER ACCESSORIES (INK JET)

PRODUCT	QUANTITY	iM PRODUCT CODE	NHS CODE
Printer ribbon	1	800002	Enquire with iM Med
Printer roll for integral printer	20	800003	Enquire with iM Med



PRINTER ROLLS (THERMAL)

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
For ED flow	20	800010	Enquire with iM Med
For Fibro-cleaner	20	800001	Enquire with iM Med
For FD8E printer	20	800011	Enquire with iM Med



BARCODE P-TOUCH PRINTER FOR FD8E

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	1	900150	FYY063



CONSUMABLES



BLACK AND WHITE PRINTER TAPES
for P-Touch printer

SIZE	QUANTITY	iM PRODUCT CODE	NHS CODE
9mm	1	900151	FYY061
12mm	1	900152	FYY064
24mm	1	900153	FYY060



6 CHANNEL SURROGATE DEVICE

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fibro-cleaner	1	800006	Enquire with iM Med
Ed flow	1	800007	Enquire with iM Med



TVC KIT (TOTAL VIABLE COUNT)


Used for taking TVC water samples.
Kit is fully autoclavable.

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	6	900154	Enquire with iM Med



EDS8 IRRIGATOR DISINFECTION KIT

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
5 Channel	1	900087	Enquire with iM Med
6 Channel	1	900086	Enquire with iM Med



FIBRO-CLEANER LINK BASKETS

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
With centre basket	1	121020	Enquire with iM Med
Without centre basket	1	125022	Enquire with iM Med



PLASTIC FEET FOR EWD BASKETS

Bolt-through

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	1	121021	Enquire with iM Med



ENT ENDOSCOPE BASKET PARTS

PRODUCT	QUANTITY	iM PRODUCT CODE	NHS CODE
Base for Fibro-cleaner	1	121023	Enquire with iM Med
Base for ED flow	1	125024	Enquire with iM Med
Support rail for Fibro-cleaner	1	121025	Enquire with iM Med
Support rail for ED flow	1	125026	Enquire with iM Med
Silicon matting (autoclavable)	1	121027	Enquire with iM Med



ENT LEAK TEST MANIFOLD

L/T connection and tubing not supplied

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	1	121028	Enquire with iM Med



OPERATOR RFID CARD

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
N/A	10	500051	Enquire with iM Med

CONNECTORS



The iM Med connector range is manufactured in the UK for the broad range of endoscopes used clinically.

Connectors and separators are produced by precision computer controlled machining methods to ensure they consistently meet the necessary quality controlled manufacturing standards.



ALL CONNECTORS AND SEPARATORS CAN BE SOLD INDIVIDUALLY.

LEAK TEST CONNECTORS



ED FLOW LEAK TEST CONNECTOR
(BEIGE FEMALE)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
All scopes	N/A	900009	Enquire with iM Med



FIBRO-CLEANER LEAK TEST CONNECTOR
(BEIGE MALE)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
All scopes	900119	900005	Enquire with iM Med



LEAK TEST CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon	900108	900030	Enquire with iM Med



ETO LEAK TEST CAP

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon	N/A	900164	Enquire with iM Med



LEAK TEST CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Storz	900109	900034	Enquire with iM Med



LEAK TEST PLUG FOR FIBRO CLEANER / STEELCO

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	900119	900017	Enquire with iM Med



LEAK TEST CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	900106	900052	Enquire with iM Med



ETO LEAK TEST CAP

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus / Storz	N/A	900162	Enquire with iM Med



LEAK TEST CONNECTORS



LEAK TEST CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	900107	900059	Enquire with iM Med



ETO LEAK TEST CAP

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	N/A	900163	Enquire with iM Med



ETO LEAK TEST CAP

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Storz	N/A	900165	Enquire with iM Med

3



CHANNEL SEPARATOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon 580	N/A	900168	Enquire with iM Med



CHANNEL SEPARATOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus 1	900092	900035	Enquire with iM Med



CHANNEL SEPARATOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus 2	900093	900036	Enquire with iM Med



VALVE HOUSE CONNECTORS + CHANNEL SEPARATORS



CHANNEL SEPARATOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax 70/80	900097	900040	Enquire with iM Med



CHANNEL SEPARATOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax 90	900115	900041	Enquire with iM Med



CHANNEL SEPARATOR ENDOCHOICE GI
SCOPES

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	900120	900043	Enquire with iM Med



CHANNEL CONNECTOR DRYER

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon	N/A	900137	Enquire with iM Med



BRONCHOSCOPE CONNECTOR (NEW)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon	N/A	900033	Enquire with iM Med

3



CHANNEL SEPARATOR V2

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon G5	N/A	900042	Enquire with iM Med



CHANNEL SEPARATOR (FIBRO-CLEANER)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	900101	900023	Enquire with iM Med



CHANNEL SEPARATOR (FD8)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	900099	900025	Enquire with iM Med



CHANNEL SHUTTER

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	900105	900044	Enquire with iM Med



AIR, WATER, SUCTION, BIOPSY, AUXILIARY + ELEVATOR CHANNEL CONNECTORS



WATER CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon	N/A	900127	Enquire with iM Med



FIBRO-CLEANER AND ED FLOW AUX/BRIDGE
CHANNEL

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	N/A	900002	Enquire with iM Med



WATER CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon	N/A	900031	Enquire with iM Med



AIR WATER CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	N/A	900138	Enquire with iM Med



SUCTION CONNECTOR (THREE PIN)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Storz	N/A	900126	Enquire with iM Med

3



SUCTION CONNECTOR (TWO PIN)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon	900104	900045	Enquire with iM Med



MALE LUER LOCK CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	N/A	900067	Enquire with iM Med



AUXILIARY CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	N/A	900053	Enquire with iM Med



BRIDGE CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	N/A	900054	Enquire with iM Med



BRON BIOPSY CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	900090	900004	Enquire with iM Med



AIR, WATER, SUCTION, BIOPSY, AUXILIARY + ELEVATOR CHANNEL CONNECTORS



BRONCHO SUCTION CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	900112	900003	Enquire with iM Med



CHANNEL CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	900103	900075	Enquire with iM Med



WATER CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Olympus	N/A	900049	Enquire with iM Med



AIR CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon / Olympus	900121	900047	Enquire with iM Med



BIOPSY CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Fujinon / Olympus	900091	900001	Enquire with iM Med

3




BUTTON SHUTTER

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	900116	900058	Enquire with iM Med



IRRIGATION CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	N/A	900060	Enquire with iM Med



AIR, WATER, SUCTION, BIOPSY, AUXILIARY + ELEVATOR CHANNEL CONNECTORS



WATER JET CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	N/A	900056	Enquire with iM Med



SUCTION CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	N/A	900048	Enquire with iM Med



CHANNEL CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	900102	900055	Enquire with iM Med



SHUTTER OF-B153

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	900132	900132	Enquire with iM Med



SHUTTER OF-B212

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
Pentax	N/A	900133	Enquire with iM Med

3



CLEAR
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
Small hose	No	900046	Enquire with iM Med
Steelco hose	No	900175	Enquire with iM Med
600mm tubing 4x2 mm	No	900010	Enquire with iM Med
100mm tubing (small ID)	No	900084	Enquire with iM Med
600mm tubing	Yes	900069	Enquire with iM Med



BEIGE
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
600mm hose 4x2mm	No	900011	Enquire with iM Med



GREY
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
Steelco Hose	No	900179	Enquire with iM Med



RED
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
Steelco Hose	No	900180	Enquire with iM Med



BLUE
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
Steelco Hose	No	900178	Enquire with iM Med



SCOPE CHANNEL CONNECTORS



BLACK
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
Steelco hose	No	900182	Enquire with iM Med
600mm tubing 4x2 mm	No	900016	Enquire with iM Med
600mm tubing	Yes	900074	Enquire with iM Med



BROWN
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
Steelco hose	No	900176	Enquire with iM Med
600mm tubing 4x2 mm	No	900013	Enquire with iM Med
600mm tubing	Yes	900071	Enquire with iM Med



GREEN
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
Steelco hose	No	900177	Enquire with iM Med
600mm tubing 4x2 mm	No	900015	Enquire with iM Med
600mm tubing	Yes	900073	Enquire with iM Med



YELLOW
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
Steelco hose	No	900181	Enquire with iM Med
600mm tubing 4x2 mm	No	900012	Enquire with iM Med
600mm tubing	Yes	900070	Enquire with iM Med




ORANGE
HOSES AND TUBING

SIZE	M/C CONNECTOR	iM PRODUCT CODE	NHS CODE
600mm hose 4x2mm	No	900014	Enquire with iM Med
600mm tubing	Yes	900072	Enquire with iM Med



One-Time Connection System
for all processing and storage



Miele Group Member

PRODUCT	COMPATIBILITY	iM PRODUCT CODE	NHS CODE
OCS Hub	Steelco	900170	Enquire with iM Med
OCS Socket	Steelco	900171	Enquire with iM Med
OCS Black Flow Restrictor	Steelco	900183	Enquire with iM Med
OCS White Flow Reducer	Steelco	900128	Enquire with iM Med

SCOPE CHANNEL CONNECTORS



ED FLOW CHANNEL TUBE CPC CONNECTOR
(GREY 90°)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	N/A	900008	Enquire with iM Med



FIBRO-CLEANER OR FD8 CHANNEL TUBE
CPC CONNECTOR (GREY STRAIGHT)

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	N/A	900006	Enquire with iM Med

3



Y-PIECE POLYPROPYLENE

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	N/A	900018	Enquire with iM Med



CPC INLINE FEMALE CONNECTOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
EDS8 Cleaning Kit	N/A	900185	Enquire with iM Med



5 CHANNEL IRRIGATOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
EDS8	N/A	900085	Enquire with iM Med



6 CHANNEL IRRIGATOR

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
EDS8	N/A	900086	Enquire with iM Med



ED FLOW M/C AUX/BRIDGE CHANNEL

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
N/A	N/A	900158	Enquire with iM Med



IRRIGATOR DISINFECTION KIT

COMPATIBILITY	O-RING KIT	iM PRODUCT CODE	NHS CODE
EDS8	N/A	900087	Enquire with iM Med



SCOPE CHANNEL CONNECTORS



FEMALE PLASTIC LUER LOCK

COMPATIBILITY	O-RING KIT	PRODUCT CODE	NHS CODE
N/A	N/A	900122	Enquire with iM Med



MALE LUER TO HOSE BARB ADAPTER 1/4"

COMPATIBILITY	O-RING KIT	PRODUCT CODE	NHS CODE
N/A	N/A	900131	Enquire with iM Med



MALE PLASTIC LUER LOCK

COMPATIBILITY	O-RING KIT	PRODUCT CODE	NHS CODE
N/A	N/A	900123	Enquire with iM Med



POLYPROPYLENE T PIECE

COMPATIBILITY	O-RING KIT	PRODUCT CODE	NHS CODE
N/A	N/A	900136	Enquire with iM Med



NEW ELEVATOR MACHINE CONNECTOR SS

COMPATIBILITY	O-RING KIT	PRODUCT CODE	NHS CODE
N/A	N/A	900088	Enquire with iM Med

3



COMPLIANT, COMPATIBLE, TYPE TESTED
DECONTAMINATION CHEMISTRIES

iM-Med.com/consumables

MANUFACTURED IN THE UK



DRYING, STORAGE, TRANSPORTATION



Ensuring the drying,
storage and
transportation of
decontaminated medical
devices is maintained for
optimised patient safety.



CLEAN STAGE 4

FLEXIBLE ENDOSCOPE TROLLEYS

Safe and secure transportation.



Endoscope transportation trolleys manufactured in robust high quality stainless steel and to the highest standards, to ensure equipment and staff safety.

- Lockable wheels
- Lockable doors
- Ergonomic height
- Additional space
- Tray protection
- Improved handling
- Easy to manoeuvre
- Easily cleanable with a wide range of standard products

DESCRIPTION	SUITABLE FOR THERMAL DISINFECTION IN A TROLLEY WASHER	iM PRODUCT CODE
304 stainless steel 5 level trolley	NO	126040
316L stainless steel 5 level trolley	YES	126041
316L stainless steel 5 level trolley with two handles	YES	126042

4



PROLONG THE ASEPTIC STORAGE AND SAFE TRANSPORTATION OF ENDOSCOPES FOR UP TO 30 DAYS.

ENDOSCOPE PACKING, STORAGE AND TRANSPORTATION SYSTEM

A solution for the aseptic storage and safe transportation of endoscopes. This unique process safely contains the endoscopes for up to 30 days.

Following decontamination and compliant drying in an endoscope drying cabinet, thereby ensuring all residual bacteria are fully desiccated, an endoscope is packed using the space-saving, free-standing iM EndoPax™ system.

The system is based on a preservation process widely used in the food industry to reliably eliminate the risk of microbial growth.

With a rapid, user-friendly process, the system delivers a process ensuring the internal and external lumens are fully conditioned using an inert gas mix, to inhibit aerobic and anaerobic bacterial growth.

Without the addition of chemistry to the terminal packing process, potential risks are reduced to the endoscope, patient and user.

The medical device bag is heat sealed ensuring a reliable seal integrity until the user requires the endoscope.

The iM EndoPax™ system has an independent monitor to provide a fully validated and traceable process.

iM EndoPax™ enables endoscopes to be prepared and conveniently ready for routine lists, emergency situations, out of hours and off-site endoscopy work. Endoscopes can be transported in a controlled and tamper evident environment using transportation cases or lockable trolleys.

iM-Med.com/endopax



iM PRODUCT CODE	NHS CODE
600020	FYY004

CONSUMABLES



VACUUM BAGS

For use with iM EndoPax™ packing system

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
iM EndoPax™	80	123005	FYY008



WHITE SELF-ADHESIVE LABELS

For iM EndoVac™ and iM EndoPax™ IMS printer

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
iM EndoVac™ / iM EndoPax™	500	800005	Enquire with iM Med



CLEAR, GAMMA-IRRADIATED TRANSPORT BAGS

COMPATIBILITY	QUANTITY	iM PRODUCT CODE	NHS CODE
Fibro-cleaner	100	121001	Enquire with iM Med
ED flow	80	125005	Enquire with iM Med



iM ENDOVAC VACUUM BAGS FOR VAC-TRAY

SIZE	QUANTITY	iM PRODUCT CODE	NHS CODE
Extra large 600x800mm	100	123004	FAL18546
Large 600x750mm	100	123003	FAL18545
Medium 600x400mm	100	123002	FAL18393

4



TRANSPORTATION TRAYS

Compatible with a range of transport

SIZE	QUANTITY	iM PRODUCT CODE	NHS CODE
Extra-large 650x470mm	1	123014	Enquire with iM Med
Large 530x370mm	1	123013	Enquire with iM Med
Small 520x420mm	1	123011	Enquire with iM Med



CONTOURED TRANSPORTATION TRAY

SIZE	QUANTITY	iM PRODUCT CODE	NHS CODE
520x420mm	1	123015	Enquire with iM Med



PLASTIC FLAT BASE ENDOSCOPE TRAY (GREY)

SIZE	QUANTITY	iM PRODUCT CODE	NHS CODE
520X420mm	1	123021	Enquire with iM Med





Choice of Steelco
CSSD and endoscopy
cutting-edge
decontamination
solutions.



EQUIPMENT


GIVING YOU REAL CHOICE.

DISCOVER THE RANGE:

iM-Med.com/steelco



one connection.



STEELCO OCS: ONE-TIME CONNECTION SYSTEM

One-time endoscope handling between prewashing, washing/disinfection and drying/storage devices.

Reducing contamination risk

Error-proof connection

Three-times faster connection

DECONTAMINATION EQUIPMENT THAT FITS YOUR NEEDS.

DON'T ADAPT YOUR NEEDS TO FIT THE EQUIPMENT.



CHOICE OF

ENDOSCOPY DECONTAMINATION EQUIPMENT

One scope per chamber EWD, single or double options

One, two or three scopes per chamber EWD

Vertical and horizontal drying storage cabinets

CHOICE OF

SURGICAL INSTRUMENT DECONTAMINATION EQUIPMENT

Full automation

Steam and/or electric functionality

Single or double door options

Large choice of capacity

Fast cycle times



Part of the Miele Group, Steelco S.p.A. is leading manufacturer of decontamination equipment and systems for washing, disinfection and sterilisation of medical instruments in hospitals and other healthcare facilities.

Maximising safety
Maximising capacity
Minimising footprint

STEELCO ENDOSCOPY + CSSD DECONTAMINATION EQUIPMENT

iM Med partnered with Steelco to supply the full Steelco Endoscopy and CSSD range of equipment.

iM Med are the nominated provider of the Steelco range on the NHS Supply Chain Decontamination Equipment Capital Framework, with a team of decontamination specialists to discuss your requirements and design a range of solutions to meet your needs.


DISCOVER THE RANGE:

iM-Med.com/steelco





Delivering
decontamination training
and compliance for
continuous patient safety.



ACADEMY



TRAINING

A range of compliance training courses designed to match the exact requirements of your department and staff, delivered on-site or at an off-site location.

CPD certified training.

CONFERENCES

Aimed at teams involved in decontamination, from front-line decontamination technicians up to senior management supporting Trusts, to guarantee a team approach to patient safety.



PASSIONATE ABOUT PATIENT SAFETY
EXCELLENCE

Find out more at
iM-Med.com/compliance



“

The effective training of my staff is so crucial for the ongoing safe standards of my department.

The iM Academy has been fantastic. They developed a training plan with me in my role as decontamination lead and then delivered the training in an expert and engaging manner with my Team.

They stayed with us until all of the Team felt confident and iM are now our go to for any training and support we need. They understand what's needed and have provided peace of mind for the ongoing competence of my staff.

”



Specialist CSSD and
endoscopy technical
services: guaranteeing
compliance, optimising
performance.



TECHNICAL



iM MED TRULY UNDERSTANDS
WHAT OUR VALUED CUSTOMERS
REQUIRE TO DELIVER SAFE PATIENT
CARE.

iM-Med.com/technical

SPECIALIST TECHNICAL
DECONTAMINATION ENGINEERS

Dedicated engineers

Comprehensive contracts

Rapid response times

Guaranteed equipment uptime





BESPOKE REPORTING

Quarterly reporting

Customer portal for reports
and iM Academy resources

Microbiological: testing,
analysis and technical support

PLANNED SERVICE MAINTENANCE

Machine health checks

Prolong equipment life

Customer service excellence



ORDERS AND ENQUIRIES

T. +44 (0)1223 440 475

E. enquiries@im-med.com

WEBSITE

iM-Med.com

HEAD OFFICE

iM Med Ltd

1 Pembroke Avenue

Waterbeach

Cambridge CB25 9QP





DECONTAMINATION CATALOGUE

