



### EOS-H-01

- Works with the Full range of Innovex Single Use Endoscopes
- High capacity processor: Advanced image algorithm to improve image quality
- USB connection for external big capacity storage device
- Internal battery for time clock
- USB connection for external keyboard to simplify input operation
- Signal output: 2 DVI + 1 SDI + 1 VGA
- Accurate time counting
- Memorise last time settings before new settings



### EA100

- Portable DVP to allow for easy transportation
- TouchScreen 13.3"
- High capacity processor: advanced image algorithm to improve image quality
- USB connection for external storage device (flash drive)
- Signal output: DVI and SDI
- Battery 7800mAh
- IV pole and desk stand mounting clamp
- Works with the Full range of Innovex Single Use Endoscopes

**Shanghai AnQing  
Medical Instrument Co., Ltd.**  
(Subsidiary of INNOVEX)

3&4 Floor, No.2 Building, 366 Huiqing Rd, East  
Zhangjiang High-Tech Park, 201201 Shanghai, China  
Tel: +86 21 61117375  
Fax: +86 21 61117374  
Email: [global@innovexmed.com.cn](mailto:global@innovexmed.com.cn)




## SINGLE USE FLEXIBLE BRONCHOSCOPE

SINGLE USE STILL NEEDS PERFORMANCE

# WHEN VISION MATTERS

INNOVEX Single-Use Flexible Bronchoscope provides every patient a clean scope with high quality view – every time. Eliminate the cost and risk associated with reprocessing and cross-contamination.



1 Lever lock

3 Photo & video button

5 Aspiration connector

7 Dual deflection:  
210°up / 180°down

9 Excellent maneuverability  
Working length 620mm

2 Deflection Lever

4 Aspiration button

6 Working channel luer port

8 Field of view: 110°±10%


10 Excellent maneuverability  
Working length 670mm

● 640K pixels


**BS41H-12EU/US**  
4.1mm/1.2mm



**BS50H-20EU/US**  
5.0mm/2.0mm




**BS59H-28EU/US**  
5.9mm/2.8mm




● 160K pixels

**BS27U-12EU/US**  
2.8mm/1.2mm



**BS38U-20EU/US**  
3.8mm/2.0mm



DMR-RC2-005 A/1