


SIMAI GENIUS MORCELLATOR SYSTEM



Revolutionising Urology Surgery with Precision & Efficiency

ADVANCED CUTTING TECHNOLOGY

- Advanced cutting design ensures fast rotary morcellation
- Consistent Suction for superior tissue retention
- Adjustable blade speed (up to 20g/min)
- High quality and durable cutting blades
- Reduced operating time with enhanced safety and efficiency.

ERGONOMIC, LIGHT & COMPACT DESIGN

- Reduces fatigue during long procedures
- Easy-to-detach motorised handle for sterilisation
- Compact: Fast Setup and easy moveability

LARGE CAPACITY COLLECTION

- 8 litre capacity for greater efficiency during procedures.



SIMAI

GENIUS

MORCELLATOR SYSTEM

CLINICAL BENEFITS



Shorter Operating Times

Improved efficiency reduces the duration of surgeries, decreasing anaesthesia time and patient recovery period



Precision Surgery

High speed and effective morcellation ensures a clear surgical field and visibility, reducing potential trauma to surrounding tissue

WHY SIMAI GENIUS MORCELLATOR?

- ✓ The Simai Genius Morcellator is a game changer for BPH enucleation procedures, offering unmatched speed and efficiency
- ✓ Provides fast and precise morcellation, processing 13-20 grams per minute
- ✓ Features a lightweight handle design and single foot pedal for easy setup
- ✓ Compact dedicated trolley, making it ideal for all operating theatres.

KEY SPECIFICATIONS

Morcellator Blade	Working length 360mm / Diameter 4.8mm
Fast Suction Pressure	0-700 mmGh (0-90 Kpa)
Foot Switch	Single foot pedal offering 2 stage operation – suction, suction & morcellation
Collection Canisters	2 x 4 Lt collection canisters
Input Power	300VA
Safety Class	I
Sterilisation	Motor Handle is easy to detach and can be sterilised by auto-clave or low temperature plasma
Morcellator Blades	Single use and Reusable options