









## Crown & Bridge Material Selector - Which is Best for Your Needs

<b>BENEFITS &amp; FEATURES</b>	 <b>PFM</b>	 <b>Full Cast semi-precious</b>	 <b>Full Cast high-noble</b>	 <b>Zirconia Solid full contour</b>	 <b>Multi Layer Zirconia</b>	 <b>HT Zirconia</b>	 <b>E-Max</b>	 <b>E-Max Empress</b>
Strength	Strong	Very high	Excellent	1200 MPa+, superior	Strong	600-800 MPa	400 MPa	160 MPa
Durability	Durable	Very high	Superior	Highly Durable	Durable	Durable	Durable	Durable
Longevity	20-30 years*	30+ years*	20-30 years*	30+ years*	20-30 years*	20 years*	20 years*	20 years*
Tissue Response	Healthy, good	Healthy	Excellent	Healthy	Healthy	Healthy	Excellent	Healthy
Free From Metal	✗	✗	✗	✓	✓	✓	✓	✓
Translucence	Moderate	None	None	Opaque	Life-like	Natural	Very Natural	Very Natural
Versatility	Anywhere	Posterior	Posterior	Anywhere	Ant/Bridges	Ant/Premolars	Corwns/Veneers	Veneers Anterior Crown
Integrity Level	Good	High	High	High	Good	Good	High	Good
Shade Matching	Fair	Poor	Poor	Good Opaque	Very good	Exceptional	Exceptional	Exceptional
Plaque Resistance	High	High	High	High	High	High	High	High
Biocompatibility	Good	Exceptional	Exceptional	Exceptional	Exceptional	Exceptional	Exceptional	Exceptional
Oplaescence	Low	None	None	None	Moderate	High	Excellent	Excellent
Esthetic Level	Moderate	Low	Low	Low	High	Very high	Excellent	Excellent
Additional Info	Reliable, esthetic, affordable	Great balance of cost & quality	Posterior Use - Not for Bruxers	Strongest ceramic, ideal for bruxers & cost effective alt. to full cast crown	Esthetic bridge alternative, cosmetic	Esthetics w/enhanced translucency	Best for veneers, anteriors	Cosmetic

Call or chat with our technical support team: **(877) 785-9710**

\* Lifespan may vary based on individual patient care and maintenance.

Learn more at [summit-horizon.com](https://summit-horizon.com)

