

## Defined Perimeter Mattress Cover

HCPCS:



- The Defined Perimeter Mattress Cover creates a raised rail, defined perimeter for enhanced fall prevention, without using patient restraints
- Open in mid-section for safe and easy patient egress and ingress
- Built-in medium density foam bolsters safely and comfortably define edges of bed for patient and prevent falling
- Easy-to-install and secure
- Breathable Nylon makes this a truly universal cover, compatible with both foam and powered air surfaces

PRODUCTS	Description
14333	UniCover,DefinedPerimeter, 1/EA
14333-42	UniCover, Defined Perimeter, 42", 1/ea
SPECIFICATIONS	Description
Cover Material	Nylon
Fits Mattresses	34-43" (W) X 80" (L) X 6"-10" (H)
Foam Density	Medium
Foam Sections	5.5" (W) x 29.5" (L) x 4" (H)
Shipping Weight	5 lbs
Warranty	Limited 1 Year