

Deluxe Aluminum Bath Chair

HCPCS: E0240



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Easy, tool-free assembly of back seat and legs (Figure A - B)
- Aluminum frame is lightweight, durable and corrosion proof
- Legs are height adjustable in 1" increments, and angled with suction-style tips to provide additional stability
- Made using recycled materials

PRODUCTS

Description

RTL12202KDR	Bath Bench, w/Back, KD,RTL,1 EA
-------------	---------------------------------

SPECIFICATIONS

Description

Carton Shipping Weight	8.5 lbs
Outside Legs	19.5" (W) x 19.5" (D)
Seat Dimensions	19.25" (W) x 11.5" (D)
Seat Height	14 21
Warranty	Limited Lifetime
Weight	7 lbs
Weight Capacity	300 lbs