

Lightweight Bariatric Homecare Bed

HCPCS: 15300LW: E0301 | 15302LW and 15303LW: E0302



- Lightweight frame design offers providers easier and more manageable bed deliveries and assemblies
- Allows for additional space for patient comfort and recovery as well as caregiver efficiency
- Product Safety & Performance Standards: Meets IEC 60601-2-52; 60601-1; 60601-1-2; 60601-1-11
- Safe working load capacity of 600 lbs
- Full electric bariatric homecare bed with dual motor assembly for head and frame adjustability
- Ultra lightweight frame is easily set up by 1 individual; no links to connect and no tools required
- Square tube steel frame slat deck provides durable strength, no-sag mattress support, and is easy to clean
- Scratch resistant wood grain panel bed ends include hi lo motors for frame height adjustability
- UL approved motor is self contained for reduced weight and noise and can be easily installed and removed with patient in bed
- Reduce service calls with large glow in the dark hand pendant that is easy-to-use for adjustments; LED indicators available for power and lockout included
- Battery Back up accessory available (#15302LWB) and sold separately
- Half Rails accessory available (15021BV-LW) and sold separately
- Rail accessories meet FDA and ASTM entrapment guidelines and standards
-
-

PRODUCTS	Description
15300LW	Lightweight Full Elec 42" Bariatric Bed
15302LW	Lightweight Full Elec 48" Bariatric Bed
15303LW	Lightweight Full Elec 54" Bariatric Bed

SPECIFICATIONS	Description
Battery	Backup available but not included
Bed Height	15300LW: 34" - 42" with casters 15302LW and 15303LW: 35" - 44" with casters
Carton Shipping Weight	15300LW: 255 lbs (2 boxes) 15302LW: 294 lbs (3 boxes) 15303LW: 312 lbs (3 boxes)
Casters	5" casters (2 locking and 2 unlocking)
Color	Brown
Ground Clearance	15300LW: 18 - 26" 15302LW and 15303LW: 16 - 24"
Motor/Drive System	15300LW: Dual Full Electric w/Hi Lo 15302LW and 15303LW: Linear Full Electric w/Hi Lo
Sleep Surface	15300LW: 42"W x 80" L 15302LW: 48"W x 80" L 15303LW: 54"W x 80" L
Warranty	5 years on Motor and Hand Control Limited Lifetime on Welds and Frame 1 year on all other parts and components
Weight Capacity	15300LW: Safe working load 600lbs (to include occupant and any accessories) 15302LW: Safe working load 750lbs (to include occupant and any accessories) 15303LW: Safe working load 1000lbs (to include occupant and any accessories)

SPECIFICATIONS**Description**

Material Composition	Bed Frame: Steel Bed Ends: Steel base & Wood grain panels Motors: electrical components
----------------------	---

ACCESSORIES**Description**

15302LWB	LW Bari Bed Accessory battery pk,48/54, BX
----------	--

15021BV-LW	Half Side Rail,42in/48/54 LW Bari bed,1pr, BX
------------	---