




Molift RgoSling 4-point

Molift RgoSlings offer a wide range of models and sizes ensuring we have a solution, whatever the users requirements. Molift RgoSlings are designed for a 4-point sling bar.





Molift RgoSling Repositioning sheet

An in-bed lift and positioning sheet

- Turning and transferring patients higher in the bed.
- Lifting from bed to another supine surface, stretcher etc.
- Stays in bed - for comfortable and dignified care
- Facilitates care situations - saves time for caregivers.
- Colour coded straps facilitate communication when attaching the straps equally on both sides.

Size	Quantity	Repositioning sheet
L 1340 x W 1300 mm	1 pc	1721570
L 1340 x W 1300 mm	10 pcs	1721570-C

Maximum user mass: 500 kg