

Molift Stretcher

A scoop stretcher, ideal for horizontal hoisting situations

Easy to use

Molift Stretcher is a scoop stretcher for hospitals and institutions. It is designed for horizontal transfers and when sling hoisting is considered inappropriate due to medical conditions. Molift Stretcher is an optimal solution for transfer of clients who are comatose or sedated, as well as clients with suspected spinal injuries and/or multiple fractures.

Ideal for emergency care

In combination with Molift standard 4-point sling bar, it is ideal for emergency care needs, when time is critical and the patient has limited mobility. Molift Stretcher can be divided in two parts in length and is therefore easy to place under the client.

Molift Stretcher has five handles along each side, which provides a manual hoisting alternative from the floor and allows horizontal and lateral transfer between beds and surfaces.

Features

- A unique patented scoop design, for safe and efficient moving.
- Low total weight and easy to carry.
- A safe locking system - easily opened when unloaded - completely locked when loaded.
- Many possible uses and compatible with X-ray, CT and MRI. Designed for minimal shadow in X-ray.



Description	Item no.	Dimensions	
Molift Stretcher MRI	2150108	Height Length outer Width outer	126 mm 1980 mm 630 mm
Material Plastic composite body Stainless steel locks		Compatibility Molift Stretcher MRI: Compatibility with MRI. Approved for X-ray and CT-scan	<div> </div> <div> </div> <div> </div>
Weight Total weight	14.5 kg		<p>* The SWL always dependent on the hoists SWL when less than the stretcher.</p> <p>For more information visit etac.com</p>