



# PRODUCT OVERVIEW



## Manufacturing Certifications with full certification.

- FDA Registered
- FDA 510(k) Premarket Notification
- ISO13485 Quality System – Certificate
- CE Mark
- Certified ISO Class 8 Clean Rooms
- CPP (Critical Product Production) Room
- ISO 10993 Biocompatible
- ISO 11135 Sterilization Validation
- 16CFR1610:2007 – Compliant Flammability Rated



[www.onyxmedical.com](http://www.onyxmedical.com)



## Product Groups:

1. Microscope Drapes
2. Dome Bags and Band Bags
3. C-Arm Drapes
4. Camera Drapes
5. Probe Covers
6. Miscellaneous
  - Back Table Covers
  - Surgical Gowns
7. Face Shields

## 1. Microscope Drapes



## Major Surgical Microscope Brands



## 1. Microscope Drapes



[www.onyxmedical.com](http://www.onyxmedical.com)



1. Wide range of Microscope drapes to fit specific microscope needs.
2. The Microscope drapes are made from high quality polyethylene film.
3. Easy to drape with a special anti-blocking formula to prevent the PE film from sticking.
4. Straps to secure the drape allow a quick and easy draping.
5. All the drapes come with an angled lens to prevent glare during the procedure.
6. The lens covers used are made of high grade material and offer a high optical quality, no distortion and clear visibility.



## 1. Microscope Drapes

# Microscope Draping




[www.onyxmedical.com](http://www.onyxmedical.com)

### Questions

- Lens holder diameter?
- Length of the drape?
- Width of the drape?
- Number of Oculars?
- When available code used for cross reference?

Lens holder diameter:

Inner diameter of the lens holder =  
outer diameter of the Microscope  
lens housing




## 1. Microscope Drapes



[www.onyxmedical.com](http://www.onyxmedical.com)

# Unique Features









## 1. Microscope Drapes



[www.onyxmedical.com](http://www.onyxmedical.com)

## Instruction for use:






## 2. Band Bags and Dome Bags



[www.onyxmedical.com](http://www.onyxmedical.com)




A wide range of Dome and Band bags with latex free elastic, packed with a wrap so picking and using a bag is a simple job. All our bags are made of high grade special anti-blocking PE Film



Standard positioning  
tapes on every Band Bag.

High Quality sewn-in elastic  
in Dome Bags.





Latex free elastic (blue)  
bands in all band bags

## 2. Band Bags and Dome Bags



[www.onyxmedical.com](http://www.onyxmedical.com)



Dome Bags are round, Band Bags are Square with an elastic opening.

These Bags are the most versatile equipment drape. Can be used to cover virtually everything in the OR/Cath Lab.

Computers, monitors, keyboards, chairs, foot switches, etc.

Also used to cover imaging equipment such as C-Arms, Fluoroscopes.

Offers a wide range in sizes.

C-Arm Drapes and C-Arm Kits

### 3. C-Arm Drapes and C-Arm kits




[www.onyxmedical.com](http://www.onyxmedical.com)

## Major C-arm Brands

**PHILIPS**

**SIEMENS**



### 3. C-Arm Drapes and C-Arm kits





[www.onyxmedical.com](http://www.onyxmedical.com)





#### C-Arm Drape Types:

- One piece straight type
- One piece C-Shape



#### C-Arm Kit Drape types:

- Clip drape with upper and lower part (Band or Dome Bag)
- 2-part C-Shape
- Foot switch cover (optional)




### 3. C-Arm Drapes and C-Arm kits



[www.onyxmedical.com](http://www.onyxmedical.com)



The C-Arm kit uses an organizer pouch in packaging.  
End-user can easily locate the right part.





### 3. C-Arm Drapes and C-Arm kits



[www.onyxmedical.com](http://www.onyxmedical.com)

C-Pouch to preserve sterility when arm moves from non-sterile to sterile field.



C-arm rotation





Up/Down travel



Swiveling



Axial rotation







### 3. C-Arm Drapes and C-Arm kits



[www.onyxmedical.com](http://www.onyxmedical.com)

#### O-Shape drape




#### 4. Camera Drapes

A surgical video camera is used during endoscopic applications for example in laparoscopic or arthroscopic procedures.



[www.onyxmedical.com](http://www.onyxmedical.com)

Camera with cable



Scope



processor




Only part draped is the camera and cable.  
The Scope is cleaned and disinfected.

In a Laparoscopic procedure use of Scope Warmer is recommended


# Key in choice of camera drapes:

- ☐ Width of the tube film
- ☐ Length of the tube film
- ☐ Insert assistance, optional with Cardboard, or Ring
- ☐ Use of elastic bands and tape
- ☐ Tip of the camera Drape:
  - Open
  - Closed (tear-off)
  - Elastic White or Blue Kraton
  - Coupler system (allows change of Scope without re-draping)

### Camera Drape with Cardboard insert assistance




Camera Drape with ring insert assistance design insert assistance.




## 4. Camera Drapes




White Elastic Tip




Blue Elastic Tip



Perforated Tip




Open Tip




## 4. Camera Drapes

Closed Camera Drape system.  
Change scope without re-draping





## 5. Probe Covers



You can find a huge variety in Ultrasound probe covers.  
This results in a wide offering in sizes.  
On top of that Material, tip construction  
and use of Gel adds to the number of codes

## 5. Probe Covers

### **Material Options offers options in material used:**

- Polyethylene (PE)


Commonly used basic material

- Thermoplastic elastomer (TPE)

When customer wants to have a softer material

### **Box Tip Construction:**

The tip should not have a seam that interferes with the ultrasonic image



## 5. Probe Covers





[www.onyxmedical.com](http://www.onyxmedical.com)

### High Grade Ultrasonic Gel

Purpose is to use inside the probe cover and enhances the Ultrasound image quality.

Can however also be used outside the probe cover and even intra-operative.



## 5. Probe Covers



[www.onyxmedical.com](http://www.onyxmedical.com)



Probe Covers can have tape for a tight fit and latex free rubber bands.



No Gel needed, probe can be fixed with tape.

This tape is easy to remove and guarantees a clear ultrasonic image.



6. Miscellaneous

Back Table Covers



Product Number	Description	Size (W x L)	Cross References		
			Cardinal	Kimberly Clark	Medline
E6001	Back Table Cover	44" x 75" (24" x 75")	A8376N	79445 / 79448	DYNJP2309
E6001FF	Back Table Cover	44" x 75" (24" x 75")	A8376N	79445 / 79448	DYNJP2309
E6002	Back Table Cover	50" x 90" (30" x 90")	A8378N	79056 / 79058	DYNJP2316 / DYNJP2319
E6002FF	Back Table Cover	50" x 90" (30" x 90")	A8378N	79056 / 79058	DYNJP2316 / DYNJP2319

E6001 Features

~ Absorbent Reinforcement 24" x 75"

E6001FF Features

~ Fan Folded

~ Absorbent Reinforcement 24" x 75"

E6002 Features

~ Absorbent Reinforcement 30" x 90"


E6002FF Features

~ Fan Folded

~ Absorbent Reinforcement 30" x 90"

6. Miscellaneous

Back Table Covers



Product Number	Description	Size (W x L)	Cross References		
			Cardinal	Kimberly Clark	Medline
E6003	Back Table Cover	44" x 90" (24" x 90")	A8377N	79498	DYNJP2311
E6004	Back Table Cover	44" x 44" (24" x 44")		79440	

E6003 Features

~ Absorbent Reinforcement 24" x 90"

(Also comes Fan Folded)

E6004 Features

~ Absorbent Reinforcement 24" x 44"

(Also comes Fan Folded)



## 6. Miscellaneous



[www.onyxmedical.com](http://www.onyxmedical.com)

### **Surgical Gowns**

- Sterile packed and disposable
- Easily activated
- Can be used to warm two endoscopes at a time
- Scope stays warm for 45 minutes to one hour
- Easy to use different size scopes.
- No alcohol
- No need for costly decontamination, sterilization, and additional sterile water for Thermos-type warmers
- Disposable & environmentally friendly
- Sterile
- Latex and DEHP free
- CE marked

## 6. Miscellaneous

# Surgical Gowns – Standard Non-Reinforced Gown


AAMI Level	1	2	3	4
Sleeve Types:	Set-In / Raglan (by special request)			

### Features

~ adjustable velcro & pile on gown neck

### Use

~ for lower fluid level procedures



Regular SMS /  
Soft SMS



Fabric  
Reinforcement



Poly  
Reinforcement

Product Number	Description	Size (W x L)	Cross References		
			Cardinal	Kimberly Clark	Medline
EP0121	Unreinforced SMS Gown	L	9515N	72878 / 90018	DYNJP2001S / SPT2001SMS
EP0131	Unreinforced SMS Gown	XL	9545N	72888 / 90048	DYNJP2002S / SPT2002SMS
EP0141	Unreinforced SMS Gown	XXL	9575N		DYNJP2003S / SPT2003SMS
E0121	Unreinforced Soft SMS Gown	L	9518		DYNJP2701 / SPT2701
E0131	Unreinforced Soft SMS Gown	XL	9548		DYNJP2702 / SPT2702
E0141	Unreinforced Soft SMS Gown	XXL	9578		DYNJP2703 / SPT2703

## 6. Miscellaneous

# Surgical Gowns – Fabric Reinforced Gown

AAMI Level	1	2	3	4
Sleeve Types:	Set-In / Raglan (by special request)			


### Features

~ adjustable velcro & pile on gown neck

~ fabric reinforced

### Use

~ for medium fluid level procedures



Regular SMS /  
Soft SMS



Fabric  
Reinforcement



Poly  
Reinforcement

Product Number	Description	Size (W x L)	Cross References		
			Cardinal	Kimberly Clark	Medline
EP0221	Fabric Reinforced SMS Gown	L	9511N	78911 / 78917 / 90112	DYNJP2101S / SPT2101SMS
EP0231	Fabric Reinforced SMS Gown	XL	9541N	78921 / 78927 / 90142	DYNJP2102S / SPT2102SMS
EP0241	Fabric Reinforced SMS Gown	XXL	9571N		DYNJP2103S / SPT2103SMS
E0221	Soft SMS Fabric Reinforced Gown	L			DYNJP2704 / SPT2704
E0231	Soft SMS Fabric Reinforced Gown	XL			DYNJP2705 / SPT2705
E0241	Soft SMS Fabric Reinforced Gown	XXL			DYNJP2706 / SPT2706

## 6. Miscellaneous

# Surgical Gowns – Poly Reinforced Gown

AAMI Level	1	2	3	4
Sleeve Types:	Set-In / Raglan (by special request)			


### Features

~ adjustable velcro & pile on gown neck

~ poly reinforced

### Use

~ for heavier fluid level procedures



Regular SMS /  
Soft SMS



Fabric  
Reinforcement



Poly  
Reinforcement

Product Number	Description	Size (W x L)	Cross References		
			Cardinal	Kimberly Clark	Medline
EP0321	Poly Reinforced SMS Gown	L	9010 / 9010N	72311	
EP0331	Poly Reinforced SMS Gown	XL	9040 / 9040N	72321	
EP0341	Poly Reinforced SMS Gown	XXL	9070 / 9070N		
E0321	Soft SMS Poly Reinforced Gown	L			
E0331	Soft SMS Poly Reinforced Gown	XL			
E0341	Soft SMS Poly Reinforced Gown	XXL			

## 6. Miscellaneous

# Surgical Gowns – AAMI 4 Gown

AAMI Level	1	2	3	4
Sleeve Types:	Set-In / Raglan (by special request)			

### Features

~ adjustable velcro & pile on gown neck


~ breathable material

### Use

~ for heavier fluid level procedures



AAMI 4 Breathable Material



Product Number	Description	Size (W x L)	Cross References		
			Cardinal	Kimberly Clark	Medline
E40121	AAMI 4 Gown	L	9011 / 9011N	72438	DYNJP2306P
E40131	AAMI 4 Gown	XL	9041 / 9041N	72445	DYNJP2307P
E40141	AAMI 4 Gown	XXL	9071 / 9071N	72447	DYNJP2308P
E40151	AAMI 4 Gown	XL, X-Long	9041EL / 9041ELN	73041	DYNJP2303P

## 7. Face Shields



[www.onyxmedical.com](http://www.onyxmedical.com)



### **Standard Shield (Model 16000)**

9 x 13 inch full face  
protection



### **Drape Shield (Model 18000) Short (7 inch)**

fluid barrier fabric  
below Standard Shield  
Fabric stays outside of  
protective garments.  
Can cup under chin.



### **Drape-U Shield (Model 20000-100)**

Full face, head and neck protection with  
AAAM Level III barrier fabric  
Adjustable side fabric flaps open/close





[www.onyxmedical.com](http://www.onyxmedical.com)



[www.onyxmedical.com](http://www.onyxmedical.com)

**Onyx Medical Inc.  
16 Digital Drive, Suite 120  
Novato, CA 94949**

**Toll-Free:  
800-333-5773  
Phone  
415-884-4412**

**EMAIL US at  
[info@onyxmedical.com](mailto:info@onyxmedical.com)**