

PRODUCT OVERVIEW

Implants



Selected plates, screws and instruments for osteosynthesis

„KÖNIGSEE IMPLANTATE is your specialist for plate systems with standard or (variable) angle-stable screw connections for all extremities.

The products we manufacture in Germany and sell worldwide are specially tailored to the needs adapted to an osteosynthesis and guarantee doctor and patient a smooth operation with the aim of improving the patient's quality of life.“

The product overview gives you an overview of our range. If you are missing a specific implant, please contact our service team or your local sales representative for advice on our products.

STERILE IMPLANTS

All implants are available both non-sterile packed and sterile packed. To order sterile packed implants, please complete the Article number with a leading „S“.

Example: 3.133.14 not packed sterile → S3.133.14 packed sterile.

SETS OF IMPLANTS

Our implant sets usually consist of implants, the associated screws and instruments. Based on these complete sets, we also offer sets with sterile packed implants. You can recognize the sets with sterile plates and screws analogous to the original set by a „PS +“ in the article number. Some sets are also available in a reduced version, with a minimum number of sterile plates and screws. You can recognize this by the „PS“ at the end.

Example:

19.104.00 non-sterile →

19.104.00PS+ (plates + screws packed sterile)

19.045.00I (associated instruments)

INSTRUMENTS

The instrument sets for using the PS and PS + variants can be found on page 43 or in the appropriate places in the product overview.

You can usually recognize sets by the „I“ at the end of the article number.

MATERIAL

The article numbers at Königsee Implantate differentiate between implant steel and titanium, starting with an even or an odd number.

Example screws: 3.133.14 titanium → 2.133.14 steel

Example plates: 5.799.06 titanium → 4.799.06 steel

Our plates and nails are made of pure titanium, screws are made of titanium alloy. Deviations in plates or nails can be found directly with the product.

INFORMATIONS

WiKo®

„The system name WiKo® is derived from the German designation „Winkelstabile-Kombination“, which translates into English as „angle-stable combination“. This system combines conventional plating technologies with angle-stable locking.

WiKo® is a registered trademark of Königsee Implantate GmbH.“

VaWiKo®

„The system name VaWiKo® is derived from the German designation „Variabel winkelstabile Kombination“, which translates into English as „variable angle-stable combination“. This system combines conventional plating technologies with angle-stable locking in variable angles.

VaWiKo® is a registered trademark of Königsee Implantate GmbH.“

TiOB®

TiOB® is a surface treatment that brings optimized properties.

ABBREVIATIONS

st	self tapping
sb	self drilling
bt	back tapping
IH	Internal Hexagon
SL	Six Lobe
as	angle-stable
vas	variable angle-stable
T	titanium
S	steel
I	instruments
PS+	plates + screws original set
PS	plates + screws minimum equipment

	4	Clavicle
	6	Humerus
	9	Elbow
	10	Radius
	15	Ulna
	15	Hand
	18	Femur
	21	Patella
	21	Tibia
	25	Fibula
	25	Foot
	31	Spine
	31	Periprosthetic
	33	Children
	35	Wires
	35	Cannulated
	38	Screws
	39	Systems
	43	Instruments

Clavicle shaft plate 3.5 as anterior



Clavicle shaft plate 3.5 as anterior med.

- Material: titanium alloy
- Surface: TiOB®
- Holes: 6 – 8
- Lengths: 80 – 103 mm
- Item no.: 5.3802.XX
- Set Instruments: 19.035.00I



Clavicle shaft plate 3.5 as anterior lat.

- Material: titanium alloy
- Surface: TiOB®
- Holes: 8 – 10
- Lengths: 105 – 128 mm
- Item no.: 5.3812.XX
- Set Instruments: 19.035.00I

Clavicle plates 3.5 as



Lateral clavicle plate 3.5 as

- Head holes: 3
- Shaft holes: 4 | 6
- Lengths: 65 | 82 mm
- Item no. left: 5.3512.XX
- Item no. right: 5.3502.XX
- Set Instruments: 19.035.00I




Clavicle shaft plate 3.5 as

- Holes: 6 – 10
- Lengths: 99 | 124 mm
- Item no. left: 5.3532.XX
- Item no. right: 5.3522.XX
- Set Instruments: 19.035.00I



ACJ hook plate 3.5 as

- Head holes: 2
- Shaft holes: 3
- Length: 78 mm
- Offset: 14 | 18 mm
- Hook: 18 mm
- Item no. left: 5.3552.XX18
- Item no. right: 5.3542.XX18
- Set Instruments: 19.035.00I



Clavicle plate 3.5 as IH T SET

- angle-stable, Small fragment
- Set no.: 19.478.00

Internal hexagon



Cortical screw 3.5 as st

- Item no.: 3.133.XX



Standard Cortical screw 3.5 st

- Item no.: 3.132.XX

Lateral clavicle plate 2.7 vas

Lateral clavicle plate 2.7 vas

- Head holes: 5
- Shaft holes: 4 | 6
- Lengths: 65 and 82 mm
- Item no. left: 5.3613.XX
- Item no. right: 5.3603.XX
- Set Instruments: 19.027.00I



Lateral clavicle plate 3.5 vas

Lateral clavicle plate 3.5 vas

- Head holes: 5
- Shaft holes: 4 | 6
- Lengths: 65 and 82 mm
- Item no. left: 5.3513.XX
- Item no. right: 5.3503.XX
- Set Instruments: 19.035.00I



Clavicle shaft plate 3.5 vas R95

- Holes: 6 – 10
- Lengths: 74; 99 and 124 mm
- Item no. left: 5.3533.XX
- Item no. right: 5.3523.XX
- Set Instruments: 19.035.00I



Clavicle shaft plate 3.5 vas R70

- Holes: 9
- Length: 107 mm
- Item no. left: 5.3533.09/70
- Item no. right: 5.3523.09/70
- Set Instruments: 19.035.00I



Clavicle shaft plate 3.5 vas R105

- Holes: 9
- Length: 111 mm
- Item no. left: 5.3533.09/105
- Item no. right: 5.3523.09/105
- Set Instruments: 19.035.00I



ACJ hook plate 3.5 vas

- Head holes: 2
- Shaft holes: 2 – 6
- Lengths: 66; 74 and 97 mm
- Offset: 14 mm
- Hook: 16 mm
- Item no. left: 5.3573.XX/1416
- Item no. right: 5.3563.XX/1416
- Set Instruments: 19.035.00I



Cortical screw 3.5 as st

- Item no.: 3.133.XX




Standard Cortical screw 3.5 st

- Item no.: 3.132.XX

Clavicle plate 3.5 vas IH T SET

- variable angle-stable, Small fragment
- Set no.: 19.479.00



Internal hexagon

Ligamentfix



Ligamentfix holes 4 length 12
· Item no.: 5.990.01



Ligamentfix holes 4 length 15
· Item no.: 5.990.02



Ligamentfix holes 5 length 15
· Item no.: 5.990.05



Ligamentfix holes 6 length 18
· Item no.: 5.990.06

Screw Anchor

· Material: titanium alloy



Screw LASA-DR 6.0
· Lengths: 30 | 35 | 40 mm
· Item no.: 3.734.XX



Screw TESA 5.0
· Item no.: 3.7302.30



Screw TRANSFIX 5.0 AC joint
· Item no.: 3.735.22



Screw TRANSFIX 5.0 elbow
· Item no.: 3.7351.22



Screw Anchor IH T SET
· Aiming device system
· Set no.: 19.503.00

for shoulder- and
elbow surgery

Proximal humerus plate FIVE 4.0/3.5 as/vas



Proximal humerus plate Five 3.5/4.0 as
· Head holes: 5
· Shaft holes: 3 – 14
· Lengths: 74 – 218 mm
· Item no. titanium: 5.799.XX
· Item no. steel: 4.799.XX
· Set Instruments: 19.035.00I



Proximal humerus plate Five 3.5/4.0 vas
· Head holes: 5
· Shaft holes: 3 – 14
· Lengths: 74 – 218 mm
· Item no.: 5.797.XX
· Set Instruments: 19.035.00I



**Guide block steel for
prox. humerus pl. FIVE**
· Item no.: 10.173.00

Proximal humerus plate FIVE 4.0/3.5 as/vas



Standard Cortical screw 3.5 st
· Item no.: 3.132.XX



Cortical screw 3.5 as st
· Item no.: 3.133.XX




Cancellous screw 4.0 as st
· Item no.: 3.147.XX

Prox.Humeruspl. FIVE as IH SET

- Set no.: titanium: 19.364.00
- Set no.: steel: 19.365.00

Prox.Humeruspl. FIVE vas IH T SET

- Set no.: titanium: 19.367.00



Internal hexagon

Proximal humerus plate EIGHT 4.0/3.5 vas



Guide block PEEK for
prox. humerus pl. EIGHT
· Item no.: 10.179.10

Proximal humerus plate EIGHT 4.0/3.5 vas

- Surface: TiOB®
- Head holes: 8
- Shaft holes: 2 – 14
- Lengths: 61 – 253 mm
- Item no.: 5.790.XX
- Set Instruments: 19.035.00I



Standard Cortical screw 3.5 st
· Item no.: 3.132.XX




Cortical screw 3.5 as st
· Item no.: 3.133.XX



Cancellous screw 4.0 as st
· Item no.: 3.147.XX

Prox.Humerus plate EIGHT vas IH T SET

- Set no.: 19.3673.00



Internal hexagon



Standard Cortical screw 3.5 st
· Item no.: 3.132.XX




Cortical screw 3.5 as st
· Item no.: 3.133.XX



Cortical screw 4.0 as st
· Item no.: 3.148.XX

Prox.Humerus plate EIGHT vas Cortical IH T SET

- Set no.: 19.3683.00



Internal hexagon

VaWiKo® proximal humerus plate 4.0/3.5




Guide block for prox. humerus pl.
VaWiKo®
· Item no. 10.179.15

VaWiKo® proximal humerus plate 4.0/3.5

- Surface: TiOB®
- Head holes: 8
- Shaft holes: 3 – 8
- Lengths: 77 – 157 mm
- Item no.: 5.07913.8X
- Set Instruments: 19.035.00I



VaWiKo® proximal humerus plate 4.0/3.5



Internal hexagon

VaWiKo® Prox.humerus plate 3.5 IH SET

· Set no.: 19.3751.00




Standard Cortical screw 3.5 st

· Item no.: 3.132.XX



Cortical screw 3.5 as st

· Item no.: 3.133.XX



Internal hexagon

VaWiKo® Prox.humerus 4.0/3.5 IH Cortical SET

· Set no.: 19.3752.00



Standard Cortical screw 3.5 st

· Item no.: 3.132.XX




Cortical screw 3.5 as st

· Item no.: 3.133.XX



Cortical screw 4.0 as st

· Item no.: 3.148.XX



Internal hexagon

VaWiKo® Prox.humerus 4.0/3.5 IH Cancellous SET

· Set no.: 19.3753.00



Standard Cortical screw 3.5 st

· Item no.: 3.132.XX



Cortical screw 3.5 as st


· Item no.: 3.133.XX



Cancellous screw 4.0 as st

· Item no.: 3.147.XX

VaWiKo® Prox.humerus aiming device INSTRUMENTS



VaWiKo® Prox.humerus aiming device INSTRUMENTS

· Set no.: 19.375.00

Proximal humerus Plate 6.0/4.5 as



Proximal humerus Plate 6.0/4.5 as

- Head holes: 3
- Shaft holes: 5 – 16
- Lengths: 73 – 217 mm
- Item no. steel: 4.795.XX
- Item no. titanium: 5.795.XX
- Set Instruments: 19.045.00I

Proximal humerus Plate 6.0/4.5 as



Standard Cortical screw 4.5 st
· Item no.: 3.152.XX




Cortical screw 4.5 as st
· Item no.: 3.545.XX



Cancellous screw 6.0 as st
· Item no.: 3.163.XX

Prox.Humerus plate 6.0/4.5 as IH SET

- Set no. titanium: 19.104.00
- Set no. steel: 19.103.00



Internal hexagon

Distal humerus plate Y 3.5/4.5 as



Standard Cortical screw 3.5 st
· Item no.: 3.132.XX



Cortical screw 3.5 as st
· Item no.: 3.133.XX



Cortical screw 4.5 as st
· Item no.: 3.545.XX



Standard Cortical screw 4.5 st
· Item no.: 3.152.XX


Distal humerus plate Y 3.5/4.5 as

- Head holes: 8
- Shaft holes: 5 – 8
- Lengths: 141 – 199 mm
- Item no.: 5.8591.8X



Distal Humerus plate 3.5/4.5 as IH T SET

- Set no.: 19.675.010



Internal hexagon

Distal humerus plates 3.5 as/vas



Distal humerus Plate 3.5 as radial without extension

- Head holes: 4
- Shaft holes: 4 – 10
- Lengths: 78 – 150 mm
- Item no. left: 5.8552.4X
- Item no. right: 5.8542.4X
- Set Instruments: 19.035.001

Distal humerus Plate 3.5 as radial with extension

- Head holes: 5
- Shaft holes: 4 – 10
- Lengths: 78 – 150 mm
- Item no. left: 5.8532.5X
- Item no. right: 5.8522.5X
- Set Instruments: 19.035.001



Distal humerus plate 3.5 vas ulnar

- Head holes: 4
- Shaft holes: 4 – 10
- Lengths: 78 – 150 mm
- Item no. left: 5.8573.4X
- Item no. right: 5.8563.4X
- Set Instruments: 19.035.001



Olecranon plates




Olecranon plate 3.5 vas without claws

- Holes: 8 | 10 | 14
- Lengths: 76 | 99 | 151 mm
- Item no.: 5.08753.XX
- Set Instruments: 19.035.00I



Olecranon plate 3.5 vas with claws

- Holes: 8 | 10 | 14
- Lengths: 76 | 99 | 151 mm
- Item no.: 5.08743.XX
- Set Instruments: 19.035.00I



Elbow 3.5 as IH T SET

Distal humerus plates 3.5 as, vas and Olecranon plates

- Set no.: 19.551.00

Internal hexagon



Standard Cortical screw 3.5 st

- Item no.: 3.132.XX



Cortical screw 3.5 as st


- Item no.: 3.133.XX

Proximal radius Plate 2.0 as



Proximal radius Plate 2.0 as

- Head holes as: 2
- Shaft holes: 3 | 4
- Lengths: 32 | 40 mm
- Item no.: 5.884.2X



Dist. & prox. Radius plate 2.0/3.5 as IH T SET

- Set no.: 19.436.10

Internal hexagon



Cortical screw 2.0 as st

- Depth stop
- Item no.: 3.114.XX



Standard Cortical screw 2.0 st

- Item no.: 3.110.XX



Standard Cortical screw 3.5 st

- Item no.: 3.132.XX



Distal Radius plate 2.0/3.5 as IH T SET

- Set no.: 19.436.00
- Set Instruments: 19.035.00I

Internal hexagon



Cortical screw 2.0 as st

- Depth stopp
- Item no.: 3.114.XX



Standard Cortical screw 3.5 st

- Item no.: 3.132.XX

Distal radius plates 2.0/3.5 as

Screw connection:

Head angle-stable
 Shaft standard

palmar plate

Distal palmar radius plate 2.0/3.5 as

- Head holes: 6
- Shaft holes: 3 | 5
- Lengths: 50 | 67
- Item no. left: 5.899.6X
- Item no. right: 5.899.6X



Distal palmar radius plate 2.0/3.5 as LH middle

- Head holes: 6
- Shaft holes: 3
- Length: 50 mm
- Item no. left: 5.8995.63
- Item no. right: 5.8985.63



dorsal plate

Distal radius plate 2.0/3.5 as dorsal broad

- Head holes: 5
- Shaft holes: 3
- Length: 63 mm
- Item no. left: 5.893.53
- Item no. right: 5.892.53



Distal radius plate 2.0/3.5 as dorsal

- Head holes: 4
- Shaft holes: 3
- Length: 63 mm
- Item no. left: 5.893.43
- Item no. right: 5.892.43



Distal radius plate 2.4 vas

Distal radius plate 2.4 vas narrow

- Head holes: 7
- Shaft holes: 2–12
- Width: 23 mm
- Lengths: 46–131 mm
- Item no. left: 5.08173.XX
- Item no. right: 5.08273.XX
- Set Instruments: 19.024.00I



Distal radius plate 2.4 vas broad

- Head holes: 9
- Shaft holes: 2–12
- Width: 28 mm
- Lengths: 46–131 mm
- Item no. left: 5.08193.XX
- Item no. right: 5.08293.XX
- Set Instruments: 19.024.00I



Distal radius plates 2.4 vas



★ Six lobe

Distal radius plate 2.4 vas SL T SET
Set no.: 19.117.00



Cortical screw 2.4 SL8 st
Item no.: 3.0224.XX



Cortical screw 2.4 as SL8 st TiOB®
Item no.: 3.0124.XX



Pin 1.8 as SL8; TiOB®
Item no.: 3.0324.XX

Distal radius plates 2.4 Watershed vas



Distal radius plate 2.4 Watershed vas

- Head holes: 7
- Shaft holes: 3–4
- Width: 23 mm
- Lengths: 50 – 59 mm
- Item no. left: 5.08373.XX
- Item no. right: 5.08473.XX
- Set Instruments: 19.024.00I



★ Six lobe

Dist. radius plate 2.4 Watershed vas
SL T Set
Set no.: 19.117.10



Cortical screw 2.4 SL8 st
Item no.: 3.0224.XX



Cortical screw 2.4 as SL8 st TiOB®
Item no.: 3.0124.XX



Pin 1.8 as SL8; TiOB®
Item no.: 3.0324.XX

Distal radius plates 2.7 as



Distal radius plate 2.7 vas broad

- Head holes: 7
- Shaft holes: 2 – 13
- Width: 27 mm
- Lengths: 43 – 136 mm
- Item no. left: 5.9852.XX
- Item no. right: 5.9842.XX
- Set Instruments: 19.027.00I



Guide block PEEK Radius pl. 2.7 broad

- Item no. left: 10.181.07
- Item no. right: 10.180.07

Distal radius plates 2.7 as



Guide block PEEK Radius pl. 2.7 schmal

- Item no. left: 10.181.05
- Item no. right: 10.180.05

Distal radius plate 2.7 as narrow

- Head holes: 5
- Shaft holes: 2 – 13
- Width: 23 mm
- Lengths: 43 – 136 mm
- Item no. left: 5.9832.XX
- Item no. right: 5.9822.XX
- Set Instruments: 19.027.00I



Standard Cortical screw 2.7 st

- Item no.: 3.120.XX



Cortical screw 2.7 as st

- Item no.: 3.125.XX

Dist. radius plate 2.7 as IH T SET

- Set no.: 19.116.00



Internal hexagon

Distal radius plates 2.7 vas



Guide block PEEK Radius pl. 2.7 broad

- Item no. left: 10.181.07
- Item no. right: 10.180.07

Distal radius plate 2.7 vas broad

- Head holes: 7
- Shaft holes: 2 – 13
- Width: 27 mm
- Lengths: 43 – 136 mm
- Item no. left: 5.9853.XX
- Item no. right: 5.9843.XX
- Set Instruments IH: 19.027.00I
- Set Instruments SL: 19.027.10I



Guide block PEEK Radius pl. 2.7 narrow

- Item no. left: 10.181.05
- Item no. right: 10.180.05

Distal radius plate 2.7 vas narrow

- Head holes: 5
- Shaft holes: 2 – 13
- Width: 23 mm
- Lengths: 43 – 136 mm
- Item no. left: 5.9833.XX
- Item no. right: 5.9823.XX
- Set Instruments IH: 19.027.00I
- Set Instruments SL: 19.027.10I



Standard Cortical screw 2.7 st

- Item no.: 3.120.XX



Cortical screw 2.7 as st

- Item no.: 3.125.XX

Dist. radius plate 2.7 vas IH T SET

- Set no.: 19.1163.00



Internal hexagon

Distal radius plates 2.7 vas



Dist. radius plate 2.7 vas SL T SET
· Set no.: 19.1161.00

★ Six lobe



Cortical screw 2.7 flat head SL10 st
· Item no.: 3.1281.XX



Cortical screw 2.7 as SL10 st
· Item no.: 3.1251.XX

Distal radius plates 3.5 as/vas



Distal radius plate 3.5 Vario as
· Surface: TiOB®
· Head holes: 3
· Shaft holes: 3 – 13
· Lengths: 51 – 136 mm
· Item no. left: 5.887.3X
· Item no. right: 5.886.3X
· Set Instruments: 19.035.00I



Distal radius plate 3.5 Classic as
· Head holes: 3
· Shaft holes: 2 – 13
· Lengths: 42 – 136 mm
· Item no. left: 5.877.4X
· Item no. right: 5.876.4X
· Set Instruments: 19.035.00I



Distal radius plate 3.5 Minor as
· Head holes: 3
· Shaft holes: 3 – 8
· Lengths: 51 – 93 mm
· Item no. left: 5.883.3X
· Item no. right: 5.882.3X
· Set Instruments: 19.035.00I




Distal radius plate 3.5 Quadro as
· Head holes: 4
· Shaft holes: 3 – 8
· Lengths: 51 – 93 mm
· Item no. left: 5.877.4X
· Item no. right: 5.876.4X
· Set Instruments: 19.035.00I



Distal radius plate 3.5 Classic vas
· Head holes: 3
· Shaft holes: 3 – 13
· Lengths: 51 – 136 mm
· Item no. left: 5.8773.3X
· Item no. right: 5.8763.3X
· Set Instruments: 19.035.00I



Distal radius plate 3.5 Minor vas
· Head holes: 3
· Shaft holes: 3 – 13
· Lengths: 51 – 136 mm
· Item no. left: 5.8833.3X
· Item no. right: 5.8823.3X
· Set Instruments: 19.035.00I



Dist. radius plate 3.5 as IH T SET
· Set no.: 19.112.00



Standard Cortical screw 3.5 st
· Item no.: 3.132.XX



Cortical screw 3.5 as st
· Item no.: 3.133.XX

Internal hexagon

Distal ulna plate 2.0 as

Distal ulna plate 2.0 as

- Head holes: 4
- Shaft holes: 3 – 6
- Lengths: 38 – 60 mm
- Item no. left: 5.889.4X
- Item no. right: 5.888.4X


Distal Ulna plate 2.0 as IH T SET
 Set no.: 19.506.00


Internal hexagon


Standard Cortical screw 2.0 st
 · Item no.: 3.110.XX

Cortical screw 2.0 as st
 · Depth stop
 · Item no.: 3.114.XX

Hand plates Mini 2.0

Straight plate Mini 2.0 as

- Holes: 5 | 13
- Lengths: 29 | 77 mm
- Item no.: 5.9012.XX



Straight plate Mini 2.0 as round bridge

- Holes: 4
- Length: 29 mm
- Item no.: 5.9022.04



K plate Mini 2.0 as

- Holes: 6
- Length: 37 mm
- Item no.: 5.9032.06



Double bridge plate 2.0 as skew

- Holes: 4
- Length: 15 mm
- Item no.: 5.9052.04S



Double bridge plate 2.0 as skew

- Holes: 7
- Length: 25 mm
- Item no.: 5.9052.07S



Double bridge plate 2.0 as

- Holes: 10
- Length: 45 mm
- Item no.: 5.9052.10



Hand plates Mini 2.0

**T plate Mini 2.0 as broad**

- Head holes: 3
- Shaft holes: 4 – 8
- Lengths: 30; 42; 54 mm
- Item no.: 5.9332.XX

**T plate Mini 2.0 as narrow**

- Head holes: 2
- Shaft holes: 4 – 8
- Lengths: 30; 42; 54 mm
- Item no.: 5.9312.XX

**T plate Mini 2.0 as**

- Head holes: 3
- Shaft holes: 4
- Length: 33 mm
- Item no.: 5.9342.34

**Y-plate Mini 2.0 as**

- Head holes: 3
- Shaft holes: 4 – 8
- Lengths: 33; 45; 57 mm
- Item no.: 5.9322.XX

**MF plate as IH T SET**

- Set no.: 19.508.00

 Internal hexagon

**Cortical screw 2.0 as st**

- Depth stop
- Item no.: 3.114.XX

Ring plates vas

**Ring plate 2.0 vas**

- Holes: 7
- Ø 14,0 mm
- Item no.: 5.9953.14

**Ring plate 2.7 vas**

- Holes: 5
- Ø 16,5 mm
- Item no.: 5.9953.16

for arthrodesis in the
carpal area

Ring plates vas



Standard Cortical screw 2.0 st
· Item no.: 3.110.XX



Cortical screw 2.0 as st
· Depth stop
· Item no.: 3.114.XX



Cortical screw 2.7 as st
· Item no.: 3.125.XX

Ring plate vas IH T SET
· Set no.: 19.439.00



Internal hexagon

Hand Arthrodesis plate vas

Hand Arthrodesis plate 2.7/3.5 vas

- Offset 0 | 2
4° | 5°
- Head holes: 3
- Shaft holes: 3
- Lengths: 96 | 97 mm
- Item no. 5.9963.33XX



Hand Arthrodesis plate 2.7/3.5 vas

- Offset 0
4°
- Head holes: 4
- Shaft holes: 5
- Lengths: 133 mm
- Item no. 5.9963.4540



for revisions



Standard Cortical screw 2.7
flat head st
· Item no.: 3.128.XX



Cortical screw 2.7 as st
· Item no.: 3.125.XX




Standard Cortical screw 3.5 st head 5
· TiOB®
· Item no.: 3.138.XX



Cortical screw 3.5 as st head 5
· TiOB®
· Item no.: 3.136.XX

Hand Arthrodesis plate 2.7/3.5 vas
IH T SET

- Set no.: 19.671.020



Internal hexagon

RoSA® – Rotationally stable screw anchor



RoSA screw Ø10

- Lengths: 50 | 70 | 90 mm
- Item no.: 7.306.XX



RoSA blade

- Lengths: 60 | 80 | 100 mm
- Item no.: 7.307.XX



RoSA sliding sleeve plate


- Angle: 122° | 129° | 136° | 143°
- Lengths: 28 | 60 | 92 mm
- Item no.: 7.308.XXXPX



RoSA® connecting screw

- Item no.: 7.309.01 (cannulated)
- Item no.: 7.309.02 (not cannulated)

Angular stable implant system for bone head-conserving osteosynthesis



RoSA IH T SET

- Set no.: 19.606.00

Internal hexagon



Standard Cortical screw 4.5 st

- Item no.: 3.152.XX



Cortical screw 4.5 as st

- Item no.: 3.545.XX




Internal hexagon



RoSA IH INSTRUMENTS

- Set no.: 19.608.00 (2 Levels)



Internal hexagon

RoSA IH INSTRUMENTS Removal

- Set no.: 19.610.00

RoSA trochanter stabilization plate 4.5 vas

- Surface: TiOB®
- Lengths neck: 0 | 4 | 8



- 2 Head holes:
- Length: 64 mm
- Item no.: 7.304.XX



- 4 Head holes:
- Length: 117 mm
- Item no.: 7.305.XX

RoSA IH T Trochanter stabilization
plate SET
Set no.: 19.612.00 (2 Levels)



Cortical screw 4.5 as st
· Item no.: 3.545.XX



Internal hexagon

DHS/DCS

DHS flange

- Angle: 135° | 140° | 145° | 150°
- Lengths: 58 – 250 mm
- Item no. titanium 135°: 7.315.XX
- Item no. steel 135°: 6.315.XX
- Smaller angles available in steel



DHS/DCS screw Ø12.5

- Lengths: 50 – 130 mm
- Item no. titanium: 7.310.XX
- Item no. steel: 6.310.XX



DHS compression screw

- Item no. titanium: 7.311.01
- Item no. steel: 6.311.01




Support plate 3.5/4.5

- Item no. titanium: 7.395.00
- Item no. steel: 6.395.00



DHS IH SET

- Set no. titanium: 19.122.00
- Set no. steel: 19.121.00



Internal hexagon




Internal hexagon

DHS/DCS IH INSTRUMENTS Removal


Set no.: 19.125.00

DHS/DCS IH INSTRUMENTS

Set no.: 19.120.00



Internal hexagon



Internal hexagon

DHS Screws IH SET

- Set no. steel: 19.168.10
- Set no. titanium: 19.168.30

Femur nail from Königsee (KFN)

**KFN 125°**

- Length: 210 mm
- Ø: 10 – 13 mm
- Inside-Ø: 4.7 mm
- Item no.: 7.3KF.XX/210

**KFN 130°**


- Length: 210 mm
- Ø: 10 – 13 mm
- Inside-Ø: 4.7 mm
- Item no.: 7.1KF.XX/210

**KFN locking screw 5.5 st**

- Surface: TiOB®
- Length: 25 – 65 mm
- Item no.: 3.157.XX


**KFN rotation lock screw 6.0**

- Surface: TiOB®
- Length: 65 – 120 mm
- Item no.: 3.428.XX

**KFN IH T SET**

- Set no.: 19.344.00

Internal hexagon

**KFN INSTRUMENTS Removal**

- Set no.: 19.344.10

Internal hexagon

KFN 125° Ø10.5

- Length: 340 – 380 mm
- Ø: 10.5 mm
- Inside-Ø: 4.7 mm
- Item no. left: 7.3KF10,5L/XXX
- Item no. right: 7.3KF10,5R/XXX

KFN 130° Ø10.5

- Length: 340 – 420 mm
- Ø: 10.5 mm
- Inside-Ø: 4.7 mm
- Item no. left: 7.1KF10,5L/XXX
- Item no. right: 7.1KF10,5R/XXX

**KFN clasp screw**

- Item no.: 7.311.12

**KFN clasp screw with locking function**

- Item no.: 7.311.30

**Compression screw KFN**

- Item no.: 7.311.01

**KFN sleeve for support screw 12.0**

- Item no.: 3.310.52

**KFN support screw 10.0**

- Length: 70 – 130 mm
- Item no.: 7.313.XX

Distal femur plate

Distal femur plate 4.5-6.0 as

- Head holes: 7
- Shaft holes: 4 – 16
- Lengths: 118 – 309 mm

**titanium**

- Item no. left: 5.9492.XX
- Item no. right: 5.9482.XX
- Set Instruments: 19.045.00I

steel

- Item no. left: 4.9492.XX
- Item no. right: 4.9482.XX
- Set Instruments: 19.045.00I

Distal femur plate 4.5-6.0 as



Standard Cortical screw 4.5 st
· Item no.: 3.152.XX



Cortical screw 4.5 as st
· Item no.: 3.545.XX



Cortical screw 5.5 as st
· Item no.: 3.549.XX

Distal Femur plate 4.5/5.5 as IH
SET


- Set no. titanium: 19.602.00
- Set no. steel: 19.600.00



Internal hexagon

Distal Femur plate 4.5/5.5 as IH
aiming device INSTRUMENTS

- Set no.: 19.604.00



Internal hexagon

Patella plate 3.5 holeØ5 vas

Patella plate 3.5 holeØ5 vas


- Holes: 11
- Item no.: 5.8513.11
- Set Instruments: 19.035.10I



Cortical screw 3.5 as st head 5
· TiOB®
· Item no.: 3.136.XX

Patella plate 3.5 vas IH T SET

- Set no.: 19.671.010



Internal hexagon

VaWiKo® tibia head plate 4.0

VaWiKo® tibia head plate 4.0 broad

- Head holes: 7
- Width: 40 mm
- 4 Shaft holes; Length: 88 mm
· Item no.: 5.7083.704
- 8 Shaft holes; Length: 140 mm
· Item no.: 5.7083.708
- Set Instruments: 19.493.20I



VaWiKo® tibia head plate 4.0




VaWiKo® tibia head plate 4.0 narrow

- Head holes: 6
- Width: 35 mm

4 Shaft holes; Length: 88 mm
· Item no.: 5.7073.604

8 Shaft holes; Length: 140 mm
· Item no.: 5.7073.608

· Set Instruments: 19.493.20I



VaWiKo® Tibia head plate 4.0 SL T SET

Set no.: 19.493.20

★ Six lobe



Pestle INSTRUMENTS

Set no.: 19.493.30



Cortical screw 4.0 as ISR15 st
Item no. 3.0140.XX

HTO* plate



HTO plate 6.0/4.5 as medial

- Head holes: 3
- Shaft holes: 4 | 7
- Lengths: 63 | 93 mm
- Item no. left: 5.7252.3X
- Item no. right: 5.7242.3X




HTO spacer

- Heights: 5 – 15 mm
- Item no.: 5.7221.XX



HTO screw for spacer

- Size: M4 x 9 mm
- Item no.: 3.721.09



◆ Internal hexagon

HTO plate 4.5/6.0 as IH T SET

· Item no.: 19.357.00

*High Tibia Osteotomie



Cortical screw 4.5 as st
· Item no.: 3.545.XX



Cancellous screw 6.0 as st
· Item no.: 3.163.XX

T tibia head plate 6.0/4.5 as


T tibia head plate 6.0/4.5 as

- Head holes: 4
- Shaft holes: 2 – 16
- Lengths: 79 – 287 mm
- Item no. left: 5.777.XX*¹
- Item no. right: 5.778.XX*¹
- Set Instruments: 19.045.00I



T Tibia head plate 4.5/6.0 as IH SET

- Set no. steel: 19.404.00
- Set no. titanium: 19.408.00



 Internal hexagon



Standard Cortical screw 4.5 st
· Item no.: 3.152.XX*¹



Cortical screw 4.5 as st
· Item no.: 3.545.XX*¹



Cancellous screw 6.0 as st
· Item no.: 3.163.XX*¹

*¹ These items are also available in steel

L tibia head plate 6.0/4.5 as


L tibia head plate 6.0/4.5 as

- Head holes: 4
- Shaft holes: 3 – 16
- Lengths: 79 – 287 mm
- Item no. left: 5.773.XX*¹
- Item no. right: 5.774.XX*¹
- Set Instruments: 19.045.00I



L Tibia head plate 4.5/6.0 as IH SET

- Set no. steel: 19.402.00
- Set no. titanium: 19.406.00



 Internal hexagon



Standard Cortical screw 4.5 st
· Item no.: 3.152.XX*¹



Cortical screw 4.5 as st
· Item no.: 3.545.XX*¹




Cancellous screw 6.0 as st
· Item no.: 3.163.XX*¹

T & L Tibia head plate 4.5/6.0 as IH T SET

T & L Tibia head plate 4.5/6.0 as IH T SET

- Set no. titanium: 19.407.00



 Internal hexagon



Standard Cortical screw 4.5 st
· Item no.: 3.152.XX



Cortical screw 4.5 as st
· Item no.: 3.545.XX




Cancellous screw 6.0 as st
· Item no.: 3.163.XX

Medial tibia head plate 4.5/6.0 as



Medial tibia head plate 4.5/6.0 as

- Head holes: 3
- Shaft holes: 4 – 18
- Lengths: 77 – 301 mm
- Item no. left: 5.7192.XX
- Item no. right: 5.7182.XX
- Set Instruments: 19.045.00I



Med. Tibia head plate 4.5/6.0 as IH T SET

Set no.: 19.409.00

Internal hexagon



Standard Cortical screw 4.5 st

· Item no.: 3.152.XX



Cortical screw 4.5 as st

· Item no.: 3.545.XX



Cancellous screw 6.0 as st

· Item no.: 3.163.XX

Distal tibia plate distal 3.5 as/vas medial




Distal tibia plate distal 3.5 as medial

- Holes: 3 – 15
- Lengths: 76 – 268 mm
- Item no. left: 5.7332.XX
- Item no. right: 5.7322.XX
- Set Instruments: 19.035.00I



Distal tibia plate distal 3.5 vas medial

- Head holes: variable angle-stable
- Holes: 3 – 15
- Lengths: 76 – 268 mm
- Item no. right: 5.732.XX
- Item no. left: 5.733.XX
- Set Instruments: 19.035.00I



Distal Tibia plate 3.5 vas IH T SET

· Set no.: 19.491.00

Internal hexagon



Standard Cortical screw 3.5 st

· Item no.: 3.132.XX



Cortical screw 3.5 as st

· Item no.: 3.133.XX



Distal Tibia plate 3.5/4.0 vas IH T SET

· Set no.: 19.493.00

Internal hexagon



Standard Cortical screw 3.5 st

· Item no.: 3.132.XX



Cortical screw 4.0 as st

· Item no.: 3.148.XX

VaWiKo® distal fibula plate 3.5 holeØ5

VaWiKo® distal fibula plate 3.5 holeØ5

- Head holes: 4
- Shaft holes: 2 – 7
- Lengths: 83 – 148 mm
- Item no.: 5.07813.XX
- Set Instruments: 19.035.10I + 19.035.00I



Standard Cortical screw 3.5 st head 5

- TiOB®
- Item no.: 3.138.XX




Cortical screw 3.5 as st head 5

- TiOB®
- Item no.: 3.136.XX

VaWiKo® Distal Fibula plate 3.5 IH T SET

- Set no.: 19.474.00/2



Internal hexagon

Distal fibula plate 3.5 vas

Distal fibula plate 3.5 vas

- Holes: 7 | 9 | 11
- Lengths: 96 – 148 mm
- Item no.: 5.07803.XX
- Set Instruments: 19.035.00I



Standard Cortical screw 3.5 st

- Item no.: 3.132.XX



Cortical screw 3.5 as st

- Item no.: 3.133.XX

Distal Fibula plate 3.5 vas IH T SET

- Set no.: 19.474.00




Internal hexagon


Calcaneus

Calcaneus plate 3.5/4.0 vas with extension

- Holes: 10
- Size: 72 mm
- Item no. left: 5.895.00
- Item no. right: 5.894.00
- Set Instruments: 19.035.00I



Calcaneus




Calcaneus plate 3.5/4.0 vas with bridges

- Holes: 9
- Size: 65 mm
- Item no. left: 5.8793.09
- Item no. right: 5.8783.09
- Set Instruments: 19.035.00I



Calcaneus plate 3.5/4.0 vas small

- Holes: 8
- Can be ordered optionally to the set
- Item no. left: 5.8813.01
- Item no. right: 5.8803.01
- Set Instruments: 19.035.00I



Calcaneus plate 3.5/4.0 vas IH T SET

- Set no.: 19.326.00

Internal hexagon



Cortical screw 3.5 as st

- Item no.: 3.133.XX



Standard Cancellous screw 4.0

- Item no.: 3.145.XX



Cancellous screw 4.0 as st

- Item no.: 3.147.XX



Cannulated screw 4.0 partial thread st bt sd; TiOB®

- Item no.: 3.440.XX




Calcaneus plate 3.5/4.0 as large

- Holes: 10
- Size: 73 mm
- Item no. left: 5.879.01
- Item no. right: 5.878.01
- Set Instruments: 19.035.00I



Calcaneus plate 3.5/4.0 as small

- Holes: 8
- Size: 55 mm
- Item no. left: 5.881.01
- Item no. right: 5.880.01
- Set Instruments: 19.035.00I



Calcaneus plate 3.5/4.0 as IH T SET

- Set no.: 19.320.00

Internal hexagon



Standard Cortical screw 3.5 st

- Item no.: 3.132.XX



Cortical screw 3.5 as st

- Item no.: 3.133.XX



Standard Cancellous screw 4.0

- Item no.: 3.145.XX



Cancellous screw 4.0 as st

- Item no.: 3.147.XX

**VaWiKo® Lapidus Plate 3.5 holeØ5
plantar offset 2**

- Holes: 4
- Item no. right: 5.2923.02
- Item no. left: 5.2933.02


FOOT plate AD 3.5 holeØ5 vas

- Holes: 4
- Lengths: 23 | 28 mm
- Item no.: 5.2813.23
5.2813.28d


**FOOT plate AD interposition 3.5
holeØ5 vas**

- Holes: 5
- Lengths: 28 | 33 mm
- Item no.: 5.2813.28
5.2813.33


Multi-AD plate 3.5 holeØ5 vas

- 6 Loch with long bridge
- Length: 43 mm
- Item no.: 5.2823.43

Multi-AD plate 3.5 holeØ5 vas

- Holes: 6 – 12
- Lengths: 38 – 83 mm
- Item no.: 5.2823.XX


T-AD plate 3.5 holeØ5 vas

- Holes: 2/2
- Lengths: 33 | 38 mm
- Item no.: 5.2833.XX

T-AD plate 3.5 holeØ5 vas DC

- Holes: 2/2 | 2/3
- Lengths: 33–38 | 43 mm
- Item no.: 5.2833.XXK


L-AD plate 3.5 holeØ5 vas skew

- Holes: 2/2 | 2/3
- Lengths: 33 | 43 mm
- Item no. left: 5.2873.XX
- Item no. right: 5.2863.XX


LCF plate 3.5 holeØ5 vas

- Holes: 8
- Length: 74 mm
- Item no. right: 5.2883.74
- Item no. left: 5.2893.74



Set Hindfoot plates 3.5
holeØ5 vas can be
combined with Set Forefoot
plates 2.7 vas.


**Standard Cortical screw 3.5 st
head 5**

- TiOB®
- Item no.: 3.138.XX


Cortical screw 3.5 as st head 5

- TiOB®
- Item no.: 3.136.XX

Hindfoot plate 3.5 vas IH T SET

- Set no.: 19.524.00
- Set Instruments: 19.035.10I



Forefoot plates 2.7



VaWiKo® T plate 2.7
 · Holes: 2/7
 · Length: 73 mm
 · Item no.: 5.2623.27



Straight plate Mini 2.7 vas
 · Holes: 2
 · Length: 18 mm
 · Item no.: 5.2843.18



Straight plate Mini 2.7 vas DC
 · Holes: 4
 · Length: 37 mm
 · Item no.: 5.2843.37K



Straight plate Mini 2.7 as DC2
 · Holes: 5 | 6
 · Lengths: 46 | 55 mm
 · Item no.: 5.2753.46K
 5.2753.55K



L-AD plate 2.7 vas skew
 · Holes: 2/2 – 2/3
 · Lengths: 28 | 37 mm
 · Item no. left: 5.2773.XX
 · Item no. right: 5.2763.XX



T-AD plate 2.7 vas
 · Holes: 2/2
 · Length: 28 mm
 · Item no.: 5.2783.28



T-AD plate 2.7 vas DC
 · Holes: 2/3
 · Length: 37 mm
 · Item no.: 5.2783.37K



MTP arthrodesis plate 2.7 vas
 · Holes: 6
 · Length: 41 mm
 · Item no. left: 5.2793.41L
 · Item no. right: 5.2793.41R



· Holes: 6
 · Length: 46 mm
 · Item no. left: 5.2793.46L
 · Item no. right: 5.2793.46R



· Holes: 8
 · Length: 51 mm
 · Item no. left: 5.2793.51L
 · Item no. right: 5.2793.51R



Double bridge plate 3.5 hole Ø5 vas bridge thin
 · bridge thickness: 1 mm
 · Holes: 4
 · Size: 18 x 20 mm
 · Item no.: 5.2803.20d

for revision

Double bridge plate 3.5 hole Ø5 vas

- Holes: 4
- Size: 16 x 18 mm
- Item no.: 5.2803.18



- Holes: 4
- Size: 18 x 20 mm
- Item no.: 5.2803.20



- Holes: 4
- Size: 16 x 22 mm
- Item no.: 5.2803.22



- Holes: 4
- Size: 18 x 24 mm
- Item no.: 5.2803.24



- Holes: 8
- Size: 16 x 37 mm
- Item no.: 5.2803.37



Cortical screw 2.7
flat head st
· Item no.: 3.128.XX



Cortical screw 2.7 as st
· Item no.: 3.125.XX

Forefoot plate 2.7 vas IH T SET

- Set no.: 19.523.00
- Set Instruments: 19.027.00I



Internal hexagon



Cortical screw 2.7 flat head
SL10 st
· Item no.: 3.1281.XX



Cortical screw 2.7 as SL10 st
· Item no.: 3.1251.XX

Forefoot plate 2.7 vas SL T SET

- Set no.: 19.5231.00
- Set Instruments: 19.027.10I



Six lobe

Step-AD plate 2.7 vas

- Holes: 5
- Steps: 2 | 3 | 4 | 5 mm
- Colours: gold | pink | light blue | dark violett
- Item no.: 5.2913.22/X
- Set Instruments: 19.027.00I (IH)
- Set Instruments: 19.027.10I (SL)



Foot plate Y 2.7 vas




Foot plate Y 2.7 vas

- Holes: 2/2
- Size: 34 x 18 mm
- Item no.: 5.2903.34



Spacer foot plate Y 2.7

- Size: 8.0 x 3.5 mm
- Heights: 2 | 3 | 4 | 5 | 6 mm
- Item no.: 5.290.XX



Forefoot plate 2.7 Y vas IH T SET

- Set no.: 19.523.20
- Set Instruments: 19.027.00I (IH)
- Set Instruments: 19.027.10I (SL)


 Internal hexagon

Forefoot surgery



Bone clamp 8x10

- intended
- Width: 8 | 10 mm
- Item no.: 7.399.XX



Forefoot surgery T SET

- Set no.: 19.529.00

FRAKO 3.0/4.0 DrillFix 2.0



Drill Fix screw 2.0

- Item no.: 3.127.XX
- Set Instruments: 19.335.00I



Fracture compression screw 3.0/4.0

- Item no.: 3.123.XX
- Set Instruments: 19.335.00I



FRAKO 3.0/4.0 DrillFix 2.0 SET

- Set no.: 19.335.00

for Scarf-, Chevron-,
Weil-osteotomies and
fractures

Königsee Implantate Cervical Variable Anterior Plate System

KI CVA plate 4.5/5.2

- Surface: TiOB®
- Holes: 4 – 14
- Lengths: 24 – 93 mm
- Item no.: 5.555.XXXX




Spine screw 4.5
· Item no.: 3.031.XX



Spine screw 5.2
· Item no.: 3.033.XX

KI CVA plate 4.5/5.2 SL T SET

Set no.: 19.651.00



★ Six lobe

Fracture compression screw 3.0/4.0 st sd can

Fracture compression screw 3.0/4.0 st sd can

- Lengths: 32 – 44 mm
- Inner Ø: 1.3 mm
- Item no.: 3.123.XX
- Set Instruments: 19.380.10I



DENS Axis screwing IH T SET

Set no.: 19.380.10



⬡ Internal hexagon

LC titanium strip


LC titanium strip

- with clamp
- Size: 8.0 x 0.4 mm
- Length: 260 mm
- Item no.: 7.440.01



LC titanium strip T SET

· Set no.: 19.519.00



CPL plate 4.5 as




CPL plate 4.5 as narrow

- Notches: 6 – 12
- Lengths: 124 – 214 mm
- Width: 12 mm
- Item no.: 4.2341.XX
- Set Instruments: 19.045.00I


CPL plate 4.5 as broad

- Notches: 8 – 14
- Lengths: 174 – 270 mm
- Width: 16 mm
- Item no.: 4.2461.XX
- Set Instruments: 19.045.00I




Cable osteosynthesis 1.0 S SET

- Cables: Ø 1.0
- Set no.: 19.518.00



Cable osteosynthesis 1.5 S SET

- Cables: Ø 1.5
- Set no.: 19.518.00/1



CPL plate 4.5 as IH S SET

- Set no.: 19.516.00

Internal hexagon



Cortical screw 4.5 as st

- Item no.: 2.545.XX



Standard Cortical screw 4.5 st

- Item no.: 2.152.XX

Distal femur plate 4.5-6.0 as notch



Distal femur plate 4.5-6.0 as notch

- Holes: 12
- Length: 246 mm
- Width: 18 mm
- Item no. left: T.292.12L
- Item no. right: T.292.12R
- Combinable: 19.518.00/1
19.600.00
- Set Instruments: 19.045.00I

Distal femur plate 4.5-6.0 as notch

- Holes: 16
- Length: 309 mm
- Width: 18 mm
- Item no. left: T.292.16L
- Item no. right: T.292.16R
- Combinable: 19.518.00/1
19.600.00
- Set Instruments: 19.045.00I

for rope cerclages Ø 1.5

VaWiKo® EPI plate 4.0

for Hemi-/Epiphyseodesis

VaWiKo® EPI plate hole distance 12;
TiOB®

· Item no. 5.6983.2120

without spikes

hole spacing 12

VaWiKo® EPI plate hole distance 12
with spikes; TiOB®

· Item no. 5.6983.2122



with spikes



hole spacing 16

VaWiKo® EPI plate hole distance 16
without spikes; TiOB®

· Item no. 5.6983.2160

VaWiKo® EPI plate hole distance 16
with spikes; TiOB®

· Item no. 5.6983.2162



Cortical screw 4.0 as st

· Item no.: 3.148.XX



Cannulated screw 4.0 st sb

· TiOB®

· Item no.: 3.445.XX

VaWiKo® EPI plate 4.0 vas IH T SET

Set no.: 19.563.00

Set Instruments: 19.563.00I



Internal hexagon

DET screw

Dynamic Epiphyseal
Teleskopic screw for
Epiphysiolysis captis
femoris

DET screw 7.0 st sb

· Lengths: 40 – 80 mm

· Inner Ø: 2.2 mm

· Item no.: 6.285.XX



DET screw sleeve 7.5


· Lengths: 40 – 80 mm

· Inner Ø: 5.6 mm

· Item no. : 6.286.XX

DET 7.0 IH S SET

· Set no.: 19.501.00



Internal hexagon

Gliding nail




Gliding nail 300

- Length: 300 mm
- Ø: 1.2 – 5.0
- Item no. titanium: 7.282.XX
- Item no. steel: 6.282.XX

Gliding nail 450

- Length: 450 mm
- Ø: 1.2 – 5.0
- Item no. titanium: 7.281.XX
- Item no. steel: 6.281.XX



Gliding nail 300 SET

- Set no.: titanium: 19.290.00
- Set no.: steel: 19.292.00

Gliding nail 450 SET

- Set no. titanium: 19.290.10
- Set no. steel: 19.292.10

for shaft fractures

Lengths 320 and 350 mm ebenfalls erhältlich, aber nicht auf Set.

Telescope Nail Bailey-Dubow



Telescope Nail Bailey-Dubow 3.2

- Lengths: 70 – 250 mm
- Item no.: 7.288.00/XXX

Telescope Nail Bailey-Dubow 4.0

- Lengths: 120 – 250 mm
- Item no.: 7.289.00/XXX

Telescope Nail Bailey-Dubow 4.8

- Lengths: 70 – 300 mm
- Item no.: 7.290.00/XXX

Telescope Nail Bailey-Dubow 5.5

- Lengths: 130 – 300 mm
- Item no.: 7.291.00/XXX

Telescope Nail Bailey-Dubow 6.4

- Lengths: 250 – 350 mm
- Item no.: 7.292.00/XXX

for Osteogenesis imperfecta



Drilling point Telescope Nail Bailey-Dubow

- Item no.: 7.2XX.01



Cancellous screw 4.0 Telescope Nail Bailey-Dubow

- Item no.: 7.2XX.02

Kid plate femur corrective prox. 3.5/6.0



Kid plate femur corrective prox. 3.5/6.0 AD4

- Holes: 3 | 4
- Lengths: 70 | 81 mm
- Item no.: 4.6402.04XX



Kid plate femur corrective prox. 3.5/6.0 AD6

- Holes: 3 | 4
- Lengths: 73 | 84 mm
- Item no.: 4.6402.06XX

„for varization and / or derotation osteotomy at the pediatric coxal end of the femur“

Kid plate femur corrective prox. 3.5/6.0

Kid plate femur corrective prox. 3.5/6.0 AD8

- Holes: 3 | 4
- Lengths: 77 | 88 mm
- Item no.: 4.6402.08XX



Standard Cortical screw 3.5 st

- Item no.: 2.132.XX



Cortical screw 3.5 as st

- Item no.: 2.133.XX




Cannulated screw 6.0 as st sb

- Item no.: 2.540.XX

WiKo® DVO children 3.5/6.0 as IH S SET

- Set no.: 19.561.10
- Set Instruments: 19.561.10I



Compression wire Ø1.2/1.8



Compression wire Ø1.2/1.8

- Thread distance: 5 | 10 | 15 mm
- Length: 150 mm
- Item no.: 2.919.18/XX



TIS screw 7.5 st sd

TIS screw 7.5 st sd

- Surface: TiOB®
- Lengths: 50 – 200 mm
- Item no.: 3.530.XX




TIS screw 7.5 thread 32 st bt sd

- Surface: TiOB®
- Lengths: 50 – 200 mm
- Item no.: 3.531.XX



TIS screws 7.5 st IH T SET

- Set no.: 19.681.00
- Set Instruments: 19.075.00I



Internal hexagon

CHILDREN

WIRES

CANNULATED

TIS screw 7.5 st sd




Internal hexagon

TIS screws 7.5 st IH T upgrade
SCREWS

- Set no.: 19.681.00/1
- Set Instruments: 19.075.00I

SPASY™ Bone Cement Filler

- Item no. D.71309.00I



Cannulated Screws



Internal hexagon

Cann. screws 2.7 st sd IH T SET

- Set no.: 19.170.00
- Set Instruments: 19.170.00I

Cannulated screw 2.7
partial thread st sb; TiOB®

- Item no. 3.448.XX



Cannulated screw 2.7 st sb; TiOB®

- Item no. 3.450.XX



Internal hexagon

Cann. screws 3.5 st sd IH T SET

- Set no.: 19.184.00
- Set Instruments: 19.184.00I

Cannulated screw 3.5
partial thread st bt sd

- Item no. 3.504.XX^{*1}



Cannulated screw 3.5 st sb

- Item no. 3.604.XX^{*1}



Internal hexagon

Cann. screws 4.0 st sd IH T SET

- Set no.: 19.190.00
- Set Instruments: 19.190.00I

Cannulated screw 4.0
partial thread st bt sd; TiOB®

- Item no. 3.440.XX



Cannulated screw 4.0 st sb; TiOB®

- Item no. 3.445.XX

^{*1} These items are also available in steel




**Cannulated screw 4.5
partial thread st bt sd**
· Item no. 3.434.XX^{*1}



Cannulated screw 4.5 st sb
· Item no. 3.439.XX

Cann. screws 4.5 st sd IH T SET

- Set no.: 19.203.00
- Set Instruments: 19.203.00I



Internal hexagon




**Cannulated screw 6.5
thread 16 st bt sd; TiOB®**
· Item no. 3.520.XX^{*1}



**Cannulated screw 6.5
thread 32 st bt sd; TiOB®**
· Item no. 3.523.XX^{*1}

Cann. screws 6.5 st sd IH T SET

- Set no.: 19.210.00
- Set Instruments: 19.210.00I



Internal hexagon

available separately:

Cann. screws 7.0 st sd IH T SET

- Surface: TiOB®
- Item no. 3.525.XX^{*1}



**Cannulated screw 7.0
thread 16 st bt sd**
· Item no. 3.414.XX^{*1}



**Cannulated screw 7.0
thread 32 st bt sd**
· Item no. 3.424.XX^{*1}

DB Screws 7.0 st sb IH T SET

- Set no.: 19.223.00
- Set Instruments: 19.223.00I



Internal hexagon

available separately:

Cannulated screw 7.0 st sb

- Item no. 3.411.XX



Further cannulated screws that are not self-drilling are available on request.

*1 These items are also available in steel

SF screws as 3.5/4.0



Internal hexagon

SF screws as 3.5/4.0 IH Cancellous 4.0 T SCREWS

· Set no.: 19.533.00



Cortical screw 3.5 st
· Item no.: 3.132.XX



Cortical screw 3.5 as st
· Item no.: 3.133.XX



Cancellous screw 4.0
partial thread bt
· Item no.: 3.140.XX



Cancellous screw 4.0
· Item no.: 3.145.XX



Cancellous screw 4.0 as st
· Item no.: 3.147.XX



Internal hexagon

SF screws as 3.5/4.0 IH Cortical 4.0 T SCREWS

· Set no.: 19.533.10



Cortical screw 3.5 st
· Item no.: 3.132.XX



Cortical screw 3.5 as st
· Item no.: 3.133.XX



Cancellous screw 4.0
partial thread bt
· Item no.: 3.140.XX




Cancellous screw 4.0
· Item no.: 3.145.XX



Cortical screw 4.0 as st
· Item no.: 3.148.XX

LF screws 4.5/6.0



Internal hexagon

LF screws 4.5/6.0 IH T SET

· Set no. titanium: 19.168.40


· Set no. steel: 19.168.20




Cortical screw 4.5 st
· Item no.: 3.152.XX



Cancellous screw 6.5
thread 16
· Item no.: 3.170.XX



Cancellous screw 6.5
thread 32
· Item no.: 3.180.XX



Cancellous screw 6.5
· Item no.: 3.190.XX

Rep Ret MF SET


· Set no.: 19.426.00

Rep Ret SF SET


· Set no.: 19.424.00

Rep Ret LF SET

· Set no.: 19.425.00

**IMPLANTEX® – Metal removal****IMPLANTEX® INSTRUMENTS**

· Set no.: 19.511.00

**Mini fragment Standard****Standard Cortical screw 1.5 st**

· Item no.: 3.100.XX

**Standard Cortical screw 2.0 st**

· Item no.: 3.110.XX


**Standard Cortical screw 2.7 st**

· Item no.: 3.120.XX

MF plate IH T SET

· Set no. titanium: 19.141.00

· Set no. steel: 19.140.00

**Straight plate Mini 1.5**

· Item no. titanium: 5.404.XX

· Item no. steel: 4.404.XX

**L plate Mini 2.0 skew**

· Item no. titanium: 5.92X.XX

· Item no. steel: 4.92X.XX

**Adaption plate Mini 2.0**

· Item no. titanium: 5.307.XX

· Item no. steel: 4.307.XX

**Straight plate Mini 2.0**

· Item no. titanium: 5.900.XX

· Item no. steel: 4.900.XX

**T plate Mini 2.0 long hole**


· Item no. titanium: 5.930.XX

· Item no. steel: 4.930.XX

**Condylar plate Mini hand 1.5 pin 21**


· Item no. titanium: 5.405.XX

Small fragment Standard




SF Plate IH SET

- Set no. titanium: 19.467.00
- Set no. steel: 19.466.00
- Set Instruments: 19.035.00I



Internal hexagon



Standard Cortical screw 3.5 st

- Item no.: 3.132.XX



Standard Cancellous screw 4.0

- Item no.: 3.145.XX



Standard Cancellous screw 4.0 partial thread bt

- Item no.: 3.140.XX



Low contact plate 3.5/4.0

- Item no. titanium: 5.220.XX
- Item no. steel: 4.220.XX



Third tubular plate 3.5 bulge

- Item no. titanium: 5.808.XX
- Item no. steel: 4.808.XX



T plate 3.5

- Item no. titanium: 5.817.XX
- Item no. steel: 4.817.XX



Reco plate 3.5 strong

- Item no. titanium: 5.313.XX
- Item no. steel: 4.313.XX



T plate 3.5 skew long hole

- Item no. titanium: 5.812.XX
- Item no. steel: 4.812.XX



Calcaneus Plate 3.5 small

- Item no. titanium: 5.86X.01
- Item no. steel: 4.86X.01

Straight plate 3.5 as



Internal hexagon

Straight plate 3.5 as IH T SET

- Set no.: 19.473.00
- Set Instruments: 19.035.00I



Standard Cortical screw 3.5 st

- Item no.: 3.132.XX



Cortical screw 3.5 as st

- Item no.: 3.133.XX



Third tubular plate 3.5 as

- Item no.: 5.809.XX
- Set Instruments: 19.035.00I




Reco plate 3.5 vas strong

- Item no.: 5.3133.XX
- Set Instruments: 19.035.00I

Straight plate 3.5 T PLATES

Straight plate 3.5 T PLATES

- Set no.: 19.535.00
- Set Instruments: 19.035.001



 Internal hexagon

Third tubular plate 3.5 as

- Item no.: 5.809.XX



Reco plate 3.5 vas strong

- Item no.: 5.3133.XX



Straight plate 3.5 vas

- Item no.: 5.823.XX



Distal fibula plate 3.5 vas


- Item no.: 5.07803.XX



Kombi 2.7/3.5 IH T SCREWS INSTRUMENTS

Kombi 2.7/3.5 IH T SCREWS INSTRUMENTS

- Set no.: 19.671.00
- Set Instruments: 19.035.001



 Internal hexagon



Cortical screw 2.7
flat head st

- Item no.: 3.128.XX



Cortical screw 2.7 as st

- Item no.: 3.125.XX



Cortical screw 3.5 st head 5


- Item no.: 3.138.XX



Cortical screw 3.5 as st head 5

- Item no.: 3.136.XX

Kombi 3.5/4.5 IH T SCREWS INSTRUMENTS



Internal hexagon

Kombi 3.5/4.5 IH T SCREWS INSTRUMENTS

· Set no.: 19.675.00



Cortical screw 3.5 st
· Item no.: 3.132.XX



Cortical screw 3.5 as st
· Item no.: 3.133.XX



Cortical screw 4.5 st
· Item no.: 3.152.XX



Cortical screw 4.5 as st
· Item no.: 3.545.XX

LF PLATES



Internal hexagon


LF plate T PLATES

· Set no. titanium: 19.162.00
· Set no. steel: 19.161.00
· Set Instruments: 19.045.00I



Straight plate 4.5/6.5 DC narrow

· Item no. titanium: 5.530.XX
· Item no. steel: 4.530.XX



Straight plate 4.5/6.5 DC broad

· Item no. titanium: 5.540.XX
· Item no. steel: 4.540.XX

Low Contact plates 4.5/6.0 as




Low contact plate 4.5/6.0 as narrow

· 2 Kombi-Kompressionslöchern
· Holes: 2 – 16
· Lengths: 106 – 286 mm
· Item no. titanium: 5.231.XX
· Item no. steel: 4.231.XX



Low contact plate 4.5/6.0 as broad

· 2 Kombi-Kompressionslöchern
· Holes: 6 – 20
· Lengths: 106 – 358 mm
· Item no. titanium: 5.241.XX
· Item no. steel: 4.241.XX



Internal hexagon

Low Contact plate 4.5/6.0 as IH T SET

· Set no. titanium: 19.514.00
· Set no. steel: 19.512.00
· Set Instruments: 19.045.00I



Standard Cortical screw 4.5 st
· Item no.: 3.152.XX



Cortical screw 4.5 as st
· Item no.: 3.545.XX

INSTRUMENTS Mini fragment

MF 2.4 tools SL INSTRUMENTS


· Set no.: 19.024.00I



★ Six lobe

MF 2.7 tools IH INSTRUMENTS


· Set no.: 19.027.00I



◆ Internal hexagon

MF 2.7 tools SL INSTRUMENTS

· Set no.: 19.027.10I



★ Six lobe

INSTRUMENTS Small fragment

SF 3.5 tools IH INSTRUMENTS


· Set no.: 19.035.00I



◆ Internal hexagon

SF 3.5 tools IH head 5
INSTRUMENTS

· Set no.: 19.035.10I




◆ Internal hexagon

INSTRUMENTS Large fragment

LF tools IH INSTRUMENTS

· Set no.: 19.045.00I



◆ Internal hexagon

TIS screw 7.5 st IH INSTRUMENTS

· Set no.: 19.075.00I



◆ Internal hexagon



Königsee Implantate GmbH

OT Aschau | Am Sand 4
07426 Allendorf/Thuringia
Germany

Fon: + 49 36738 498 - 560

Fax: + 49 36738 498 - 569

export@koenigsee-implantate.de

www.koenigsee-implantate.com

PRODUCT OVERVIEW

Implants