




Tel: 1-613-258-3344
Email: inquires@tlcmedicalsupply.ca & mario@tlcmedicalsupply.ca




Baffin Automatic

The Baffin Automatic multifunctional device is designed for uprighting, stabilization in sitting and lying positions, and patient transport. The product is designed for children and adolescents.

The Baffin Automatic allows the patient to build a correct sitting posture pattern. The patient sitting in the device can participate in educational activities, occupational therapy and play. Baffin promotes patient activation and the building of cognitive functions. The device is controlled by electric actuators. Changing the position between sitting, lying and standing is done using a remote control.



COLOURS



IMPORTANT FEATURES



The patented "spine" system allows anatomical or corrective adjustment of the seat back.



Adjustable hip supports and armrests.




Table - the large working space of the table tray provides convenience during therapeutic activities.



System of supports, vest, stabilizing belts.



Split foot platforms. It is possible to set a preset angle of flexion in the knee joint and position the foot in dorsiflexion in the sole, as well as the possibility of supination and pronation.



Adjustable, flexible headrest that adapts to the shape of the head.

TECHNICAL DATA

	SIZE S	SIZE M	SIZE L
Back length	44 cm / 17"	52 cm / 20.4"	60 cm / 23.6"
Seat depth	18÷28 cm / 7"-11"	21÷31 cm / 8.2"- 12.2"	30÷39 cm / 11.8"- 15.3"
Seat width	18÷29 cm / 7"-11.4"	22÷33 cm / 8.6"- 13"	22÷33 cm / 8.6"- 13"
Backrest angle	0°÷90°	0°÷90°	0°÷90°
Range of verticalization	0°÷86°	0°÷86°	0°÷86°
Maximum Patient weight	40 kg / 88 lb	60 kg / 132 lb	60 kg / 132 lb