





Tel: 1-613-258-3344
Email: inquires@tlcmedicalsupply.ca &
mario@tlcmedicalsupply.ca



Baffin TRIO 4W

The Baffin TRIO 4W multifunctional device is designed for upright standing, stabilization in sitting and lying positions, and patient transport. The product is tailored for adolescents and adults.

A patient equipped with the device can participate in educational activities, occupational therapy and play while maintaining a proper sitting position. The Baffin TRIO is an excellent tool for physiotherapists to facilitate the start of the uprighting process even in the most demanding Patients with knee and hip contractures.



COLOURS



IMPORTANT FEATURES



The patented "spine" system allows anatomical or corrective adjustment of the seat back.



Systems of removable bolsters, hip caps, and footrests to facilitate transfer of the patient to the device



The large table-now as large as 2200 cm²-will make playing and learning even easier.



Smooth change of position. Great possibilities for passive posture correction in sitting and lying positions.



Stable, subtle frame with 4 swivel wheels.



Adjustable, flexible headrest that adapts to the shape of the head.

TECHNICAL DATA

Back length

SIZE I

61 cm / 24"

SIZE II

61 or 77 cm / 24" or 30.3"

Seat depth

36÷50 cm / 14"-19.5"

36÷50 cm lub 50÷64 cm / 14"-19.5" or 19.5"-25.2"

Seat width

27÷40 cm / 10.6"-15.7"

34÷45 cm

Backrest angle

80°÷115°

80°÷115°

Range of verticalization

0°÷85°

0°÷85°

Maximum Patient weight

90 kg / 198 lb

90 kg / 198 lb