

Price List valid from 01.04.2024

**BAFFIN TRIO 4W**

LIW's order number:	Date:
Dealer name:	Customer's order number:
Invoice information:	
Street:	
Postal Code:	
Shipping information:	
Street:	
Postal Code:	
Tel.	
Ordered For:	
Quantity:	
<b>BAFFIN TRIO 4W SIZE 1 BASE CONFIGURATION (STD): size 1 3 285 CAD</b>	
<b>BAFFIN TRIO 4W SIZE 2 BASE CONFIGURATION (STD): size 2 3 285 CAD</b>	



**COLOUR**



Blue – standard colour for Size I



Grey – standard colour for Size II

**COLOURS available on request.  
Lead time will vary.**

BASIC PARAMETERS	SIZE	Baffin TRIO 4W SIZE I	Baffin TRIO 4W SIZE II
Seat depth	A	36 -50 cm / 14"-19.5"	50-64 cm / 19.5"-25.2" size I optional/available on request
Back length	B	61 cm / 24"	77 cm / 30.3" size I optional/available on request
Seat width	C	27-40 cm / 10.6" – 15.7"	34-45 cm / 13.4"-17.7"
User's height	D	135-170 cm / 53"-67"	160-190 cm / 64"-74.8"
Maximum User's weight		90 kg / 198 lb	90 kg / 198 lb



Blue




Red



Grey

DIMENSIONS		
A		cm / inches
B		cm / inches
C		cm / inches
D		cm / inches
E		cm / inches
F		cm / inches
G		cm / inches
Weight		kg / lb
Additional notes		



Code no.	Pict. no.	Name	Comments	Qty.	Price (CAD)
BAFT-1 , BAFT-2 BAFW-1 , BAFW-2	-	Multifunctional electrical seat / standing /lying frame consisting of a seat including upholstery		1	STD
BAFT-L-04	6	Lateral supports with elbow blocks – size L F = min 36,5 cm / 14.3” (from seat to top edge of the lateral’s pad) <i>Other sizes available – see pictures below</i>	2 pcs.	1	STD
BAFT-C-01	8	Armrests	2 pcs.	1	STD
BAFT-C-02	9	Hip supports	2 pcs.	1	STD
BAFT-B-01-L BAFT-B-01-XL BAFT-B-01-XXL	10	Pelvic belts Size L      XL      XXL		1	STD
n/a	11	Legs abduction belts	2 pcs.	1	STD
BAFT-K-01-1 BAFT-K-01-2	13	Kneepads 14 cm / 5.5 ” wide For TRIO I      For TRIO II	2 pcs.	1	STD
BAFT-K-02-1 BAFT-K-02-2	13	Kneepads 17 cm / 6.7 ” wide For TRIO I      For TRIO II			
BAFT-F-01-L BAFT-F-01-XL BAFT-F-01-XXL	14	3D rotational foot platforms Size L      XL      XXL	2 pcs.	1	STD
BAFT-F-04-1 BAFT-F-04-2 BAFT-F-04-3	16	Foot platforms height adjustment Range 1      Range 2      Range 3	2 pcs.	1	STD
BAFT-F-02	15	Foot straps	4 pcs.	1	STD
BAFT-T-01	19	Tray (acrylic)		1	STD

**HEADRESTS – Minimum one option below required**

BAFT-H-01-1	1	Butterfly headrest with multi adjustable headrest hardware		1	680
BAFT-H-02-1	2	Headrest enlarged on both sides	upcharge	1	120
BAFT-H-03-1	3	Headrest right side enlarged	upcharge	1	100
BAFT-H-04-1	4	Headrest left side enlarged	upcharge	1	100

**LATERAL SUPPORTS – one option required**

BAFT-L-91	-	Additional foam for lateral supports	upcharge	1	100
BAFT-L-92	-	Additional foam for hip supports	upcharge	1	100

**TRUNK SUPPORTS – required indicated by\***

BAFT-B-02*	18	Second Spine Backrest for TRIO I		1	2 000
BAFT-B-03*	18	Second Spine Backrest for TRIO II		1	2 500
BAFT-V-01-L BAFT-V-01-XL BAFT-V-01-XXL	7	Vest Size L      XL      XXL		1	150


**FOOT SUPPORT ACCESSORIES - optional**

BAFT-F-03-L BAFT-F-03-XL BAFT-F-03-XXL	17	3D rotational food platforms <b>with heel supports</b> Size L      XL      XXL	Upcharge option	1	100
--	----	---	-----------------	---	-----

**ADDITIONAL ACCESSORIES – required indicated by\***

BAFT-A-06 *	-	Electric actuator for TRIO I		1	2 335
BAFT-A-07 *	-	Electric actuator for TRIO II		1	2 835
BAFT-T-02	20	Tray padding		1	120
BAFT-A-03 *	-	Battery and charger	1 set.	1	1 500
BAFT-A-04	21	Respiratory shelf		1	200
BAFT-A-02	12	Additional seat filler pad		1	150
BAFT-A-05-01 BAFT-A-05-02	-	Extra set of standard upholstery For TRIO I      TRIO II		1	400

**PICTURES**



**1** Butterfly headrest

A = 40 cm / 15.7"  
B = 25 cm / 9.8"  
H = 15 cm / 5.9"

**2**

Headrest enlarged on **both sides**

Up horn + 5 cm / + 2"  
Bottom horn + 7 cm / + 2.7"

**3**


Headrest **right side** enlarged

Up horn + 5 cm / + 2"  
Bottom horn + 7 cm / + 2.7"

**4**

Headrest **left side** enlarged

Up horn + 5 cm / + 2"  
Bottom horn + 7 cm / + 2.7"



**5** Lateral supports

Spacing: 27 - 40 cm / 10.6" - 15.7"

**Size L**

Pad: 16 cm x 15 cm / 6.2" x 5.9"  
F = min 32 cm / 12.6"

**Size XL**

Pad: 22 cm x 12 cm / 8.6" x 4.7"  
F = min 34 cm / 13.3"

**6**

Lateral supports with elbow blocks

Spacing: 27 - 40 cm / 10.6" - 15.7"

**Size L**

Pad: 16 cm x 15 cm / 6.2" x 5.9"  
F = min 36.5 cm / 14.3"

**Size XL**

Pad: 22 cm x 12 cm / 8.6" x 4.7"  
F = min 39 cm / 15.3"

**7**

**Vest L**


A - 39 cm / 15.3"  
B - 25 cm / 9.8"  
C - 20 cm / 7.8"  
D - 8.5 cm / 3.3"  
E - 29 cm / 11.4"

**Vest XL**

A - 46 cm / 18.1"  
B - 32 cm / 12.5"  
C - 21 cm / 8.3"  
D - 17 cm / 6.7"  
E - 26 cm / 10.2"

**Vest XXL**

A - 46 cm / 18.1"  
B - 32 cm / 12.5"  
C - 21 cm / 8.3"  
D - 17 cm / 6.7"  
E - 26 cm / 10.2"



**8** Armrests

**9** Hip supports

Spacing: 27 - 40 cm / 10.6" - 15.7"  
Pad dimension: 23 cm x 18 cm / 9.0" x 7.0"

**10**

**Pelvic belts**

**L:**  
A - 45 cm / 17.7"  
B - 10 cm / 3.9"

**Pelvic belts XL:**

A - 52 cm / 20.5"  
B - 12 cm / 4.7"

**Pelvic belts XXL:**

A - 59 cm / 23.2"  
B - 14 cm / 5.5"

**11** Legs abduction belts



**12**

**Additional seat filler pad**

**13**


**Kneepads:** 14 cm wide / 5.5" / 17 cm wide / 5.5"

**14**

**3D rotational foot platforms size L**  
Length: 24 cm / 6.4"  
Width: 12 cm / 4.7"  
Height: 5.5 cm / 2.1"


**3D rotational foot platforms size XL**  
Length: 28 cm / 11.0"  
Width: 12.5 cm / 4.9"  
Height: 5.5 cm / 2.1"

**3D rotational foot platforms size XXL**  
Length: 32 cm / 12.6"  
Width: 16 cm / 6.2"  
Height: 5.5 cm / 2.1"



**15**


**Feet straps** (set of 4 pcs.)



**16**

**Feet platforms height adjustment**


Dimension <b>G</b>	3D Feet platforms ( <b>14</b> )
Platform height range 1	25,5-38 cm / 10.0" – 14.9"
Platform height range 2	31,5-43,5 cm / 12.4" – 17.1"
Platform height range 3	41,5-51,5 cm / 16.3"- 20.2"



**17**

**3D rotational foot platforms with heel supports**

Heel supports height: 12 cm / 4,7"



**18**

**Second Spine Backrest**




**19**

**Tray**



**20**

**Tray padding**



**21**

**Respiratory shelf**