


Price List valid from 01.04.2024


BAFFIN AUTOMATIC

The multifunctional standing frame Baffin Automatic is designed to set the patient in multiple sitting positions or the standing position. It is available in 3 sizes: S, M, L for users between 85 cm 33" and 130 cm 53" tall.




LIW's order number:	Date:
Dealer name:	Customer's order number:
Invoice information:	
Street:	
Postal Code:	
Shipping information:	
Street:	
Postal Code:	
Tel.	
Ordered For:	
Quantity:	
PRICE FOR STANDARD CONFIGURATION (STD): Size S CAD 3 285	
PRICE FOR STANDARD CONFIGURATION (STD): Size M CAD 3 285	
PRICE FOR STANDARD CONFIGURATION (STD): Size L CAD 3 285	

COLOUR

 Red – standard colour for Size M

 Blue – standard colour for Size L


COLOURS available on request.
Lead time will vary.

 Grey  Red  Blue




BASIC PARAMETERS (cm)	SIZE	Baffin Automatic SIZE S	Baffin Automatic SIZE M	Baffin Automatic SIZE L
Seat depth	A	18-28 / 7"-11"	21-31 / 8.2"- 12.2"	30-39 / 11.8"- 15.3"
Back length	B	44/ 17"	52 / 20.4"	60/ 23.6"
Seat width	C	18-29 / 18"-11.4"	22-33 / 8.6"- 13"	22-33/ 8.6"- 13"
User's height	D	85-105 / 33.4"-41.3"	100-120 / 39.3"-47.2"	115-135 / 45.2"- 53.1"
Maximum User's weight		40 kg / 88 lb	60 kg / 132 lb	60 kg / 132 lb

DIMENSIONS		
A		cm / inches
B		cm / inches
C		cm / inches
D		cm / inches
E		cm / inches
F		cm / inches
G		cm / inches
Weight		kg / lb
Additional notes		




Code no.	Pict. no.	Name	Comments	Qty.	Price (CAD)
BAF-S BAF-M BAF-L	-	Multifunctional electrical seat / standing /lying frame consisting of a seat including upholstery		1	STD
BAF-L-02-S BAF-L-02-M BAF-L-02-L	6	Lateral supports with Elbow blocks Size S M L	2 pcs.	1	STD
BAFARS-C-01	-	Armrests (model RS)	2 pcs.	1	STD
BAF-C-02	-	Hip supports	2 pcs.	1	STD
BAF-B-01-S BAF-B-01-M BAF-B-01-L	9	Pelvic belts Size S M L		1	STD
BAF-K-03	10	Legs abduction belts	2 pcs.	1	STD
BAF-K-01	11	Kneepads	2 pcs.	1	STD
BAF-F-01-S BAF-F-01-M BAF-F-01-L	12	3D Foot platforms Size S M L	2 pcs.	1	STD
BAF-F-03-S BAF-F-03-M BAF-F-03-L	14	Foot platforms height adjustment Range 1 2 3	2 pcs.	1	STD
BAF-F-02	13	Foot straps	4 pcs.	1	STD
BAF-T-01	16	Tray (acrylic)		1	STD
HEADRESTS - optional					
BAF-H-01-1	1	Butterfly headrest with multi adjustable headrest hardware		1	680
BAF-H-02-1	2	Headrest enlarged on both sides	upcharge	1	120
BAF-H-03-1	3	Headrest right side enlarged	upcharge	1	100
BAF-H-04-1	4	Headrest left side enlarged	upcharge	1	100
LATERAL SUPPORTS - optional					
BAF-L-91	-	Additional foam for lateral supports	upcharge	1	100
BAF-L-92	-	Additional foam for hip supports	upcharge	1	100
TRUNK SUPPORTS - optional					
BAF-B-02	7	Second Spine Backrest	Required	1	2 000
BAF-V-01-S BAF-V-01-M BAF-V-01-L	8	Vest Size S M L		1	150
FOOT SUPPORT ACCESSORIES - optional					
BAF-F-03-L BAF-F-03-XL BAF-F-03-XXL	15	3D Foot platforms with heel supports Size S M L	Upgrade option.	1	100
ADDITIONAL ACCESSORIES – required indicated by*					
BAF-A-06 *	-	Electric actuator		1	2 335
BAF-T-02	17	Tray padding		1	120
BAF-A-03 *	-	Battery and charger	1 set.	1	1 500
BAF-A-04	18	Respiratory shelf		1	200
BAF -A-05-01 BAF -A-05-02	-	Extra set of standard upholstery Size S M L		1	400


PICTURES




1 Butterfly headrest 22 cm x 18 cm / 8.6" x 7"	2 Headrest enlarged on both sides Up horn + 5 cm / +1.9" Bottom horn + 7 cm / +2.7"	3 Headrest right side enlarged Up horn + 5 cm / +1.9" Bottom horn + 7 cm / +2.7"	4 Headrest left side enlarged Up horn + 5 cm / +1.9" Bottom horn + 7 cm / +2.7"
---	--	---	--




5 Lateral supports Spacing: 13 - 30 cm / 5.1" x 11.8 "			6 Lateral supports with elbow blocks Spacing: 13 - 30 cm / 5.1" x 11.8 "			7 Second Spine Backrest
Size S Pad: 8 x 9 cm / 3.1" x 6.5"	Size M Pad: 13 x 11 cm / 5.11" x 4.3"	Size L Pad: 16 x 15 cm / 6.2" x 5.9"	Size S Pad: 8 x 9 cm / 3.1" x 6.5"	Size M Pad: 13 x 11 cm / 5.11" x 4.3"	Size L Pad: 16 x 15 cm / 6.2" x 5.9"	



8 Vest S A - 26 cm / 10.2" B - 17 cm / 6.6" C - 14 cm / 5.5" D - 8,5 cm / 3.3" E - 15 cm / 5.9"	8 Vest M A - 31 cm / 12.2" B - 20 cm / 7.8" C - 16 cm / 6.2 " D - 8,5 cm / 3.3" E - 19 cm / 7.4"	8 Vest L A - 39 cm / 15.3" B - 25 cm / 9.8" C - 20 cm / 7.8" D - 8,5 cm / 3.3" E - 29 cm / 11.4"	9 Pelvic belts S: A - 31 cm / 12.2" B - 8 cm / 3.1"	9 Pelvic belts M: A - 35 cm / 13.77" B - 8 cm / 3.1"	9 Pelvic belts L: A - 45 cm / 17.7" B - 10 cm / 3.9"	10 Legs abduction belts
--	---	---	---	--	--	-----------------------------------



11 Kneepad supports	12 Foot platforms S Length: 18 cm / 7" Width: 9 cm / 3.5"	12 Foot platforms M Length: 21,5 cm / 8.4" Width: 12 cm / 4.7"	12 Foot platforms L Length: 24 cm / 9.4" Width: 12 cm / 4.7"
-------------------------------	---	--	--



13 Feet straps (set of 4 pcs.)	14 Foot platforms height adjustment Range 1 Range 2 Range 3	14 Dimension G 12 - 30 / 4.7"- 11.8" 16 - 33 / 16"- 6.3" 26 - 44 / 10.23" - 17.3"	15 3D Foot platforms with heel supports Heel supports height: 12 cm
--	--	--	---



16 Tray	17 Tray padding	18 Respiratory shelf
-------------------	---------------------------	--------------------------------