


Essentials Rollator Parts List

Model #: 700-1500



Call out	Part #	Description	UoM
1	700-1600	Brake handle	EA
2	700-1103	Backstrap	EA
3	700-1605	Seat	EA
4	700-1610	Underseat basket	EA
5	700-1110	C-Clip	EA
6	700-1620	Front fork tube w. wheel	EA
7	700-1625	Rear leg with brake wheel (Left)	EA
8	700-1630	Rear leg with brake wheel (Right)	EA
9	700-1125	Plastic seat support	EA
10	700-1127	Handle bar - Lock	EA
-	520-1100	Cane holder	EA
-	520-1200	Bottle holder	EA

Click [part number](#) to see parts images