

Price List valid from 01.04.2024

LORI STANDER

3 in 1, prone, upright and supine standing up to 60° abduction.

LIW's order number:	Date:
Dealer name:	Customer's order number:
Invoice information:	
Street:	
Postal Code:	
Shipping information:	
Street:	
Postal Code:	
Tel.	
Ordered For:	
Quantity:	
PRICE FOR STANDARD CONFIGURATION (STD): size 1 3 285 CAD	
PRICE FOR STANDARD CONFIGURATION (STD): size 2 3 285 CAD	



COLOUR




Sky blue - standard



Pink - available on request.
Lead time will vary.

BASIC PARAMETERS	Lori Size 1. baby	Lori Size 2. junior
Height	65-100 cm / 25.6"-39.4"	90-135 cm / 35.4"-53.1"
Max. Weight	35 kg / 77 lb	50 kg / 110 lb

DIMENSIONS		
A		cm / inches
B		cm / inches
C		cm / inches
D		cm / inches
E		cm / inches
Weight		Kg / lb
Additional notes		



Code no.	Pict. no.	Name	Comments	Qty.	Price (CAD)
LO1 or LO2	-	Frame with wheels	Prone, upright and supine standing frame up to 60° abduction	1	STD
LO-C-01-1 LO-C-01-2	9	Trunk support size 1 size 2		1	STD
LO-B-01-1 LO-B-01-2	12	Trunk positioning strap size 1 size 2		1	STD
LO-L-01-S LO-L-01-M	4	Trunk laterals size S size M	Left and right	2	STD
LO-C-02-1 LO-C-02-2	11	Pelvic support size 1 size 2		1	STD
LO-B-01-1 LO-B-01-2	12	Pelvic positioning strap size 1 size 2		1	STD
LO-L-01-S LO-L-01-M	4	Hip laterals size S size M	Left and right	2	STD
LO-K-01-1 LO-K-01-2	14	Knee supports size 1 size 2	Left and right	2	STD
LO-F-03-S LO-F-03-M LO-F-03-L	15,16	3D rotational foot platforms with straps size S M L	Left and right	2	STD
LO-T-01	19	Tray (acrylic)		1	STD

HEADRESTS - optional

LO-H-01-1 LO-H-01-2	1,2	Butterfly headrest for supine positioning with mounting system size 1 size 2		1	360
LO-H-02-1 LO-H-02-2	1,3	Butterfly headrest for prone positioning with mounting system size 1 size 2	Available only together with 3D back rotational support with pelvic supports (picture no. 11)	1	450
LO-H-03-1 LO-H-03-2	1,3	Butterfly headrest for prone positioning with mounting system size 1 size 2	Available without 3D back rotational support with pelvic supports (picture no. 11)	1	550

LATERAL SUPPORTS - optional

LO-L-02	5	Elbow blocks	For use with Trunk laterals	2	175
LO-L-06-1 LO-L-06-2	6	Trunk laterals Swing Away size 1 size 2	Upgrade option	2	375
LO-L-04	8	Elbow blocks for swing away laterals	For use with swing away laterals	2	175




















TRUNK AND HIP SUPPORT ACCESSORIES - optional

LO-V-01-S LO-V-01-M LO-V-01-L	7	Vest with attachment size S M L		1	150
LO-L-05-1 LO-L-05-2	11	3D back rotational support with pelvic supports size 1 size 2		1	550
LO-B-03-1 LO-B-03-2	13	Corrective hip belt size 1 size 2		1	200

ADDITIONAL ACCESSORIES - optional

LO-T-02	17	Tray with bowl and cover	Upgrade option	1	300
LO-T-03	18	Padding for elbows	For use with both trays	2	80

PICTURES

<p>1</p> 	<p>2</p> 	<p>3</p> 	<p>4</p> 
<p>Butterfly headrest Size 1 Size 2 A- 23 cm / 9" A- 38 cm / 15" B- 18 cm / 7" B- 25 cm / 9,8" H- 12 cm / 4,7" H- 16 cm / 6,2"</p>	<p>Butterfly headrest for supine positioning</p>	<p>Butterfly headrest for prone positioning</p>	<p>Trunk or hip laterals Pad in size S: 8,5 cm x 9 cm / 3,35" x 3,5" Pad in size M: 15,5 cm x 14 cm / 6,1" x 5,5"</p>
<p>5</p> 	<p>6</p> 	<p>7</p> 	
<p>Elbow blocks Size: 9 cm x 8,5 cm / 3,5" x 3,35"</p>	<p>Trunk laterals Swing Away Pad in size S: 8,5 cm x 9 cm / 3,35" x 3,5" Pad in size M: 15,5 cm x 14 cm / 6,1" x 5,5"</p>	<p>Vest Size S Size M Size L A - 26 cm / 10,2" A - 31 cm / 12,2" A - 39 cm / 15,35" B - 17 cm / 6,7" B - 20 cm / 7,9" B - 25 cm / 9,8" C - 14 cm / 5,5" C - 15-16 cm / 6,1" C - 19-20 cm / 7,7" D - 8,5 cm / 3,35" D - 8,5 cm / 3,35" D - 8,5 cm / 3,35" E - 15 cm / 5,9" E - 19 cm / 7,5" E - 29 cm / 11,4"</p>	
<p>8</p> 	<p>9</p> 	<p>10</p> 	<p>11</p> 
<p>Elbow blocks for trunk laterals swing away Size: 9 cm x 8,5 cm / 3,5" x 3,35"</p>	<p>Trunk support Size 1: 20 cm x 17 cm / 7,9" x 6,7" Size 2: 25 cm x 22 cm / 9,8" x 8,7"</p>	<p>Pelvic support Size 1: 20 cm x 14,5 cm / 7,9" x 2,6" Size 2: 25 cm x 20,5 cm / 9,8" x 8,1"</p>	<p>Back rotational support with pelvic supports Size 1: 20 cm x 14,5 cm / 7,9" x 2,6" Size 2: 25 cm x 20,5 cm / 9,8" x 8,1"</p>
<p>12</p> 	<p>13</p> 	<p>14</p> 	<p>15</p> 
<p>Trunk or Pelvic positioning strap Size 1: (perimeter 25-65 cm / 9,8"-25,6" Size 2: (perimeter 54-71cm / 21,2"-28"</p>	<p>Corrective hip belt Size 1: perimeter max 52cm / 20,5" Size 2: perimeter max 60cm / 23,6"</p>	<p>Knee supports Size S: 7,5 x 5,5 cm / 2,95" x 2,2" Size M: 9,5 x 7,5 cm / 3,7" x 2,95"</p>	<p>3D rotational foot platforms Size S: 18x9 cm / 7" x 3,5" Size M: 21,5 x 12 cm / 8,5" x 4,7" Size L: 24 x 12 cm / 9,4" x 4,7"</p>
<p>16</p> 	<p>17</p> 	<p>18</p> 	<p>19</p> 
<p>Foot straps</p>	<p>Tray with bowl and cover</p>	<p>Padding for elbows</p>	<p>Tray (acrylic)</p>