

Specifications

| Colour | Model Numbers | |
|---|-----------------------------------|------------------------------------|
| ■ Sky Blue ■ Midnight Blue | 600-210006 600-210007 | 600-190006 600-190007 |
| Overall width | 24.75" / 63 cm | 24.75" / 63 cm |
| Folded width | 10.5" / 26.7 cm | 10.5" / 26.7 cm |
| Length | 28.25" / 71.75 cm | 28.25" / 71.75 cm |
| Handle height range | 32.4" – 40.2" 82.3 cm – 102 cm | 29.4" – 37.2" 74.7 cm – 94.5 cm |
| Distance between handles | 17.75" / 45 cm | 17.75" / 45 cm |
| Seat Height | 21" / 54 cm | 19" / 48 cm |
| Seat width | 16" / 40.6 cm | 16" / 40.6 cm |
| Seat depth (rollator) | 12.5" / 31.75 cm | 12.5" / 31.75 cm |
| Seat depth (transport chair) | 16" / 40.6 cm | 16" / 40.6 cm |
| Weight (Transport) | 30 lbs (13.6 kg) | 30 lbs (13.6 kg) |
| Weight (Total) | 32 lbs (14.5 kg) | 32 lbs (14.5 kg) |
| Weight Capacity | 275 lbs (125 kg) | 275 lbs (125 kg) |



SKY BLUE



MIDNIGHT BLUE




Limited lifetime warranty

Please visit:

triumphmobility.com

for additional information



SCAN ME!

Contact us

Triumph Mobility

151 Bentley Ave, Unit 4

Nepean, Ontario, K2E 6T7

Tel: +1 (855) 546 0711

Fax: +1 (866) 225 3059

Email: info@triumphmobility.com

www.triumphmobility.com



TriumphMobility

©2020Triumph Mobility Inc. All rights reserved.
 Triumph™ and Prestige™ are trademarks of
 Triumph Mobility Inc. The Triumph Prestige
 rollator and transport chair is covered by pending
 patent and industrial design protection.

Printed in Canada

File name: R_202204_TriumphPrestige_brochure_EN

Proudly designed in Canada



TRIUMPH
prestige

The Prestige Lifestyle

- ✓ No tools needed
- ✓ No fussy instructions

Finally, a support to keep pace with your busy life.



ROLLATOR MODE

TRANSPORT CHAIR MODE



Perfect Flexibility

10.5 in




The Prestige locks closed at **10.5 in** at its widest point – perfect for tucking it away in a restaurant or at home. Its economic size also makes it easy to fit in a car trunk, weighing only **13.6 kg**



The Prestige transforms from rollator to transport chair in seconds: **swivel the flexible padded back rest, flip the handgrips around**, you're ready to move!



Don't worry about wearing it out: the frame has a **lifetime warranty** and all other parts are guaranteed for 2 years.



≤ 11 lbs (5 kg)

The basket underneath makes the Prestige perfect for shopping, securely attached at six different points and strong enough to **hold up to 11 lbs.**

What else?

Comfortable



- The padded backrest adjusts to the natural angle in both modes
- A padded seat cushion with a roomy seat depth of 12.5" in rollator mode and 16" in transport chair mode
- Ergonomic, angled and adjustable handgrips
- Padded armrest at a natural position of 9.5" from seat to arm

Stylish



- Two fashionable colours to choose from
- Wires all internal for a sleek look with nothing to snag

Safe



- Soft tires provide a natural suspension and balance
- Drum brakes provide reliable all-weather braking
- Sturdy aluminum frame design