

TRAVEL SIT


Travel SIT is an orthopedic seat for children.

LIW's order number:	Date:
Dealer name:	Customer's order number:
Invoice information:	
Street:	
Postal Code:	
Shipping information:	
Street:	
Postal Code:	
Tel.	
Ordered For:	
Quantity:	
PRICE FOR STANDARD CONFIGURATION (STD):	1 750 CAD




BASIC PARAMETERS	SIZE	Travel SIT S	Travel SIT M
Seat depth	A	26 cm / 10.2"	31 cm / 12.2"
Back length	B	52 cm / 20.5"	62 cm / 24.4"
Seat width	C	17-24 cm / 6.7"-9.4"	22-35 cm / 8.6"-13.8"
User's height	D	70-105 cm/27.5"-41.3"	95-125cm/37.4"-49.2"

COLOUR











DIMENSIONS		
A		cm/inches
B		cm/inches
C		cm/inches
D		cm/inches
E		cm/inches
F		cm/inches
G		cm/inches
Additional notes		



Code no.	Pict. no.	Name	Comments	Qty.	Price
TRSS-01 or TRSM-01	-	Seat		1	STD
TRS-H-01	1	Headrest		1	STD
TRS-L-01	2	Lateral and hip supports	Set of 10 pcs	1	STD
TRS-B-01	3	Chest and hip straps	Set of 2 pcs	1	STD
TRS-B-02	4	Pelvic straps		1	STD
TRS-A-01	5	Backpack		1	STD
TRS-A-02	6	Abductor		1	STD
TRS-V-01-S TRS-V-01-M	7	Vest S M	Choose your size – see picture 7 below	1	STD
TRS-A-03	8	Lounger frame		1	STD

PICTURES

<p>1</p> 	<p>2</p> 	<p>3</p> 
<p>Headrest</p>	<p>Lateral and hip supports</p>	<p>Chest and hip straps</p>
<p>4</p> 	<p>5</p> 	<p>6</p> 
<p>Pelvic straps Size S 31 cm x 8 cm / 12.2" x 3.2"</p>	<p>Pelvic straps Size M 35 cm x 8 cm / 13.8" x 3.2"</p>	<p>Backpack</p>
<p>7</p> 	<p>8</p> 	<p>Abductor</p>
<p>Vest S A - 26 cm / 10.2" B - 17 cm / 6.7" C - 14 cm / 5.5" D - 8,5 cm / 3.3" E - 15 cm / 5.9"</p>	<p>Vest M A - 31 cm / 12.2" B - 20 cm / 7.9" C - 16 cm / 6.3" D - 8,5 cm / 3.3" E - 19 cm / 7.5"</p>	<p>Lounger frame</p>