



**LIW Care**  
Medical Devices inc.

Tel: 1-613-258-3344


Email: [inquires@tlcmedicalsupply.ca](mailto:inquires@tlcmedicalsupply.ca) & [mario@tlcmedicalsupply.ca](mailto:mario@tlcmedicalsupply.ca)




## TravelSIT

TravelSIT is a lightweight and mobile system for stabilizing your baby in sitting, lying down and intermediate positions. The seat fits into a convenient, dedicated backpack, which will be your favorite travel companion from now on.

The backrest consists of four infinitely adjustable torso supports, which, combined with two seat belts and an anatomical headrest, are designed to support your child. The seat is equipped with six supports with which you stabilize your child and adjust the width of the seat. Also included is an adaptive wedge for lower limb inversion and pelvic stabilization straps. An additional accessory is the „Lounger”, which allows you to use the seat on the floor. All this and more available with the TravelSIT.



## COLOURS



**LIW Care**  
Medical Devices inc.

Go on a journey. and Travel SIT will take care of the correct position of your child.



Integration with siblings and peers.



Dinner together in a restaurant  
At one table.



Playground, and on it a swing.



Fun at the terrace table.



Staying with grandparents.



Holidays and no access to specialized equipment.

## TECHNICAL DATA

Patient's height

## SIZE S

70÷105 cm / 27.5" - 41.3"

## SIZE M

95÷125 cm / 37.4" - 49.2"

Back length

52 cm / 20.5"

62 cm / 24.4"

Seat depth

26 cm / 10.2"

31 cm / 12.2"

Seat width

17÷24 cm / 6.7" - 9.4"

22÷35 cm / 8.6" - 13.8"