

# Enhance Your Reporting Workflow

cvi42 offers customizable reporting for cardiac MR and CT evaluation. Build the report you need for research and daily clinical practice.

## Standard Reporting Features - Included with all cvi42 licenses

- Simple Clinical and Scientific data reporting
- Customizable template builder
- Infographics, diagrams, screenshots and measurements
- Export PDF, text, DICOM secondary capture and DICOM encapsulated pdf
- User configurable normal values and automatic classification
- Updated normal values sourced from the Healthy Hearts Consortium\*
- DICOM SR (Structured Reporting) - export to any configured DICOM node (e.g., third-party reporting systems)

## Direct PowerScribe One Integration\*\*

- Send report data directly to PowerScribe One to improve workflow and avoid transcription errors
- Populate your templates within PowerScribe One to customize your report content and format

## Upgrade to cvi42 | Report

To further enhance your reporting capabilities, cvi42 offers an additional license available for purchase – cvi42 | Report.\*\*

## cvi42 | Report - Advanced Features Include the following:

- All standard reporting functionalities
- Z-Scores
- Advanced Database Search
- Browser Access
- HL7 Integration\*\*\*

\*<https://www.jacc.org/doi/10.1016/j.jcmg.2024.01.009>

\*\*Additional license required

\*\*\*HL7 integration requires the purchase of HL7 Integration Services.

Note: Both the standard report and cvi42 | Report contents require end users to verify the contents (including automatically generated contents)

**Brief Summary:** Indications, contraindications, warnings and precautions can be found in the product labelling.

**Disclaimer:** Not all modules or features are available in every region. Contact your local Circle representative for all regional availability.

**CAUTION:** Federal law (USA) restricts these devices for sale by, or on the order of a physician. The system is intended for use only by trained Healthcare Professionals



At the Heart of  
**IMAGING**

**cvi42**®