

cvi42 | Report[®]

v.3.10.3 (1406)

Release Notes

What's New in Version 3.10.3

Reporting

This minor release introduces no new user features; it focuses on security improvements and bug fixes.

Resolved Issues

The following list contains resolved usability related issues, in addition to all other internal issues resolved during the product development process.

Number	Summary
CV-18333	cvi42 Report42: Patient/Study information pop-ups dont adjust to screen size
CV-19798	Report has too much whitespace and ends up taking too many pages
CV-19800	Study indication is not on its own line and when a lot of data is added it wraps
CV-19914	Regional LV Tissue section diagram is clipped at the right hand side
CV-19952	Addendum should be renamed
CV-20337	Report page edit and preview modes are sometimes missing report images
CV-20575	Referring physician handling in cvi42 report
CV-20574	Textbox for patient history not clear

Known Issues

Number	Summary
CV-6966	cvi42 Report: Importing Findings/Summary Text in Structured Findings also imports Summary Text <i>Workaround: Delete the unnecessary Summary Text from the Structured Findings section in the Admin pages</i>
CV-9429	cvi42 Report: Too many images causes extreme lag <i>Workaround: Delete report and send again</i>
CV-9854	cvi42 Report: Text appears in front instead of behind buttons <i>Workaround: scroll past the header, the cvi42 Report header text is gone</i>

CV-21082	User is logged out after switching between Worklist and Report multiple times in a row <i>Workaround: User is able to log back in immediately.</i>
CV-21261	Findings/Summary Text Canned Text Checkboxes not Functional when Creating New Text (SQL Sever ONLY) <i>Workaround: After creating the Findings/Summary Text, edit it using the pencil icon on the Findings/Summary Text page. On the edit page, input the desired availability using the canned text checkboxes.</i>

External Libraries – Copyright and Licensing Requirements

This is a list of free and open source software used in our products, including any modifications made by Circle Cardiovascular Imaging Inc., throughout the course of our development. To access the most up-to-date information on the open source software we use, visit:
https://www.circlecvi.com/docs/product-support/Open_Source_Software.pdf

DCMTK v3.6.7

Project Website
Licenses

Download

DICOM Toolkit (DCMTK) Project
DCMTK License
DCMJ2K Commercial License
DCMTK v3.6.7 Source Code (.zip)

Apache HTTP Server 2.4.54

Project Website
License
Download

HTTP Server Project
Apache License
Source Code and Binaries

OpenSSL v1.1.11

Project Website
License
Download

OpenSSL Project
OpenSSL License
OpenSSL Source Code

PHP 7.4.32

Project Website
License
Download

PHP Project
PHP License
Source Code and Binaries

jQuery v3.6.0

Project Website
License

jQuery Project
MIT License

Fancybox v1.3.4

Project Website
License

fancybox Project
Dual licensed under the MIT and GPL licenses

HTML5 Shiv v3.6.2pre

Project Website
License
Download

HTML5 Shiv Project
Dual licensed under the MIT or GPL Version 2 licenses
HTML5 Shiv Source Code

hoverIntent r6 2011.02.26

Project Website
License
Download

hoverIntent r6 Project
MIT and GPL licenses
hoverIntent Source Code

Superfish v1.4.8

Project Website
License

Superfish Project
Dual licensed under the MIT and GPL licenses

ImageMagick v7.1.0-4

Project Website

ImageMagick Project

License	Apache License
GPL GhostScript v9.54.0	
Project Website	GhostScript Project
License	GNU Public License v3
jQuery Validation Plugin v1.11.1	
Project Website	jQuery Validation plugin Project
License	MIT License
jQuery UI v1.10.3	
Project Website	jQuery UI Project
License	MIT License
jQuery bglframe v2.1.1	
Project Website	jQuery bgiFrame Project
License	Dual licensed under the MIT and GPL licenses
jQuery chronoStrength	
Project Website	jQuery chronoStrength Project
License	Distributed under a Create Commons by-sa license
jQuery chronoStren Cookie plugin	
Project Website	jQuery Cookie plugin Project
License	Dual licensed under the MIT and GPL licenses
jQuery JSON Plugin v2.6.0	
Project Website	jQuery JSON plugin Project
License	MIT license
jQuery ScrollTo Plugin v1.4.2	
Project Website	jQuery ScrollTo plugin Project
License	Dual licensed under MIT and GPL licenses
mPDF v6.1.3	
Project Website	mPDF Project
License	Licensed under GNU General Public License (GPL) v2
jscolor 1.4.2	
Project Website	jscolor Project
License	Licensed under the GNU Lesser General Public License (LGPL)
jQuery ScrollToFixed	
Project Website	ScrollToFixed Project
License	MIT license
Copyright	Copyright and permission notice
Mirth Connect v3.1.1	
Project Website	Mirth Connect Project
License	Mirth Connect is released under the Open Source Initiative (OSI) approved MPL 1.1.
AngularJS v1.8.2	
Project Website	AngularJS Project
License	MIT license
Select2 v4.0.13	
Project Website	Select2 Project
License	MIT license
Sortable v1.9.0	
Project Website	Sortable Project
License	MIT license
Sprintf v1.0.3	
Project Website	Sprintf Project
License	Released under the 3-clause BSD license

Go v1.16.0

Project Website
License

The Go Programming Language
Released under the 3-clause BSD license

JSPDF v1.5.3

Project Website
License

Parall
MIT license

Canvg v1.5.0

Project Website
License

Canvg Project
MIT license

ResumableJS v1.0.0

Project Website
License

ResumableJS Project
MIT license

odbc v0.0.0-20210605012845

Project Website
License

Github ODBC
BSD 3-Clause License

Logrus runtime formatter v0.0.0-20190729070250

Project Website
License

Github BanzaiCloud
MIT license

Dgrijalva/jwt-go v3.2.0

Project Website
License

GitHub dgrijalva/jwt-go
Apache 2.0 License

go-openapi/errors v0.20.1

Project Website
License

GitHub go-openapi/errors
Apache 2.0 License

go-openapi/loads v0.20.2

Project Website
License

GitHub go-openapi/loads
Apache 2.0 License

github.com/go-openapi/runtime v0.19.31

Project Website
License

GitHub go-openapi/runtime
Apache 2.0 License

github.com/go-openapi/spec v0.20.3

Project Website
License

GitHub go-openapi/spec
Apache 2.0 License

github.com/go-openapi/strfmt v0.20.2

Project Website
License

GitHub go-openapi/strfmt
Apache 2.0 License

github.com/go-openapi/swag v0.19.15

Project Website
License

GitHub go-openapi/swag
Apache 2.0 License

github.com/go-openapi/validate v0.20.2

Project Website
License

GitHub go-openapi/validate
Apache 2.0 License

github.com/google/go-cmp v0.5.6

Project Website
License

GitHub google/go-cmp
BSD-3-Clause License

github.com/gorilla/mux v1.8.0

Project Website
License

GitHub gorilla/mux
BSD-3-Clause License

github.com/gorilla/websocket v1.4.2

Project Website [GitHub gorilla/websocket](#)
License [BSD-2-Clause License](#)

github.com/jesevdk/go-flags v1.5.0

Project Website [GitHub jessevdK/go-flags](#)
License [BSD-3-Clause License](#)

github.com/jinzhu/gorm v1.9.16

Project Website [GitHub jinzhu/gorm](#)
License [MIT license](#)

github.com/kardianos/service v1.2.0

Project Website [GitHub kardianos/service](#)
License [Zlib License](#)

github.com/leonelquinteros/gotext v1.4.0

Project Website [GitHub leonelquinteros/gotext](#)
License [MIT license](#)

github.com/satori/go.uuid v1.2.0

Project Website [GitHub satori/go-uuid](#)
License [MIT license](#)

github.com/sirupsen/logrus v1.8.1

Project Website [GitHub sirupsen/logrus](#)
License [MIT license](#)

github.com/spf13/pflag v1.0.5

Project Website [GitHub spf13/pflag](#)
License [BSD-3-Clause License](#)

go.uber.org/zap v1.19.1

Project Website [GitHub uber-go/zap](#)
License [MIT license](#)

golang.org/x/net v0.0.0-20210928044308-7d9f5e0b762b

Project Website [cs.opensource.google/go/x/net](#)
License [Redistributable License](#)

gopkg.in/natefinch/lumberjack.v2 v2.0.0

Project Website [GitHub natefinch/lumberjack](#)
License [MIT license](#)

Support

For further information on installing/exchanging the dynamic libraries please contact Circle Cardiovascular Imaging support.

Phone: +1 (403) 338-1870

E-mail: support@circlecvi.com

Circle Cardiovascular Imaging Inc.

1100, 800 5th Avenue SW

Calgary, AB, Canada T2P 3T6

Phone: +1 (587) 747 4692

Website: circlecvi.com