



Queensland Mobility Solutions

ETAC SWIFT SHOWER CHAIR



Care Instructions

- Clean the product with a solvent-free cleaning agent with a pH level between 5–9, or with a 70% disinfectant solution.
- The product can be cleaned in a cabinet washer-disinfector at max. 85 °C for 3 minutes

- Overall height: 800 mm - 950 mm
- Seat depth: 500 mm
- Recess depth x width: 130 mm x 115 mm
- Seat area: 450 mm x 410 mm
- Seat height: 420 mm - 570 mm
- Seat width: 530 mm - 540 mm x 565 mm - 595 mm
- Width between arm supports: 450 mm
- Overall width: 545 mm
- Product weight: 4.7 kg
- Material: Polypropylene, aluminium
- Colour: Grey, White, Lagoon Green, Blue

Shop 3, 73-79 Wisers Road,
MAROOCHYDORE QLD 4558

07 54795918

info@qldmobilitysolutions.com.au

www.qldmobilitysolutions.com.au

