



# Queensland Mobility Solutions

## Rollz Motion 2-in-1 Walker Wheelchair



### Rollz Motion – Regular

- Tyres: PU
- Weight: 11.6 kg (wheelchair attachment 3.4 kg)
- Max user weight: 125 kg
- User height: 1.60 – 1.90 m
- Seat height: 55 cm
- Handle height: 87 – 98 cm
- Seat size (rollator): 20 × 45 cm
- Seat size (wheelchair): 42 × 45 cm
- Overall length: 63 cm
- Overall width: 67 cm
- Folded width: 31 cm
- Wheel diameter: front 8", rear 12"

### Rollz Motion – Small

- Tyres: PU
- Weight: 11.6 kg (wheelchair attachment 3.4 kg)
- Max user weight: 125 kg
- User height: 1.55 – 1.85 m
- Seat height: 50 cm
- Handle height: 82 – 94 cm
- Seat size (rollator): 20 × 45 cm
- Seat size (wheelchair): 42 × 45 cm
- Overall length: 63 cm
- Overall width: 67 cm
- Folded width: 31 cm
- Wheel diameter: front 8", rear 12"

Shop 3, 73-79 Wisers Road,  
MAROOCHYDORE QLD 4558

07 54795918

[info@qldmobilitysolutions.com.au](mailto:info@qldmobilitysolutions.com.au)

[www.qldmobilitysolutions.com.au](http://www.qldmobilitysolutions.com.au)

