

# 92405NO Trauma 2

# 17x17cm

## First Aid Bandage

**A compression bandage with an extra loose bandage. The bandage is specially designed for military and emergency use.**

### Advantages of our product

1. The dressing is easy to use to stop bleeding and protect the injury from further impact
2. Produced in Sweden
3. Supplier to the Swedish army since the late 1930s
4. The components are produced within the EU
5. The products do not contain any hard plastic components that risk breaking off when cold
6. Shelf life of 20 years
7. Opens easily with just one hand and teeth
8. Flat and soft packaging
9. Good absorption
10. Low adhesion so it does not stick to wounds



### Features

Bandages:	2 pcs. 3mx6cm <u>white</u>
Compress:	1 st Cellosoft 17x17 cm
Absorption layer:	100% Viscose
Surface layer:	100% LDPE-HDPE

MANUFACTURER	TRADEMARK	CATEGORY		URSPRUNGSLAND	PRODUCT CODE	REGULATION CONFORMITY
AKLA	AKLA	TRAUMA BANDAGE		SWEDEN	30059031	93/42/EEC
SHELF LIFE				MEDICAL DEVICE	QUALITY CERTIFICATION	ENVIRONMENTAL CERTIFICATION
20 YEARS				1s (sterile)	ISO 13485	ISO 14001
NSN		NIIN		GMDN	CPV CODES	UNCCS CODES
FSC 6510		257220160		47011	33141000	481930,481931,481940
6510-25-722-0160					33141110	
NSN Created on 1 Jan 1984					33141113	
LOGISTIK		AMOUNT	MEASUREMENTS	WEIGHT		EAN / GS1
CU		1 pc	100x70x30mm			7311410924057
CU / SKU						
CU / DU		75 pc	390x300x230mm	2,7 KG		7311410224058
KFP / Pallet		3600 pc	1200x800x1530mm	154,6 KG		

