

95490 Refill – plaster dispenser

Refill for AKLAs plaster dispenser - Beige


Language: SE/GB/NO/DK/FI/DE/NL

With the AKLA Plaster Dispenser, you can make plasters according to the size and location of the wound. Pull out the tape to the required length and close the lid to get a ready-to-use plaster.

With this refill cassette, it is easy to change/refill when you have run out of plasters or wound pads.

- This product is gentle to your skin
- Contains No health hazard chemical compounds
- Free of PVC, Natural rubber latex, Colophony
- No antimicrobial compounds intentionally added
- Low Adhesion which prevents the tape from being difficult to remove

Instruktion



Contents

Wound pads: Size: 20x25mm / Amount: 100 pcs
Absorption: >500%; <1 sek
100% Polyethylene / 25% Polyethylene + 75% Viscose

Carrier: Size: 9,1mx2,5cm
Backing: Rayon fiber / Binder: Acrylate
Adhesive: Acrylate / Core: Paper material

PRODUCER	BRAND	CATEGORY	COUNTRY OF ORIGIN	TARIF CODE	REGULATION CONFORMITY
AKLA	AKLA	PLASTER	SWEDEN	30051000	MDR 2017/745
SHELF LIFE	CE	MEDICAL DEVICE	QUALITY CERTIFICATION	ENVIRONMENT CERTIFICATION	
5 år	YES	Class 1	ISO 13485	ISO 14001	
NSN	STERIL	GMDN	CPV CODES	UNCCS CODES	
	NO				
ECLASS	UNSPSC				
40020303	42170000				
LOGISTIC	AMOUNT	MEASUREMENTS	WEIGHT	EAN / GS1	
CU	1 st	103x36x107MM	0,073 KG	7311410954900	
CU / DFP	9 PCS			7311410154904	
CU / SKU	180 pcs	600x410x440MM	0,648KG	7311410254901	
CU / Pallet	1440 pcs	1200x800x1030mm	30,184KG		

