

95492 Plaster Dispenser


AKLAs plaster dispenser - Beige

Language: Instruction with pictures adapted for all languages

With the AKLA Plaster Dispenser, you can make plasters according to the size and location of the wound. Pull out the tape to the required length and close the lid to get a ready-to-use plaster.

- There is a refill, item no. 95490 for this product
- This product is gentle to your skin
- Contains No health hazard chemical compounds
- Free of PVC, Natural rubber latex, Colophony
- No antimicrobial compounds intentionally added
- Low Adhesion which prevents the tape from being difficult to remove

Instruction



Contents

Wound pads: Size: 20x25mm / Amount: 100 pcs
Absorption: >500%; <1sek
100% Polyethylene / 25% Polyethylene + 75% Viscose

Carrier: Size: 9,1mx2,5cm
Backing: Rayon fiber / Binder: Acrylate
Adhesive: Acrylate / Core: Paper material

PRODUCER	BRAND	CATEGORY	COUNTRY OF ORIGIN	TARIF CODE	REGULATION CONFORMITY
AKLA	AKLA	PLASTER	SWEDEN	30051000	MDR 2017/745
SHELF LIFE	CE	MEDICAL DEVICE	QUALITY CERTIFICATION	ENVIRONMENT CERTIFICATION	
5 år	YES	Class 1	ISO 13485	ISO 14001	
NSN	STERIL	GMDN	CPV CODES	UNCCS CODES	
	NO				
ECLASS	UNSPSC				
40020303	42311505				
LOGISTIC	AMOUNT	MEASUREMENTS	WEIGHT	EAN / GS1	
CU	1 st	145x45x120MM	0,237KG	7311410954924	
CU / DFP					
CU / SKU	15 st	480x250x130MM	3,722KG	7311410254925	
CU / Pallet	720 st	1200x800x1190	203,656 KG		

