



USE AND CARE INSTRUCTIONS for your Invisalign® aligners

WEARING AND USING YOUR INVISALIGN® ALIGNERS

Always Remember To:

1. Wear your aligners per your doctor's instructions, usually 22 hours per day, not leaving out for longer than 40 minutes at a time. When they are not in your mouth, they are not working which will hamper treatment time. Relapse can occur if not worn as directed.
2. Wash your hands thoroughly with soap and water before handling your aligners.
3. Handle only ONE aligner at a time.
4. Rinse your aligners when removing them from the packaging.
5. Aligners should only be removed when brushing, flossing, or eating.
6. When aligners are removed, they should always be placed in their protective case [never in a napkin, on the table, etc.]

Attachments:

1. Attachments are made of clear composite material and are placed on the teeth to help the aligners make certain movements. It is very important the attachments or "bumps" on the teeth engage the corresponding "bubbles" in the aligner.
2. You must check attachments for proper engagement with each new set. If it is not fully engaged, the attachment can possibly work against the treatment plan.

Aligner Insertion:

1. Make sure you have the proper aligner—the upper for your top teeth and the lower for your bottom teeth.
2. To help avoid confusion, each aligner is engraved with your unique case number, a "U" for upper and an "L" for lower, followed by the stage number.
3. You may insert either the upper or lower aligner first. When inserting each aligner, gently push the aligners over your front teeth. Then apply equal pressure, using your fingertips, to the tops of your left and right molars (back teeth) until the aligner snaps into place.
4. Do NOT bite your aligners into position. This may damage them. Note: If you experience sharp pain or significant discomfort, discontinue use of the aligners and contact your doctor.

Aligner Removal:

1. On one side of your mouth, use your fingertip/nail on the inside of your back molar to slowly pull the aligner from your molars then walk the aligner out to the front and then the other side to release the aligner.
2. Immediately rinse aligner with water, shake off excess water, and store your aligners in the protective case provided with your starter kit.

Cleaning aligners:

1. Clean your aligners prior to each insertion. Use a soft bristle toothbrush with water and a small amount of toothpaste if you like. Note: Be sure to rinse each aligner thoroughly with water after each cleaning. We also recommend that you use the Invisalign aligner cleaning products once a week, or as needed. Cleaning crystals are available for purchase at www.invisalignaccessories.com.
2. There are other available accessories such as aligner removing instrument which can be purchased online on amazon or in general pharmacies.
2. Do NOT use denture cleaners to clean aligners. Do NOT soak them in mouthwash. These products can damage the surface of the aligner, causing it to become dull and more visible.

Proper Oral Hygiene:

1. Remove your aligners for eating and drinking. (You do not need to remove your aligners to drink cool water.)
2. Brush and floss your teeth after each meal or snack prior to re- inserting your aligners. If you don't have access to your cleaning system or a toothbrush, you can simply rinse your mouth, and then clean your aligners by holding them under warm running water. It's not the best way to clean, but it works in a pinch. Be sure to thoroughly clean your aligners at your earliest convenience.
3. Recommended to have cleaning appointments every 3-4 months during treatment for the continued health of your teeth and gums.

Note: Keep your previous set of aligners in a clean plastic bag, or as instructed by your doctor. Keep them out of reach of small children and pets. When traveling always bring your current set, previous set and next set in case of an emergency.

Refinement:

1. Scanning may be necessary during or at the end of treatment. This "refinement" is not an indication that something is wrong, but simply the doctor's way of making adjustments in the treatment.
2. At the end of the first series of aligners, it may be necessary for a patient to wear another series of refinement aligners before the plan is completed.

Retention:

1. Following your Invisalign treatment you will be required to wear aligners as retainers to prevent relapse or unwanted movements.
2. In most cases we request near full time wear for the first 1 month [16-20 hours/day]. Following this, the retainers should be worn at night for life. "Nighttime for Lifetime"
3. Failure to wear retainers as directed will generally result in relapse and will then require retreatment. Failure to wear aligners as directed would void any warranty associated with this treatment.

FREQUENTLY ASKED QUESTIONS AND ANSWERS:

Q: Will the treatment be painful?

A: Most people experience tooth soreness for a few days after starting each new stage. This is normal. It is a sign that the Invisalign aligners are working, moving your teeth to their final destination. This soreness should gradually go away a couple of days after inserting the new aligner in the series. If it doesn't, promptly call your doctor.

Q: Will wearing the Invisalign aligners affect my speech?

A: Like all orthodontic treatments, the aligners may temporarily affect your speech, and you may have a slight lisp for a day or two. However, as your tongue gets used to having aligners in your mouth, any lisp or minor speech impediment caused by your aligners should disappear.

Q: What should I do if my new Invisalign aligner doesn't snap onto my teeth?

A: Minor discrepancies between the new aligner and the current tooth position are normal, since the teeth need time to conform to the new aligner position. In the event of significant problems with aligner fit, inform your doctor.

Q: Are there restrictions on what I can eat?

A: In general, no. Unlike traditional orthodontics, you can usually eat and drink whatever you desire because you remove your aligners while eating. Thus, there is no need to restrict your consumption of any of your favorite foods and snacks, unless otherwise instructed by your doctor.

Q: Is it OK to drink hot or cold beverages while wearing Invisalign aligners?

A: Except for cool water, we recommend that you do not drink while wearing the aligners. This is to avoid formation of cavities and stains, or warping of the aligners with hot drinks and hot water.

Q: Can I chew gum while wearing Invisalign aligners?

A: NO. Gum will stick to the aligners. We recommend removing your aligners for all snacks and meals.

Q: Will smoking or chewing tobacco stain the aligners?

A: We discourage tobacco use while wearing aligners because of the possibility of aligner discoloration.

Q: Why do some of my aligners have bumps or ridges on them?

A: Depending on your specific treatment, some of the movements may require either "attachments" or "ridges" to help the aligner grip the teeth. These bumps, or wells, are where the aligner grips the attachment the doctor places on your teeth. The attachments are actually small pieces of composite the doctor affixes to your teeth that are then gripped by bumps on your aligners. The ridges are slim indentations in your aligners. Your doctor will use one or both of these features to attempt to achieve the desired movement.

Q: What if I lose or break an Invisalign aligner? **A:** In the event that an aligner is lost or broken, you should immediately inform your doctor. Your doctor will probably tell you to start wearing your last set or next set of aligners immediately. He or she will possibly order you a new set of aligners to replace the ones you just lost, which should arrive in a few days. The key is to not go an extended time without wearing an aligner.

Q: What if I lose or break an attachment?

A: In the event that an attachment is lost or broken, you should contact your doctor. This is not an emergency, but an appointment will be scheduled to have the attachment replaced. Continue wearing your aligners as instructed even if an attachment is missing.

Note: If you have additional questions after reading this pamphlet, please consult with your doctor.