

Lyric



SoundLync Quick Reference



ON 1 beep (remove)

SLEEP 2 beeps (remove)

OFF 2 beeps (hold)
4 beeps

PHONAK

Lyric



SoundLync Quick Reference

VOLUME

2 beeps (remove for 2 sec)

1 beep (hold)

Volume increase starts 1 step louder

Each step: louder beep

Maximum setting: double beep

Note: Lyric is muted during volume adjustment

002 0-12 20-0 3/20 14-0 4/14 Printed in U.S.A. © 2014 All rights reserved

PHONAK