



Request for Medical Clearance for Extended Wear Hearing Aid(s)

Date: _____

Patient Name: _____

Date of Birth: _____

Your patient, _____, is interested in being fit with Lyric hearing aids (please see reverse side for description). Because these hearing aids are deep fitting and are worn for 2-4 months at a time, they present unique fitting requirements. Outlined below are specific concerns for which we are requesting medical clearance to proceed with Lyric as indicated.

Request for Medical Clearance Re:

- ☐ Uncontrolled or insulin-dependent Diabetes
- ☐ Anticoagulant medications
- ☐ Compromised immune system
- ☐ Chemotherapy less than 6 months ago
- ☐ Implantable medical device

☐ Other: _____

Physician Recommendation:

- ☐ Patient can wear Lyric hearing aids
 - ☐ Left ear only
 - ☐ Right ear only
 - ☐ Both ears
- ☐ Patient can not wear Lyric hearing aids and should be assessed for an alternative type of hearing aids

Physician Signature

Date

Request for Cerumen Removal

Lyric fittings require that ALL cerumen is completely removed from the ear canal. There can be no cerumen visible and the canal must be clean and dry prior to inserting Lyric. This patient is being referred for cerumen removal for the purpose of fitting Lyric hearing aids.

Please remove cerumen from:

- ☐ Left ear only
- ☐ Right ear only
- ☐ Both ears

Please Fax Back To:

Clinic Name: _____

Hearing Professional Name: _____

Fax: _____

Hearing Professional Signature



Product Information

Lyric is the first extended wear hearing aid that can be worn in the ear 24 hrs per day, 7 days a week for up to 120 days. Due to the deep placement inside the ear canal, Lyric is 100% invisible and delivers natural sound quality. With Lyric, there are no batteries to change, no ongoing maintenance required and no daily insertion or removal. Developed by ENT consultants, engineers and audiologists, Lyric received FDA approval in the USA in 2007 and Health Canada approval in 2010. Lyric is inserted by a hearing care professional who has undergone special training to become an Authorized Lyric Provider.

Lyric at a glance

Materials

Biocompatible soft foam seals adapt to the shape of the ear canal, preventing the device from slipping while it is in the ear and ensuring a comfortable and secure fit

Size

12 mm long, seals available in various sizes

Battery

The optimized zinc-air system enables Lyric to be worn for an uninterrupted period of up to 120 days. The battery forms part of the module; Lyric is therefore disposable and cannot be reinserted once it has been removed from the ear.

Candidacy Factors

Audiometric configuration (mild to moderate HL), lifestyle, ear dimensions, shape and health. Otoscopy is performed using a microscope or head loupe.

Sizing

Length of the ear canal is measured using a depth sizer (similar to an in-situ probe). Overall fit in the ear canal is determined using a device sizer which represents the device itself.

Fitting

The fitting process is quick and painless. Insertion is performed using calibrated forceps to ensure accurate placement 4mm from the eardrum and the device is programmed wirelessly through Phonak software. The user may experience a little discomfort afterward, in most cases lasting from a couple of hours up to 5 days. If pain is encountered the device should be removed immediately. The user can easily remove the device with the help of a special tool that is provided with the device, but it is preferable to return to the hearing care professional for this wherever possible.

Subscription model

The hearing care professional will supply a sufficient number of hearing aids per year, and will provide comprehensive care and support.

For more information, please visit www.lyrichearing.com

