

# HyperLight-G

[www.eighteeth.com](http://www.eighteeth.com)



# HYPERLIGHT-G



## Receptor

Use with film, digital sensor, phosphor plate




## Large Battery Capacity

More than 300 exposures in one full charge

## Preset Parameter

Divided into adults and children, and can be set and controlled for different teeth



Smart angulation system, faster your workflow

## Locking system





## Safe

- Internal x-ray shielding
- Locking system to avoid unintentional exposure
- Backscatter shield for operator and patient protection(\*optional)

## Portable


- Cordless, handheld and extremely mobile, move freely between rooms, locations and cities

## Economical

- Use with film, digital sensor or phosphor plate
- Prolonged working life, at least 80,000 exposures

# HyperLight-G

[www.eighteeth.com](http://www.eighteeth.com)



Model	: HyperLight-G
Tube Voltage	: 65 kV
Tube Current	: 2 mA
Exposure Time Range	: 0.02-2 Sec
Weight	: 1.9 kg
Battery	: 2000mAh X4
Line Voltage	: 100V-240V, 50Hz/60Hz