

InfuLife Over/Under Fill Chart

Planned Rate		0.5	1	1.5	2	2	4	5	1	2	4	5	10	2	4	5	10
Model		C100R0005	C100R0010	C100R0015	C100R0020	C120R0020	C120R0040	C125R0050	C270R0010	C270R0020	C270R0040	C270R0050	C270R0100	C300R0020	C400R0040	C400R0050	C400R00100
Volume and Rate		100ml x 0.5ml/hr	100ml x 1ml/hr	100ml x 1.5ml/hr	100ml x 2ml/hr	120ml x 2ml/hr	120ml x 4ml/hr	125ml x 5ml/hr	270ml x 1ml/hr	270ml x 2ml/hr	270ml x 4ml/hr	270ml x 5ml/hr	270ml x 10ml/hr	300ml x 2ml/hr	400ml x 4ml/hr	400ml x 5ml/hr	400ml x 10ml/hr
Approx. Delivery Time		Fill Volume (ml)															
12 hrs																	
18 hrs							78	83									
22 hrs							93	104					226				
24 hrs	1 day						101	114					239				
30 hrs							124	145					279				355
38 hrs			41	58	78		155	187					332				388
44 hrs			47	68	90							220					413
46 hrs			49	71	94	89						228					421
48 hrs	2 days		51	74	98	92						236					429
60 hrs			63	93	122	109					249	285					
72 hrs	3 days		75	111	146	126					285	334			337	324	
96 hrs	4 days		98	148	195	160				215	355				383	462	
120 hrs	5 days	57	122							247	425			249	429		
132 hrs	5.5 days	63								264				266	452		
144 hrs	6 days	69								280				284	475		
168 hrs	7 days	81							200	312				320			
192 hrs	8 days	93							216					356			
216 hrs	9 days	105							233								
240 hrs	10 days	118							250								
264 hrs	11 days	130							267								
288 hrs	12 days	142							284								
312 hrs	13 days	154							301								

First Medical Source, LLC