


FINGER SPLINT



PRIMO (M'SIA) MALLETT FINGER SPLINT

- Sturdy, lightweight plastic
- Perforated and shaped for ventilation and comfort
- Protect the distal joint but still allow mobility of the proximal joint

PART#	SIZE	CIR. IN CM OVER FINGER	LENGTH IN CM
JK-FS06-0	Size 0	3.8 – 4.2 cm	5 cm
JK-FS06-1	Size 1	4.2 – 4.8 cm	5 cm
JK-FS06-2	Size 2	4.8 – 5.3 cm	5.6 cm
JK-FS06-3	Size 3	5.3 – 5.5 cm	5.7 cm
JK-FS06-4	Size 4	5.5 – 6 cm	5.8 cm
JK-FS06-5	Size 5	6 – 6.5 cm	6 cm



MEDIROYAL (SWEDEN) OVAL-8 FINGER SPLINT

Stabilizes, unloads and corrects the finger IP joints. The Oval 8 improves the hand function without being bulky. Made of 1,6 mm PP which makes it hygienic, easy to clean and comfortable. Can be individually adjusted with a heat gun. Uses the same sizing as normal standard rings. Sold piece by piece or in an assortment box of 44 pcs. Available sizes: 2-15

Indications: Swan neck deformities, boutonnière- and malletfingers, reduction of lateral deviation or to unload a trigger finger.

OVAL-8	SIZE FITS OVAL-8 CM	SIZE FITS OVAL-8 (inch)
2	4.3 cm	1 3/4"
3	4.5 cm	1 13/16"
4	4.95 cm	1 15/16"
5	5.27 cm	2 1/6"
6	5.7 cm	2 3/16"
7	5.86 cm	2 5/16"
8	6.14 cm	2 7/16"
9	6.35 cm	2 1/2"
10	6.82 cm	2 11/16"
11	7.09 cm	2 13/16"
12	7.4 cm	2 15/16"
13	7.65 cm	3"
14	8.09 cm	3 3/16"
15	8.2 cm	3 1/4"



OVAL-8 KIT (44 splints, size 2 thru 15)
OVAL-8 KIT FINGER SPLINT (Pkg: 5 per pack)