

# NEXTA™ SPEEDBOX S3



## NUTRIENT & MICROBIAL PLANTER BOX DELIVERY SYSTEM FOR SOYBEANS

NEXTA™ SPEEDBOX S3 is a best-in-class 80/20 talc/graphite and micronutrient blend that ensures optimum planting performance. It can replace any seed fluency agent.


The base pail includes 0.43 pounds of IONLOCK Zinc™, plus iron and manganese. IONLOCK Zinc increases singulation performance, adheres microbes to seed coats, and provides valuable early emergence nutrition.

The patented BIO-CAPSULE™ Technology delivery system allows for the addition of multiple biological solutions safely packaged for convenient, live deployment at planting.

ETHER™ Enzyme Technology kick-starts the uptake of available soil nutrients and water while maximizing microbial activity.

### PRODUCT BENEFITS

- Drives early-season root growth, nodulation and N-fixation
- Builds stronger, healthier plants fast with strong roots
- Puts soybeans in the position to manage season-long stressors and optimize genetic potential



Actual product packaging may be different than this visual

### BASE INGREDIENTS:

#### GUARANTEED ANALYSIS

Iron (Fe).....	0.70%
Manganese (Mn) .....	0.90%
Zinc (Zn) .....	1.58%
1.58% Water Soluble Zinc (Zn)	
<i>Derived from Iron Oxide, Manganese Oxide and Zinc Sulfate.</i>	

### BIO-CAPSULE INGREDIENTS:

#### ALSO CONTAINS BENEFICIAL SUBSTANCES

#### 11.71% Biological Components

<i>Azospirillum brasilense</i> .....	5.0x10 <sup>7</sup> CFU/g
<i>Azospirillum lipoferum</i> .....	5.0x10 <sup>7</sup> CFU/g
<i>Azotobacter chroococcum</i> .....	5.0x10 <sup>7</sup> CFU/g
<i>Azotobacter vinelandii</i> .....	5.0x10 <sup>7</sup> CFU/g
<i>Thiobacillus ferrooxidans</i> .....	5.0x10 <sup>7</sup> CFU/g
<i>Paenibacillus azotofixans</i> .....	5.0x10 <sup>7</sup> CFU/g
<i>Bacillus amyloliquefaciens</i> .....	5.3x10 <sup>8</sup> CFU/g
<i>Bacillus licheniformis</i> .....	1.1x10 <sup>9</sup> CFU/g
<i>Bacillus megaterium</i> .....	5.3x10 <sup>8</sup> CFU/g
<i>Bacillus pumilus</i> .....	8.8x10 <sup>8</sup> CFU/g
<i>Bacillus subtilis</i> .....	5.3x10 <sup>8</sup> CFU/g
<i>Trichoderma harzianum</i> .....	3.5x10 <sup>5</sup> CFU/g

<i>Methylobacterium hispanicum</i> .....	1x10 <sup>9</sup> CFU/g
--	-------------------------

Tannic Acid .....	0.01%
Lipase .....	2.5x10 <sup>4</sup> μUnits/g
Mannanase .....	1.25x10 <sup>9</sup> μUnits/g

#### 88.29% Total Other Ingredients (Inert)

### DIRECTIONS FOR USE & RECOMMENDED USE RATE

#### DO NOT ACTIVATE BIO-CAPSULE UNTIL READY FOR USE.

Once ready to plant, remove the safety clips from the BIO-CAPSULES. Push down on the buttons atop the pail to release the biologicals into the base. Seal the pail and shake aggressively to blend the contents. Before each application, aggressively shake contents to ensure the correct volume is applied.

The combined products (1 pail) treat 40 units of soybeans at 140k seeds per unit.

For single-row unit planters, use the enclosed scoop to measure out blended contents. One scoop treats 1 unit (140k seeds) soybeans.

### PRODUCT AVAILABILITY

Sold in 2-pail cases to treat 80 units of soybeans. Each pail is capable of treating 40 units of soybeans.

NET WEIGHT (per pail): 4.53 lbs (2.05 kg)



9090 Katy Fwy., Suite 400 Houston, TX. USA 77024  
Toll Free: 1-800-538-5283 | Tel: 1 (713) 461-1493

©2024 Stoller Enterprises, Inc. Revised 12/24

™Trademarks of Corteva Agriscience and its affiliated companies. BIO-CAPSULE™ & design, ETHER™ & design, and IONLOCK Zinc™ are trademarks of Meristem Crop Performance Group, Inc. Terrasym® is a registered trademark of NewLeaf Symbiotics®.

### NOTICE

The information contained herein is true and accurate to the best of our knowledge and belief, but is presented without guarantee since field conditions and use are beyond our control. Neither the manufacturer nor the seller make any warranty, expressed or implied, concerning the use of this product other than indicated on the label. Buyer assumes all risk of use of this material when such use is contrary to label instruction. Read and follow the label directions carefully.