

Dextrogel® is an easy to swallow dextrose gel that acts quickly. Dextrogel® is available in:

- 75 g (3x25 g) tubes with snap off closures
- 80 g soft plastic bottle

Each 25 g dose supplies a fast acting delivery of 10 g carbohydrate

For oral use only.

Not to be given to unconscious individuals.



Not recommended for children under 2 years.

Keep out of reach of children.

Do not use if product has been opened or punctured.

Store in a cool place, below 25°C. May be stored in a refrigerator. Store in the original container. Do not freeze.

Do not use after the 'Expiry date' printed on carton.



Directions for use

1 .Twist off cap and squeeze required amount inside the cheek.

2. Gently massage the outside of the cheek to aid absorption.

3. Alternatively, squeeze required amount into the mouth and swallow.

Repeat if necessary. Speak to your healthcare advisor regarding your carbohydrate requirements.

dextrogel™
FAST ACTING DEXTROSE GEL

Neoceuticals Ltd, Level 18, 40 Bank Street,
Canary Wharf, London, E14 5NR,
United Kingdom

Date of preparation May 2025

DLF03-6

Ingredients: Water, glucose monohydrate, citric acid, gelling agent (sodium carboxy methyl cellulose), preservative, E202. Contains no artificial colouring.

Nutritional information (per 100 g): Energy 150 Kcalories. Carbohydrates 37.5 g; Sugars 36.3 g; Sodium 0.1 g.