



Positive Sleep Solutions

Your CPAP Support Center

Courtesy of Positive Sleep Solutions, here's a printable luggage tag you can attach to your CPAP bag to make airline travel and TSA screening smoother and faster.

Directions:

1. Print the tag on your home printer using standard paper or cardstock for added durability.
2. Cut along the dotted line to trim the card to size.
3. Insert the card into a standard clear luggage tag holder (available at most retailers like Target, Walmart, or office supply stores).
4. Attach the tag securely to your CPAP travel bag before heading to the airport.

Trim at Dotted Lines



MEDICAL ALERT: RESPIRATORY EQUIPMENT



Name _____

Email _____

Telephone _____

Address _____

Need assistance? Call 419-534-4185 or visit www.positivesleepsolutions.com

DOT (Department of Transportation) rules state: CPAP Respiratory devices brought into the cabin by an individual with a disability shall not count toward a carry-on item. TSA (Transportation Security Administration) rules state: CPAP machines are allow through security checkpoints once the have undergone screening



MEDICAL DEVICE

CPAP / BiPAP / VPAP / Auto CPAP



Positive Sleep Solutions

This is a necessary medical respiratory device and should be kept with the patient at all times during travel. It may not be stored and does not count as carry-on baggage. This medical device provides life support for people with Obstructive Sleep Apnea, and should not be separated from the prescribed owner.