


IC30 ACTIVEX™ MATTRESS

Feel Rating	Medium
Mattress Height	200mm
Risk Category	High Risk
User Weight	30kg - 228kg
Warranty	10 Year Warranty



BED CODES & SIZES

SINGLE

Code: IC30S
Mattress Size: 1900 x 900

LONG SINGLE

Code: IC30LS
Mattress Size: 2030 x 900

KING SINGLE

Code: IC30KS
Mattress Size: 2030 x 1070

DOUBLE

Code: IC30D
Mattress Size: 1900 x 1350

LONG DOUBLE

Code: IC30LD
Mattress Size: 2030 x 1350

QUEEN

Code: IC30Q
Mattress Size: 2030 x 1520

SPLIT QUEEN

Code: IC30SQ
Mattress Size: 2030 x 750 (per side)

KEY FEATURES

- 200mm Mattress Height
- 200kg Weight Rating
- Zip off washable covers
- Anti Slip Base Layer
- Contour Latch System
- Multi layer immersion
- Compatible with 20cm Zerotec waterproof covers
- Flip free design
- Designed for use on adjustable beds
- 20+ waterlow rating
- 10 Year Warranty

SPECIFICATIONS

- ActiveX™ pressure relieving material with Charcoal and Gel infusion
- Vant ventilated sides and base fabric
- Bluetec temperature regulating cover
- Built in firm edge support
- Flip free design
- Zip off machine washable covers
- 30 - 200kg user weight
- Standard cover MVTR: NA.
- Medical Cover MVTR (g/24h·m2): 1635.0 Average
- Medical Cover BS7175: Crib 5
- ARTG Identifier: 259893

PRESSURE CARE

Risk Category High

This rating is based on our average peak pressure of <19mmHg. For some high risk clients Icare Zerotec covers are recommended.

Icare risk category ratings are based on average peak pressure, Therapist recommendation, case studies, industry comparison and the Waterlow scale. For some high risk clients Icare Zerotec covers are recommended. The best way to assess a client on a pressure reducing surface is to monitor pressure points with visual checks over a period of 5 days.

CARE & CLEANING

To launder the cover in a washing machine, set the main wash to a minimum of 65°C for no less than 10 minutes. The cover can be tumble dried on low heat and the temperature must not exceed 60°C. The foam can be autoclaved at 134°C.

